



# Commode Chair

## Purpose

Used in the bedroom as a substitute for a toilet, especially used at night time or when a person has difficulty getting to the toilet.

## Features


- Chair with toilet seat and bowl
- Arm rests assist with getting on and off
- The height of the commode can be adjusted to make it easier to sit and stand
- Non slip feet
- Splash guard which prevents spillage

appearance of the rubber stoppers, or damage/bending in the metal. If any of these are present, stop use and contact your supplier.

- This product should be cleaned regularly with warm soapy water unless otherwise stated by the manufacturer.


## Installation

- Ensure there is sufficient space for an easy transfer from the bed at night
- Adjust frame to the appropriate height by depressing the silver button and sliding the leg pieces up or down
- Ensure all four legs are in contact with the floor to avoid rocking or unsteadiness



## Safety and precautions

- Use two hands when transferring on and off the commode to avoid tipping
- This equipment is for sitting only. Do not stand on this equipment
- You should only use the equipment in the environment in which it is prescribed.
- 100kg weight limit unless otherwise stated by the manufacturer
- Ensure the silver buttons on the height adjustable legs have slotted into the correct holes before sitting



## Maintenance

- Wipe chair over after use to remove moisture and prolong the life of the chair
- The frame should be checked monthly for signs of rust, cracks, changes in the

Your therapist has recommended the appropriate height for your commode as:

## Other comments

All information contained in this sheet is a guideline for care only, supplied by qualified health professionals. Acknowledgements to Toowoomba Health Service District, Queensland Health

## Contact Us.

P 1800 684 277

E info@maxhealth.com.au

W maxhealthcareequipment.com.au