

MONARCH RECLINER

single motor / dual motor / petite



- * Seamed backrest provides additional lumbar support
- * Padded legrest is synchronised with the backrest (Single Motor)
- * Independent control of backrest and legrest (Dual Motor)
- * Full chaise pad provides additional legrest support and comfort
- * Supplied with magazine pocket and battery back-up
- * Hand held controller is simple to use
- * Rise or recline motion stops at any point when button is released

Size	Type	Seat Width	Seat Depth	Seat Height	Arm Height	Back Height	Overall Width	Max User Weight
Standard	Single Motor	550mm	460mm	500mm	190mm	680mm	820mm	130kg
Petite**	Single Motor	480mm	400mm	460mm	140mm	620mm	780mm	130kg
Standard	Dual Motor	550mm	460mm	500mm	190mm	680mm	820mm	130kg

** The Petite is available in Pacific Blue fabric and Coffee vinyl ONLY.