



RDNS EQUIPMENT CATALOGUE



EXTERNAL PRESCRIBERS

Ver 3.31

	BATHING AND TOILETING
	BEDROOM
	BRACES AND SUPPORTS
	CHAIRS AND SEATING
	CONTINENCE
	DAILY LIVING AIDS
	ELECTRIC BEDS
	HOME ACCESS
	LIFTERS AND HOISTS
	MOBILITY AIDS
	MOBILITY AIDS ACCESSORIES
	PRESSURE CARE
	REHABILITATION EQUIPMENT
	TRANSFER AIDS
	WHEELCHAIRS



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
BATHING AND TOILETING				
Bath Aids				
BTBDWH	Bathboard with handle	Provides a support surface to place across the bath which enables safe and easy showering in a seated position. Helps when transferring in and out of the bath as this can be performed in a seated position. The included handle assists with transfers in and out of the bath as well as safety whilst seated showering.	Adjustable width clamps: 41 cm - 62 cm Length: 69 cm Width: 26 cm Maximum user weight capacity: 150 kg	
BTBDBA	Bathboard - Bariatric	Designed to provide a stable seat over the bath. Constructed of reinforced PVC Anti - slip rubber pads improves safety and protects the bath surface	Adjustable width clamps: 50 cm - 66 cm Length: 66 cm Width: 34 cm Maximum user weight capacity: 190 kg	

BTBSSW	Bath Seat - Swivel	Designed for showering over bath and aiding with transfers into the bath by allowing this occur in a safe seated position.	<p>Maximum user weight: 110kg</p> <p>Adjustable width to fit most baths</p> <p>Moulded plastic seat with drainage holes</p> <p>Contoured back rest</p> <p>Side support arms for added security</p> <p>360° swivel seat</p> <p>Locking lever for set positioning</p>	
BTBTBW	Bath Transfer Bench - With Back	<p>Provides a support surface to place across the bath which enables safe and easy showering in a seated position. Helps when transferring in and out of the bath as this can be performed in a seated position.</p> <p>Plastic seat and backrest provides a non slip surface with holes for water drainage</p> <p>Adjustable backrest and handle to suit left or right facing baths without using tools</p> <p>Built-in dish and hand-held shower slot for user convenience</p>	<p>Height adjustable seat: 43 cm - 56 cm (25mm increments) Seat dimension: Width 66 cm Depth 37 cm Maximum user weight: 136kg</p>	

BTSMNS

Shower Mat / Bath Mat Non Slip With their raised ridges and unique design, the mats allow you to stand, without slipping, on a wet surface. Use on flat firm wet surface & press down first so the suction caps engage. The mats grip securely and safely due to their high-grade medical rubber content and additional suction cups.

Machine washable at 40°.
Dimension: 555mm x 555mm



Mobile Shower Commodes

BTMSTR6

Mobile Shower Commode,
Transit, 457mm/18inch Wide,
Aluminium

The attendant propelled mobile shower commode is lightweight and maneuverable to allow for easier toileting and showering whilst the user is seated. The chair features a lightweight, silver painted aluminium frame, which provides a rust free, high strength product. Brakes on all four wheels for security, the swing away, detachable footrests offer comfort and support whilst being transported and the armrests can be flipped back to facilitate side transfer. This chair can be used in a shower cubicle or wet room, used as a commode or over a standard toilet bowl. The front of the seat is cut-away to facilitate personal hygiene and a set insert is provided for comfort during transport. Supplied with a square commode pan as standard. The detachable back enables the shower commode chairs to be partially dismantled for more compact storage.

Overall Height: 958 – 1008mm (with castors)
Overall Length: 550mm (without footrests) / 810mm (with footrests)
Overall Width: 555mm
Seat Height: 495 – 570mm
Seat Width: 457mm
Seat Depth: 425mm
Aperture Size: 190mm x 360mm
Backrest Width: 462mm
Backrest Height (from seat): 436mm
Front & Rear Wheel Diameter: 130mm
Under Seat Clearance: 430 – 505mm
Weight: 12.4kg
Maximum User Weight: 125kg
Warranty: 12 months



BTMSTR1

Mobile Shower Commode,
Transit, Stainless Steel, 460mm

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Swingback arms with padded armpads assist side transfers and can be positioned at two different heights. Adjustable safety arms for additional security. Sliding footrest is textured and weight-bearing for added safety. Footrest slides back under frame to assist patient transfers. PVC slatted backrest. Mounted on 4 x 125mm stainless steel locking castors. Durable stainless steel frame. Tested and conforms to Australian Standards AS/NZS 3973: 2009I

Seat Width - 460mm
Seat Depth - 430mm
Seat Height - 595mm
Armrest Height - 200/240mm
Backrest Height - 340mm
Overall Width - 540mm
Max User Weight - 200kg

**BTMSTR3**

Etac Swift Mobile Shower
Commode, Transit, Swing Away
Footplates, Adjustable Height

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Soft, adjustable, non-stick, quick-drying nylon back support. Seat and back support easily detachable for cleaning or storage. Swing away arm and foot supports. Curved foot supports with side heel protection.

Fully height adjustable
from 450 - 600mm
Seat height (at 20mm
intervals)
Maximum user weight
135kgs



BTMSTR4

Mobile Shower Commode,
Transit, Bariatric, 560mm wide,
Stainless Steel

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove.

Wider seat widths and reinforced frame to suit the bariatric client

Swingback arms with padded armpads assist side transfers and can be positioned at two different heights

Adjustable safety arms for additional security

Sliding footrest is textured and weight-bearing for added safety

Footrest slides back under frame to assist patient transfers

PVC slatted backrest

Mounted on 4 x 125mm stainless steel locking castors

Durable stainless steel frame

Tested and conforms to Australian Standards AS/NZS 3973: 2009I

Seat Width - 560mm

Seat Depth - 430mm

Seat Height - 595mm

Armrest Height -
200/240mm

Backrest Height - 340mm

Overall Width - 640mm

Max User Weight - 300kg



BTMSTR5

Mobile Shower Commode,
Transit, Bariatric 610mm wide,
Stainless Steel

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Wider seat widths and reinforced frame to suit the bariatric client

Swingback arms with padded armpads assist side transfers and can be positioned at two different heights

Adjustable safety arms for additional security

Sliding footrest is textured and weight-bearing for added safety

Footrest slides back under frame to assist patient transfers

PVC slatted backrest

Mounted on 4 x 125mm stainless steel locking castors

Durable stainless steel frame

Tested and conforms to Australian Standards AS/NZS 3973: 2009I

Seat Width - 610mm
Seat Depth - 430mm
Seat Height - 595mm
Armrest Height -
200/240mm
Backrest Height - 340mm
Overall Width - 690mm
Max User Weight - 300kg



BTMSSP1

Mobile Shower Commode, Self
Propelled, SS, 460mm

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Swingback arms with padded arm pads for side transfers
Adjustable safety arms for additional security
Armrests can be positioned at two different heights
Removable PVC slatted backrest
Adjustable height swingaway legrests or slide out weight-bearing footrest
Mounted on 2 x 125mm stainless steel castors and 2 x 22 inch rear wheels
Durable stainless steel frame.

Seat Width - 460mm Seat
Depth - 430mm Seat
Height - 595mm Armrest
Height - 200/240mm
Backrest Height - 340mm
Overall Width - 660mm
Max User Weight - 200kg



BTMSSP2

Mobile Shower Commode, Self
Propelled, Adjustable Height

Easy to maneuver, for the independent user, the 24-inch wheels are positioned to make it easy to maneuver the chair without compromising user safety or comfort. It is easy to transfer to and from the chair thanks to the large, unrestricted seat in front of the wheels. The brake handles are also positioned below the seat so that they do not get in the way during transfer, while remaining easy to reach. Front seat recess; fully open at rear Long, removable, swing away arm supports Curved footplate with side protection for heels; slides in under seat Footplate features full height and tilt adjustment.

Width - 69 cm
Width between arm
supports - 43.5 cm
Seat height - 55 cm
Seat width - 48 cm
Seat depth - 44 cm
Space below seat - without
pan holder: 49.5 cm, with
pan holder: 48 cm Opening
- depth x width, 36 cm x
21,5 cm
Weight - 20 kg
Maximum user weight -
130 kg



Shower Chairs

BTSCST

Shower Chair, Standard,
Adjustable - With Arms

A standard no-fuss, safe and reliable chair which is height adjustable. It features a contoured, plastic moulded seat, which is both comfortable & hygienic. It has been cleverly designed with a gentle slope, so that water can flow away for easy drainage. The steel frame using 22.2mm tubing enables users to gain a firm grip making transferring safer and easier. Tested to Australian standards.

Overall Height : 775mm - 900mm
Overall Width : 545mm
Overall Depth : 515mm
Seat Width: 440mm
Seat Depth: 405mm
Seat Height Range: 410mm–535mm
Width between Armrests: 490mm
Armrest (Seat to top of handle): 200mm
Backrest Height: 350mm
Frame: Steel
Frame Diameter (armrest and leg): 22.2mm
Extension Leg Diameter: 25.8mm
Extension Leg (outer/inner): Outer
Rubber Tip (White): DIA 40mm x DIA 23.8mm
Colour: White
Net Weight: 3.7kg
Maximum User Weight: 125kg
Warranty: 12 months



BTSCB00

Shower Chair - Bariatric - 525
Wide

Reinforced wide frame – Provides stability for larger users
Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
Height adjustable legs – Easy push button height adjustment for changing needs
Additional height increments on one leg – Accommodates uneven surfaces
All frames are individually treated for corrosion resistance
Raised armrests – Assist the user transferring on and off the elevated seat
Durable non-slip feet – Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability

SWL 200kg Overall
Height: 735 - 910 mm
Overall Width: 570 mm
Overall Depth: 570 mm
Width Between Arms: 525 mm

**BTSCB03**

Shower Chair - Bariatric - 550
Wide

Height adjustable
Wider seat widths and reinforced frame to suit the bariatric client
Padded seat and curved backrest provide comfort and postural support
Non-slip rubber feet
Powder coated steel frame
Electrophoretic surface treatment to prevent rusting

Seat Width 550mm
Seat Height 470 - 640mm
Base Width 650mm
Base Depth 520mm
Max User weight 300kg



BTSCB04

Shower Chair - Bariatric - 600 Wide

Height adjustable
Wider seat widths and reinforced frame to suit the bariatric client
Padded seat and curved backrest provide comfort and postural support
Non-slip rubber feet
Powder coated steel frame
Electrophoretic surface treatment to prevent rusting

Seat Width 600mm
Seat Height 470 - 640mm
Base Width 700mm
Base Depth 520mm
Max User weight 300kg



BTSCB01

Shower Chair - Bariatric - 610 Wide

Backrest can be adjusted in depth +/- 60mm and has curved tubing to ensure comfort for any body shape.
Having the right height for the user ensures success in getting in and out of the chair and these chairs are adjustable from 420 to 550mm.
Butterfly armrests ensure the body doesn't get squeezed getting in to and out of the chair whilst providing a good grip for sitting and getting up.
The high-tech, lightweight tubing provides the highest strength at the lowest weight
Removable seat pad ensures these items are easy to clean.

Seat Width: 610mm Seat Depth: 400 to 520mm
Overall Width: 730mm
Weight: 7kg Maximum User Weight: 325kg



BTSCB05

Shower Chair - Bariatric - 650 Wide

Height adjustable
Wider seat widths and reinforced frame to suit the bariatric client
Padded seat and curved backrest provide comfort and postural support
Non-slip rubber feet
Powder coated steel frame
Electrophoretic surface treatment to prevent rusting

Seat Width 650mm
Seat Height 470-640mm
Base Width 750mm
Base Depth 520mm
Max User Weight 300kg




**BTSCB02**

Shower Chair, Commode & Toilet Aid, Bariatric (All in One)

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Total Width: 790mm
Total Depth: 530mm
Seat Depth: 450mm
Seat Height: 465 - 615mm (to over seat top)
Width Between Arms: 600mm
Height Adjustable: 430-570mm
Safe Working Limit: 295kg



BTSCSW1	Shower Chair, Swivel	<p>Height adjustable</p> <p>Swivel motion assists easy transfers from standing to showering position</p> <p>Seat locks in both transfer and showering positions</p> <p>Integrated arms assist transfers</p> <p>Drainage slots drain water away from seating area</p> <p>Non-slip rubber tips</p> <p>Powder coated steel frame</p>	<p>Width Between Arms 480mm</p> <p>Seat Height 480-560mm</p> <p>Base Width 560mm</p> <p>Base Depth 540mm</p> <p>Max User Weight 160kg</p>	
	Shower Stools			
	BTSSWA Shower Stool Adjustable - With Arms	<p>The standard shower stool is height adjustable and features side handles to assist with transfers. The stool has a small footprint and is ideal for smaller shower alcoves.</p>	<p>Overall Width : 510mm</p> <p>Overall Depth : 410mm</p> <p>Width Between Arms: 450mm</p> <p>Weight of Stool: 3.2kg</p> <p>Height Adjustable: 490 - 610mm</p> <p>Maximum User Weight: 110kg</p>	
BTSSNA	Shower Stool Adjustable - Without Arms	<p>Height adjustable</p> <p>Side handles to assist transfers</p> <p>Durable plastic moulded seat</p> <p>Non-slip rubber tips</p> <p>Electrophoretic surface treatment to prevent rusting</p> <p>Powder coated steel frame</p>	<p>Width 380mm</p> <p>Seat height 510 - 660mm</p> <p>Base width 550mm</p> <p>Base Depth 380mm</p> <p>Max user weight 125kg</p>	

BTSSAW2

Shower Stool Adjustable - With Backrest

Height adjustable Padded seat and back for added comfort Integrated arms assist transfers Non-slip rubber tips Electrophoretic surface treatment to prevent rusting Powder coated steel frame

Width Between arms 460mm Seat Height 520-670mm Base Width 550mm Base Depth 380mm Max User Weight 125kg





BTSSAHD




Shower Stool Adjustable -
Heavy Duty Padded - Aluminium




Constructed from sturdy no rust
aluminium with easy push button
adjustment legs
The internal sliding legs provide
smooth and effortless height
adjustment whilst avoiding any chance
of water filling the legs
Reinforced frame using 25.4mm tubing
enables user to gain a firmer grip
making transferring safer & easier
The seat is padded for comfort and is
removable for cleaning

Overall Height (footprint):
679mm - 829mm
Overall Width (footprint):
655mm
Overall Depth (footprint):
445mm
Seat Width: 445mm
Seat Depth: 275mm
Seat Height Range:
539mm–689mm
Width between Armrests:
538mm
Armrest (Seat to top of
handle): 140mm
Frame: Aluminium
Frame Diameter (armrest
and leg): 25.4mm
Extension Leg Diameter:
28.6mm
Extension Leg
(outer/inner): Outer
Rubber Tip (White): DIA
40mm x DIA 27mm
Colour: White
Net Weight: 3.0kg
Maximum User Weight:
150kg
Warranty: 12 months
ARTG ID: 291268



BTSSAWB	Shower Stool With Removable Back and Arms - Bariatric	<p>Aluminum construction in light weight and corrosion free</p> <p>Cross-brace to enhance the stability during sitting</p> <p>Blow molded bench and back provide comfort and strength</p> <p>Drainage in seat to reduce slipping</p> <p>Limited Lifetime Warranty</p>	<p>Seat width 51cm</p> <p>Seat depth 36cm</p> <p>Seat height 36-51cm</p> <p>Overall width 58cm</p> <p>Overall depth 38cm</p> <p>Overall height 68-84cm</p> <p>Weight Capacity 227kg</p> <p>Total weight 3.3kg</p>	
	Accessories			
	BTSH1T	Shower Hose Single Tap 2 Metre - Push on	<p>Fixes to existing tap spouts Lightweight plastic handpiece Supplied with two rubber reducing collars to suit most tap spouts Easily fitted and removed</p> <p>Length 2m</p>	
BTSH2T	Shower Hose Double Tap 2 Metre - Push on	<p>Hand Shower – Dual Tap Clamp-on. Fixes to existing tap spouts Lightweight plastic handpiece Supplied with two rubber reducing collars to suit most tap spouts Easily fitted and removed</p> <p>Length 2m</p>		
Toileting Aids				

BTTAX01	Toilet Seat Extension - With Lid 50mm	- Ergonomic contoured one piece seat with lid Tapered flange with front and rear access Stain and odour resistant Easy to secure with two hand adjustable brackets	Height 50mm. Maximum user weight 130kg	
BTTAX02	Toilet Seat Extension - With Lid 100mm	- Ergonomic contoured one piece seat with lid Tapered flange with front and rear access Stain and odour resistant Easy to secure with two hand adjustable brackets	Height 100mm. Maximum user weight 130kg	
BTTAX03	Toilet Seat Extension - 50mm	Increases toilet seat height Ideal for people with hip or knee injuries or those who require greater seat height Easily positioned and removed using location lip and adjustable clamps Front and rear cutouts for personal cleaning Contoured seat for greater comfort Moulded plastic is stain and odour resistant	Height: 50mm Max user weight: 190kg	

BTTAX04	Toilet Seat Extension - 100mm	Increases toilet seat height Ideal for people with hip or knee injuries or those who require greater seat height Easily positioned and removed using location lip and adjustable clamps Front and rear cutouts for personal cleaning Contoured seat for greater comfort Moulded plastic is stain and odour resistant	Height: 100mm Max user weight: 190kg	
BTTSA01	Toilet Surround - Adjustable	Frame to assist toilet transfers Height adjustable Non-slip rubber tips Powder coated steel	Height - 625-775 Width between Arms - 455 Base Width - 540 Base Depth - 490 Max User Weight - 125kg	
BTTSW01	Toilet Surround - Wide	Heavy duty fixed height frame to assist toilet transfers Fitted with rear cross bar for additional strength and stability Non-slip rubber tips Powder coated steel	Height - 685 Width between Arms - 620 Base Width - 670 Base Depth - 520 Max User Weight - 130kg	

BTTSRST

Toilet Seat Raiser - Adjustable
with Splash Guard (SWL 125kg)

The Over Toilet Aid assists transfers on & off the toilet. This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non slip rubber feet. Comes complete with a seat and splashguard

Seat Height 440-590mm
Seat Width 370mm
Width Between Arms 450mm
Max User Weight 125kg

**BTTSRTW**

Toilet Seat Raiser - Adjustable
with Splash Guard (SWL 175kg)

Raised armrests – Assist the user transferring on and off the Over Toilet Aid
Additional height increments on one leg – Accommodates uneven surfaces
Height adjustable legs – Easy push button height adjustment for changing needs

Seat Height 450 - 630mm
Seat Depth 435mm
Width Between Arms 460mm
Max User Weight 175kg

**BTTSR55**

Toilet Seat Raiser - Bariatric
550mm Wide

Provides a raised seat to assist in toileting and transfers.
Wider seat widths and reinforced frame to suit the bariatric client.
Ideal for people with hip or knee injuries or those who require greater seat height
Ergonomic angled arms assist transfers
Supplied with lift-back toilet seat and removable splashguard
Non-slip rubber tips
Powder coated steel

Width Between Arms 550mm
Seat height 550mm
Armrest height 230mm
Base width 680mm
Base Depth 610mm
Maximum user weight 200kg



BTTSR60

Toilet Seat Raiser - Bariatric
600mm Wide

Provides a raised seat to assist in
toileting and transfers
Wider seat widths and reinforced frame
to suit the bariatric client
Ideal for people with hip or knee
injuries or those who require greater
seat height
Ergonomic angled arms assist transfers
Supplied with lift-back toilet seat and
removable splashguard
Non-slip rubber tips
Powder coated steel

Width Between Arms
600mm
Seat height 550mm
Armrest height 230mm
Base width 730mm
Base Depth 610mm
Maximum user weight
250kg

**BTTSR65**

Toilet Seat Raiser - Bariatric
650mm Wide

Provides a raised seat to assist in
toileting and transfers
Wider seat widths and reinforced frame
to suit the bariatric client
Ideal for people with hip
or knee injuries or those who require
greater seat height
Ergonomic angled
arms assist transfers
Supplied with lift-
back toilet seat and removable
splashguard
Non-slip rubber tips
Powder coated steel

Width Between Arms
650mm
Seat height
550mm
Armrest height
230mm
Base width 780mm
Base Depth 610mm
Maximum user weight
300kg







BTSCB02**Shower Chair, Commode &
Toilet Aid, Bariatric (All in One)**

The bariatric All in One provides for being a toilet seat raise, shower chair and commode chair. The legs are height adjustable. It has a square bowl with easy access for removal from the front of the chair. The square bowl does not have a lid. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Total Width: 790mm
Total Depth: 530mm
Seat Depth: 450mm
Seat Height: 465 - 615mm
(to over seat top)
Width Between Arms:
600mm
Height Adjustable: 430-
570mm
Safe Working Limit: 295kg



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
BEDROOM				
Bedside Commode				
BRCMAH1	Commode - Adjustable Height	The Bedside Commode has a clip on toilet seat, which clips on in seconds. There are no metal fittings on toilet seat. It is also height adjustable and easy to clean. It has a light blue deluxe padded overseat, as well as a removable round bowl with lid and handle. The lid eliminate spillages and reduces risks of cross infection. It also features a powder coated dark blue frame. This bedside commode looks very smart in a bedroom.	Overall Height: 820 - 970mm Overall Width: 610mm Overall Length: 550mm Seat Height: 460 - 610mm Weight: 8kg Maximum User Weight: 110kg	
BRCMCL	Classic Bedside Commode	Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom Available in Champagne or Slate Vinyl Raised Padded Armrests Easy Clean Vinyl Surfaces Height Adjustable Removable, Padded Seat Non-Slip Rubber Tips	Seat Height: 475 - 615 mm	

BRCM58	Commode - Adjustable Height - Heavy Duty - 580mm	Steel Heavy Duty Commode With removable Backrest. Durable and rigid frame in extra wide for high weight loading Adjustable height.	Seat 58cm width x 46cm depth Max User Weight 200kg.	
BRCMB65	Commode - Bariatric - 650mm Wide	Wider seat widths and reinforced frame to suit the bariatric client Padded seat, padded armrests, lift-back toilet seat and removable commode bowl Adjustable seat height Padded seat cover slides back easily for access to commode Three sided vanity skirt conceals commode bowl Powder coated steel frame	Width between arms 650mm Seat height: 390 - 540mm Armrest Height: 300mm Max user weight 150kg	

BTSCB02

Shower Chair, Commode &
Toilet Aid, Bariatric (All in One)

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Total Width: 790mm
Total Depth: 530mm
Seat Depth: 450mm
Seat Height: 465 -
615mm (to over seat top)
Width Between Arms:
600mm
Height Adjustable: 430-
570mm
Safe Working Limit:
295kg



Bed Sticks

BRHBR01

Horizontal Bed Rail

This Bed Rail is designed to prevent the user from rolling out of bed. Once the user is in bed, simply slide the Bed Rail under the mattress. The weight of the user secures the Bed Rail. Constructed from sturdy epoxy coated steel.

Height: 500mm
Width (Base): 580mm
Length (Rail): 600mm
Weight: 2.5kg



Bed Positioning

BRBB01**Bed Bolster - Mattress Extension**

200 mm long foam mattress extension for additional pressure relief under heel zone, available in Standard and King Single widths.

Suitable for use in most standard residential care, hospital and home care beds with sleeping platform extension capabilities. Standard model suits any 150 mm high foam mattress; King Single model suits Triple Layer King surface.

King Single (180 high, 1050mm wide)

Single (150 high, 865mm wide)

Features

high-density memory foam upper layer with medium density foam core base
double-strength multi-stretch waterproof cover reduces the potential for shear injuries
seams are stitched, welded and taped for maximum infection control
Standard and King Single widths available

**BRBRA01****Back Rest - Adjustable**




Multi-angled pillow support
Provides assistance and support in bed
Wide aluminium frame fitted with strong breathable nylon back
Angle can be adjusted via a simple ratchet mechanism

Length 495mm

Width 700mm

Height 490mm



BRBW01	Foam Bed Wedge	<p>Angled foam cushion provides comfort and support to the upper body Ideal for those who suffer from respiratory and reflux problems Removable quilted cotton cover Wedge can be turned to provide two different recline angles (15 degrees or 75 degrees) Foam insert can be removed to alter the angle by approx 5 degrees</p>	<p>Length 680mm Width 710mm Height 240-600mm</p>	
	Bed / Chair Blocks			
	BRBK01 Bed/Chair Blocks Plastic 40mm	<p>Plastic block raisers designed to increase the height of a bed or chair. Increases height by 40 mm Stackable</p>	<p>Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg</p>	
BRBK02	Bed/Chair Blocks Plastic 100mm	<p>Plastic block raisers designed to increase the height of a bed or chair. Increases height by 100 mm Stackable</p>	<p>Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg</p>	
Overbed Tables				

BROBT01

Overbed Table

Eoxy coated frame with a vinyl top. Recessed edges, which enables the table to be cleaned easily. The 50mm twin wheel castors allow the table to be easily moved away from the bed when not in use.

Length 800mm
Overall Width 400mm
Minimum Height 710mm
Maximum Height 1050mm
Height of Castor Wheels 50mm



Self Help Poles/Ladders



BRBCA01


Bed Cradle - Adjustable Height

The bed cradle is designed to keep bed linen off sensitive areas of the patient's lower leg or feet by creating clearance between bed linen and the mattress. The lower frame slides underneath the mattress, securing the cradle to the bed, while the upper frame is positioned horizontally over the lower limbs. Bedding is then placed over the frame to remove any weight or pressure from the sensitive limbs below. Height adjustable frame allows the user to customise the space between mattress surface and bedding.

Height adjustable 440 - 640 mm



BRBL01	Bed Ladder	<p>A solution for clients who have physical difficulties sitting up in or getting out of bed</p> <p>Step by step they gradually pull themselves up from a lying to a sitting position</p> <p>The action involved also strengthens the stomach and arm muscles of the user</p> <p>Can also be used to assist in transfers from wheelchairs to beds</p> <p>Plastic ladder rungs give a firm and safe grip</p> <p>Cord ends tie around the bed legs</p> <p>Suitable for single or double beds</p>	<p>Rung Width 195mm</p> <p>Rung Diameter 25mm</p> <p>Cord Length 3m (to first rung)</p>	
BRSHPO1	Self Help Pole - Freestanding	<p>Free standing support handle provides assistance for the user to reposition or transfer in and out of bed</p> <p>Handle height is easily adjusted to suit the user</p> <p>Sits behind or beside the bed</p>	<p>Safe work load 150kg</p>	

ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
BRACES AND SUPPORTS				
BRACES				
BRSKB01	Knee Brace - Futuro Sport Adjustable	<p>Walk. Jog. Run. The FUTURO Sport Adjustable Knee Support keeps you moving with compression, strength and support where you need it. It also helps provide protection on the track, court, or green. With this adjustable knee wrap, there is no excuse not to be active anymore.</p> <p>Special Features: Adjustable straps for customised fit and support Use For: Sprains, strains, arthritis Product Type: Wrap Around Features and Benefits:</p> <p>Open knee cap design for added support Side ventilation for breathability and reduced bunching behind knee Durable neoprene-blend material for strong support, soothing warmth, and comfort Easy to put on and take off</p>	<p>One Size Fits Most, Measure knee circumference: 33cm - 45.7cm</p>	

BRSW01

Wrist Braces Roylan D Ring Small
Left

Wrist Circumference
14.6cm to 16.5cm
Length 16.5cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW02

Wrist Braces Roylan D Ring Small
Right

Wrist Circumference
14.6cm to 16.5cm
Length 16.5cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW03

Wrist Braces Roylan D Ring
Medium Left

Wrist Circumference
17.1cm to 19.1cm Length
17.8cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW04

Wrist Braces Roylan D Ring
Medium Right

Wrist Circumference
17.1cm to 19.1cm Length
17.8cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW05

Wrist Braces Roylan D Ring Large
Left

Wrist Circumference
19.7cm to 21.6 cm
Length 19cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW06

Wrist Braces Roylan D Ring Large
Right

Wrist Circumference
19.7cm to 21.6 cm
Length 19cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



ORTHOSSES

BRSAF01a

Rebound Foot Up Ankle Cuff S/M

Rebound Foot-Up is a lightweight ankle foot orthosis that offers dynamic and discreet support for drop foot or related disorders requiring dorsiflexion support

Anke Cuff 18-25cm Foot wrap, strap and bolt 20-28cm

Contoured to the shape of the leg, it is comfortable to wear and easy to use

Rebound Foot-Up provides visible improvement in gait by supporting the foot the moment it is raised



BRSAF02a

Rebound Foot Up Ankle Cuff L/XL

Rebound Foot-Up is a lightweight ankle foot orthosis that offers dynamic and discreet support for drop foot or related disorders requiring dorsiflexion support

Anke Cuff 26-33cm Foot wrap, strap and bolt 20-28cm

Contoured to the shape of the leg, it is comfortable to wear and easy to use

Rebound Foot-Up provides visible improvement in gait by supporting the foot the moment it is raised



BRSAF06

Ankle Foot Orthosis Foot-Up
Shoeless Wrap Universal



Rebound Foot-Up Foot Wrap can be used when barefoot or in open shoes or sandals. The adjustable elastic strap of the Rebound Foot Up is attached to the Foot Wrap using a nut and bolt fastened through the eyelets.

The Rebound Foot-Up Ankle Cuff and Rebound Foot-Up Foot Wrap are made from breathable material that can be worn comfortably for long periods of time.

The Rebound Foot-Up is hand-washable.

One size cut to size (20-28cm)



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
CHAIRS AND SEATING				
CSFS01	Foot Stools			
	Foot Stool / Leg Rest	The Legrest Stool is height adjustable, designed with a 10 degree angle for comfort and support and should be placed with the higher end against the chair.	Seat Pad 400mm x 500mm x 40mm Height Adjustable 420mm - 570mm Maximum User Weight 160kg	
	Kitchen Stools			
CSKSAWA	Kitchen Stool - Adjustable With Arms	Height adjustable Constructed with lightweight aluminium Will not rust so can be used in the shower, as well as your homecare environment A broad sloping padded white vinyl seat for maximum comfort Large non-slip feet White plastic seat and back.	Width Between Arms: 455mm Seat Depth: 300mm Seat to Top of Backrest: 300mm Leg Height Adjustment: 500mm Floor to Seat: 635mm Maximum User Weight: Up to 100kg	
Utility Chairs				

CSLBCWA**Low Back Utility Chair**

The Standard Utility Chair features a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder coated steel frame. The adjustable height assists users who struggle to get out of standard chairs.

Overall Width 580mm
Seat Width 455mm
Seat Depth 460mm
Seat Height 450 - 600mm
Back Height 465mm
Maximum User Weight 160kg

**CSHBCWA****High Back Utility Chair**

This chair can be used in a lounge room instead of a couch or regular armchair that is too low to get in and out of with ease. Features a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder coated steel frame. Available in Teal/Blue and Mocha.

Overall Width 640mm
Seat Width 490mm
Seat Depth 480mm
Seat Height 455 - 605mm
Back Height 730mm
Maximum User Weight 160kg

**CSHBCB****High Back Utility Chair - Bariatric**

This chair can be used in a lounge room instead of a couch or regular armchair that is too low to get in and out of with ease.
Wider seat widths and reinforced frames to suit the bariatric client
Adjustable height and seat depth
Contoured backrest
Powder coated steel frame

Upholstery Champagne Vinyl
Seat Width 600mm
Seat Depth 490-550mm
Seat Height 430-580mm
Arm Height 250mm
Back Height 740mm
Overall Width 800mm
Max User Weight 300kg



CSAPCL

Adult Postural Chairs

Chair Air Deluxe Bed Adult with
Footrest - Large

Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest, tray, neck support and footplate. The frame is completely constructed from sturdy stainless steel.

Features:

Adjustable calf rest

Reclinable seat

Air system for pressure relief

Full stainless steel frame

Three lockable castors and one steering castor (non-marking)

Breathable polyurethane cover (fire retardant)

Replaceable and easy to clean cover

Four folding wings

Standard dining tray

Magazine pocket

Seat Height: 63cm

Seat Width: 51cm

Seat Depth: 50cm

Overall Width: 73cm

Overall Length: 135 - 188cm

Table to Floor: 88cm

Table Slide: 23cm

Table Size: 60 x 40.5 x 2cm

Back Reclining: 102° - 136°

Seat Tilting: 0 - 18°

Leg Elevating: 70° down to 5° up

Weight of Chair: 43.5kg

Maximum User Weight: 120kg

Warranty: 1 year



Electric Recliner Chairs

CSPLRA1

Monarch Chair (SINGLE MOTOR)

Seamed backrest provides additional lumbar support
Padded legrest is synchronised with the backrest
Full chaise pad provides additional legrest support and comfort
Supplied with magazine pocket
Hand held controller is simple to use
Rise or recline motion stops at any point when button is released

Standard
Seat Width 550mm
Seat Depth 460mm
Seat Height 500mm
Arm Height 190mm
Back Height 680mm
Overall Width 820mm
Max User Weight 130Kg



CSPLRA2

Monarch Chair (DUAL MOTOR)

Seamed backrest provides additional lumbar support
Padded legrest is synchronised with the backrest
Full chaise pad provides additional legrest support and comfort
Supplied with magazine pocket
Hand held controller is simple to use
The dual motor provides independent control of backrest and legrest.

Standard
Seat Width 550mm
Seat Depth 460mm
Seat Height 500mm
Arm Height 190mm
Back Height 680mm
Overall Width 820mm
Max User Weight 130Kg



CSSMMFB

Mayfair Chair - Black (SINGLE MOTOR)

The Select Mayfair Chair features a very comfortable 3 Section Waterfall Backrest Design for Optimum Comfort and Upper, Middle and Lower back support.

Quality Foams and no Sag springs provide supreme Postural Support.

Backup Power Supply.

Featuring quality Actuators and American Branded Mechanisms means quiet motors for a Smooth & Discreet Operation.

Standard

Seat Width 520mm

Seat Depth 480mm

Seat Height 480mm

Overall Width 850mm

Max User Weight 150Kg



CSSMMFV

Mayfair Chair - Vanilla (SINGLE MOTOR)

The Select Mayfair Chair features a very comfortable 3 Section Waterfall Backrest Design for Optimum Comfort and Upper, Middle and Lower back support.

Quality Foams and no Sag springs provide supreme Postural Support.

Backup Power Supply.

Featuring quality Actuators and American Branded Mechanisms means quiet motors for a Smooth & Discreet Operation.

Standard

Seat Width 520mm




Seat Depth 480mm




Seat Height 480mm





Overall Width 850mm

Max User Weight 150Kg



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
CONTINENCE				
Bed Pans				
CNBPSL1	Bed Pan - Slipper	Assists toileting of user whilst confined to bed Wedge shape assists in correct positioning Poly Propylene construction		
BTBPP	Bed Pan Plastic	Light weight polypropylene bed pan for toileting while remaining in bed Wide opening cut-out area to allow for correct positioning Durable moulded plastic design with autoclavable construction Compatible with most Mobile Shower Commodes		
Urinals				
COMUNS1	Urinal Male - Non spill	Autoclavable for sanitising Integrated handle Patented non-spill spout prevents spillage in any position Spout unscrews easily for disposal and cleaning 100ml graduation markings Moulded plastic	Weight 240g 1 litre capacity 500ml non-spill capacity	

COMUST1	Urinal Male - Standard	Standard urinal bottle for males. Ideal for bed confinement or travelling	Graduated up to 1000 ml Secured lid to reduce spills and odours	
	COFUST	Urinal Female	Assists the user or carer in personal toileting Ideal for bed confinement or travelling Anatomically shaped for females Secured lid to reduce spills and odours	
	<hr/>			
	COASTI	Absorbent Sheet - With Bed Tuck Ins	Highly absorbent draw sheet with side tuck-ins. This item comes with built in waterproofing.	Width (not incl. tuck-ins): 1000mm Length: 1000mm Absorption Capacity: 2.5L
COWPTI	McIntoshes			
	Waterproof Sheeting / McIntosh PVC	Mack-E Mattress Protector with wings		

ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
DAILY LIVING AIDS				
Crockery				
DLACU05	Feeding Cup with Long Spout and Handles	Homcraft Feeding Cup, Standard 8mm Spout, 250ml, 2/pack		
DLAPL01	Plate Incurve Surround Homcraft	This flexible guard can be fitted to ordinary plates to assist with one-handed eating. Three clips attach it to the rim and create a high inward sloping face to help with food collection and prevent spillage.	Plate size 190 to 254mm external diameter. Dishwasher safe. Melamine material dishwasher safe to 180 degrees F (82°C).	
DLAPL03	Plate Manoy Large	Manoy Contoured Plate, Large, 279 x 197mm, White		
Cutlery				
DLACT12	Knife - Angled	The knife has a straight built-up handle with angled blade that allows cutting with a rocker action	Stainless steel utensils that are easy to clean and dishwasher safe. Length of handles: 114mm	

DLACT01

Fork Suregrips

Sure Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

Sure Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.



DLACT02

Knife Suregrips

SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. KNIFE IS NOT BENDABLE

Sure Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.



DLACT03

Tablespoon Suregrips

SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.



DLACT18

Teaspoon Goodgrips

SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.



DLACT19

GOODGRIPS Fork - Weighted and Bendable

Good Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

The Weighted Utensils range are available with 170g of extra weight in the built up handle to provide more control for people with a tremor or limited hand control, Parkinson's disease or spasticity.

The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use.

35mm handle diameter
170g weighted handle

**DLACT20**

GOODGRIPS Tablespoon - Weighted and Bendable

Good Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

The Weighted Utensils range are available with 170g of extra weight in the built up handle to provide more control for people with a tremor or limited hand control, Parkinson's disease or spasticity.

The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use.

35mm handle diameter
170g weighted handle



DLACT21

GOODGRIPS Teaspoon -
Weighted and Bendable

Good Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.

The Weighted Utensils range are available with 170g of extra weight in the built up handle to provide more control for people with a tremor or limited hand control, Parkinson's disease or spasticity.

The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use.

35mm handle diameter
170g weighted handle

**DLACT22**

GOODGRIPS Knife - Weighted

Good Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.






The Weighted Utensils range are available with 170g of extra weight in the built up handle to provide more control for people with a tremor or limited hand control, Parkinson's disease or spasticity.

The knife is weighted but not bendable.

35mm handle diameter
170g weighted handle



Dressing Aids

DLADA01	Button Hook Modular Kings Complete	Assists one-handed dressing Stainless steel wire passes through button hole, gripping the button. It is pulled through with a slight twisting motion Contoured handle	Length 230mm Weight 70g	
DLADA02	Dressing Stick	Wooden stick with epoxy coated hook at one end and rubber tip at the other. Assists in applying and removing clothing. Can also be used to push or pull objects.	66cm length 5/8"-diameter wooden dowel. Finished with lacquer with reinforced hook.	
DLADA03	Shoe Horn Long Handled Plastic	Plastic handle designed to assist in putting on shoes especially for people who have difficulty bending and reaching.	Length 60cm. Blue Colour. Includes wrist strap.	
DLADA04	Stocking Aid	Stocking aid allows the user to don socks and stockings independently. Side notches hold the sock securely, while the strap ends help pull stocking into position efficiently. Built with a flexible plastic gutter for improved maneuverability.	Available in contoured or notched. Straps 36". Gutter length 9".	
DLADA05	Sock Aid - Single Plastic	Simple aid designed to assist with putting on a stocking especially if bending and reaching the foot is difficult. Cut out side notches designed to hold stocking 1 metre length pull straps		

DLADA07

Stocking Compression
Applicator Ezy-As Small

The unique design of the Ezy-As
compression garment applicator makes
it possible to assist in the fitting of a
wide range of compression garments.

Small 28 - 35 cm



DLADA08

Stocking Compression
Applicator Ezy-As Medium

The unique design of the Ezy-As
compression garment applicator makes
it possible to assist in the fitting of a
wide range of compression garments.

Medium 32 - 40cm



DLADA09

Stocking Compression
Applicator Ezy-As Large

The unique design of the Ezy-As
compression garment applicator makes
it possible to assist in the fitting of a
wide range of compression garments.

Large 38cm - 47cm



DLADA10

Compression Stocking /
Garment Applicator - Small

Small



DLADA11

Sock Aid - Terry Towelling

Simple aid designed to assist with putting on a sock especially if bending and reaching the foot is difficult. Terry towelling outer fabric for greater hold with smooth nylon inner surface to reduce friction.

Nylon pull straps - 700mm



Household Aids

DLAHHA3

Tap Turner - Crystal

This Crystal Tap Turner will assist you to use the tap on any tap that does not have an 'X' style knob. The Turner clamps onto the tap and then you can use the handle to control the tap. This is a fantastic tool for those who struggle to grasp a tap and close their hand. Instead, this tool gives you a larger area handle to grip on to.

Length 190mm. Weight 112g



DLAHHA4

Tap Turner - Red/Blue

The Homecraft Tapturn Turners are strong moulded turners which simply hook over the top of a 'X' type tap, to make it easy to turn. They are supplied as a red and blue pair for hot and cold taps.

Length 135mm
Width 47mm



DLAHHA8

Non Slip Matting

A roll of versatile soft material that can be used on any surface to help stop items from moving, or to protect against scratches.

This non-slip matting has a slightly spongy texture which prevents items moving or slipping on almost any surface.

Comes in a roll (size 1000 mm long by 300 mm wide) and can be cut to size with household scissors.



Tray Mobile

DLATM01

Trolley Walker - Adjustable Height

Includes clip-on plastic trays that can be removed easily for cleaning. Can be used with or without trays or it can be used as a table.

Rails on tray edges prevent items from slipping off.

Lockable hand brakes mean this trolley can also be used as a walking aid and provides additional safety when standing from a seated position.

Handles are adjustable in height.

Easily accessible loop lockable brakes.

Supplied with a cream painted steel frame.

Comes in white or beige colour.

Overall Width: 585mm

Overall Depth: 560mm




Handle Height Range:

755mm - 905mm

Wheel Diameter: 150mm (6")

Maximum User Weight: 150kg



DLATM02	Traymobile (no brakes)	<p>Ideal for indoor use</p> <p>Two removable trays to assist transporting smaller household items</p> <p>Height adjustable</p> <p>Four swivel castors provide greater manoeuvrability</p>	<p>Frame Height 800-950mm</p> <p>Frame Width 490mm</p> <p>Frame Depth 540mm</p> <p>Max User Weight 125kg</p>	
<div data-bbox="398 427 586 456" data-label="Section-Header"> <p>Kitchen Aids</p> </div> <div data-bbox="85 475 275 735">DLAKA13</div>				
	Plastic Spreading Board	no info available		
DLAKA12	Bread Board	<p>Ideal for people with tremors or impaired sight, the clever design of this cutting aid makes kitchen work much easier. Food is firmly held in the device while a series of rubber stopper feet keep the board from moving around on your kitchen bench or table. Great for bread, cheese, vegetables or other food items, your knife slides between the two guides, making it easy to cut straight and safe, with no risk of injury to fingers or hands. The board can be adjusted for a range of different cutting thickness</p>		

DLAKA14

ETAC Food Prep System

This durable cutting board is ideal for those who have difficulty in gripping, or only have the use of one hand.

The clamp can fix items such as tins or bowls, and the stainless steel spikes keep vegetables still when peeling or cutting.

It has four non-slip rubber feet to prevent the board from moving.

Suction pads are also provided for use on wet surfaces.

Size: 315mm × 295mm.
Weight: 1.66kg.



DLAKA15

Cordless Kettle Tipper

Suits most cordless kettles.

There is a metal plate that secures the base of the cordless kettle and a frame with a Velcro strap to keep the upper part of the kettle in place.

A cut out at the foot of the tipper allows cups to be positioned as close to the kettle as required.

The frame has a built-in back stop to prevent spillage if the kettle is accidentally released during use.

Weight: 1kg



DLAKA05

Jar & Bottle Opener - Mighty Lever

The Jar and Bottle Opener has a flexible stainless steel band which is fitted around a lid, and then clamped tight using the screw at the end of the handle. The handle can then be used as a lever.

Suitable for tops 20-110mm diameter. Weight 88g.



DLAKA06

Jar Twister

A cone shaped rubber moulding with fluted finger grips on the outside, and ribbing on the inside to give a good grip for opening twist-off jar lids.

Suitable for tops 20 to 85mm diameter.
Weight: 100g



DLAKA10

Zyliss 5 in 1 Can Opener

This versatile opener is the ideal gadget for opening pull tabs on cans, plastic and metal twist-off caps, metal bottle caps, pop-tops and jar lids. No more struggling to get the caps and lids off bottles and jars, this handy gadget has all your needs covered. The soft touch grip makes it comfortable and simple to use.

Dimensions:
12Lx15Wx4H



Personal Care Aids

DLAPA05

Long Handled Toe Wiper

Helps people who have limited access to their feet toes to wash and dry between their toes.

Long bendable wire stem with a plastic built-up easy-grip handle. Flat wire frame covered in a white towelling pad for washing between the toes. Includes 2 pads. Length: 710mm. Weight: 120g.



DLAPA09

Reaching Aid Economy - 67cm

Trigger style mechanism pick up stick 67cm
Slip resistant jaw provides a secure hold
Hook trigger
Magnetic tip



DLAPA10

Reaching Aid Economy - 82cm

Trigger style mechanism pick up stick 82cm
Slip resistant jaw provides a secure hold
Hook trigger
Magnetic tip



DLAPA11

Reaching Aid Squeeze - short

Slip resistant rotating jaw designed to pick up objects at different angles 10204/65 (65 cm)
Simple hand trigger operation
Magnetic tip
Wheelchair/walking frame clip supplied.




DLAPA12

Reaching Aid Squeeze - long

Slip resistant rotating jaw designed to pick up objects at different angles
Simple hand trigger operation
Magnetic tip
Wheelchair/walking frame clip supplied.

10204/81 (81 cm)



DLAPA13

Bottom Wiper

The Homecraft Toilet Aid provides assistance with independent hygiene. The aid is designed with a curved metal handle and a molded-plastic head that gives the user the proper reach. The handle's recessed serrations allow it to securely grip paper.

Length 25.4cm. Weight 71g.




DLAPA14**Long Handled Sponge**

As you get older, it can be difficult to bend down to wash your feet or to lift them up so you can wash them. Doing so might even put you at risk for a slip and fall in the shower. This long-handled bath sponge has a wire that is strong, but also easily bendable to be rotated, twisted, and moved in any position that helps you accomplish optimal cleaning without bending! It can be used by seniors and people with disabilities to bathe independently, by caretakers to easily clean residents, or for household cleaning tasks.

Length 24". Sponge head 3.75" x 5".



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
ELECTRIC BEDS				
Electric Adjustable Beds				
BREBA01	Community Electric Bed	<p>A versatile and adaptable nursing bed which is easy to use and built for strength and durability for the home care market. Electric bed, with electric backrest/headraise, knee break and manual leg lift. Trendelenburg and Reverse Trendelenburg. 4 locking castors. Bed folds for transportation.</p> <ul style="list-style-type: none"> • Easy to fold and superior clearance when folded for transport or storage • Multiple positioning to suit a range of care needs • Easy to use hand controller with the ability to lockout individual functions • Linak actuators providing superior reliability • Strong and durable frame with a 5-year warranty 	<ul style="list-style-type: none"> • Safe Working Load - 220kg • Dimensions - 2324x1050mm • Height Range - 500 to 730mm • Mattress size - 1960 • Bed Frame Weight - 74 kg • Warranty - 5 year bed frame / 2 year electrical / 1 year castor 	

BREBA02

Electric Bed, Aged Care

A fully-featured, easy to operate, electric bed.

Four-part, high low profiling bed

Adjustable bed height, backrest, knee break and heel lift

Quick-release head and foot panels

Choice of panel designs

Central locking castors

In-built 200mm base extension

Electric hand controller with auto contour and a lockout function

4 Detachable mattress retainers

Business box holder

- Safe working load - 220 kg

- Maximum patient weight - 185 kg

- Warranty - 3 year

Electrical / 10 year

Frame / 1 year Back-Up

Battery

- Frame material - Powder Coated Steel

- Model - Standard, King

Overall dimensions



BREBICP

iCare King Single Complete Bed Package

Safe working load 200kg on king single.

Hand Controller Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.

Hard wearing commercial grade fabric.

Onyx range (BLACK) 82% Polyester, 18% Olefin. Stone range (GREY) 70% Polyester, 30% Acrylic.

Strong metal mesh under mattress allowing increased airflow

8x heavy duty swivel castors with one lockable castor on each corner.

Head lift angle 62° and knee break angle 42°

Max height 66cm from floor to top of base

Min height 22cm from floor to top of base.

Includes the IC20 Mattress a mid range mattress with a soft feel but excellent support and a matching headboard in STONE fabric.

Safe working load - 200kg.

Minimum height from floor to base is 220mm.



Accessories

BREBBS1


Bedstick Low - Townsend for Electric Beds

Clamped firmly to bed base railing
Acts as a support to roll over, to sit from lying position and to assist with getting in and out of bed.

Low profile 525 x 175
2.1kg
Max User Weight 136kg

For use for Hospital Style Beds

Manufactured in Adelaide,
Powdercoated mild steel. Clean using mild detergent and water and dry with soft cloth.

**BREBBS2**


BedstickStandard - Townsend for Electric Beds

Clamped firmly to bed base railing
Acts as a support to roll over, to sit from lying position and to assist with getting in and out of bed.

Standard Height 715x175
2.4kg
Max User Weight 136kg


For use for Hospital Style Beds

Manufactured in Adelaide,
Powdercoated mild steel. Clean using mild detergent and water and dry with soft cloth.

**BREBBR3**

Bedrail for Elec Bed - Fold Down - 3/4 length

Clamp on, fold down metal side rail system
attaches firmly to bed frame. Three quarter length rails for increased access.



BREBBR1

Bed Rail - Full Length for Electric Bed

Drop down metal side rail system attaches firmly to bed frame. Full length rails for added safety and security.

**BREBBR2**

Bed Rail Cover - Padded for Electric Bed







Soft cushions to encase metal bed rails Standard 670 x 1260mm protects limbs from hitting rails or getting caught. Soft, medical grade foam cushions, stretch cover, easy secure Velcro closure, easy wipe clean. Two sizes available, fits majority of bed rails.



**BREBSHP**

Pole Self Help For Community Bed

Support pole to aid patient repositioning while lying in bed. Mounts into inbuilt brackets at either side of bed head.



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
HOME ACCESS				
Portable Ramps				
HARMP1	Multipurpose Ramp 870mm	Fibreglass ramp that has a concertina folding action for storage and transport. 875 mm length.	835mm width 875mm length 5.6kg	
HARMP2	DecPac Portable Ramp 700mm	Fibreglass ramp that has a concertina folding action for storage and transport. 700mm length.	745mm width 700mm length 4.7kg	
HARMP3	DecPac Portable Ramp 1200mm	Fibreglass ramp that has a concertina folding action for storage and transport. 1200mm length.	750mm width 1200mm length 8.0kg	
HARMP3E	DecPac Portable Ramp 1200mm (with EDGE)	Fibreglass ramp that has a concertina folding action for storage and transport. 1200mm length with EDGE.	820mm width 1200mm length 9.7kg	
HARMP5	DecPac Portable Ramp 1650mm	Fibreglass ramp that has a concertina folding action for storage and transport. 1650mm length.	750mm width 1650mm length 10.8kg	
HARMP5E	DecPac Portable Ramp 1650mm (with EDGE)	Fibreglass ramp that has a concertina folding action for storage and transport. 1650mm length with EDGE.	820mm width 1650mm length 13kg	

HARMP6	DecPac Portable Ramp 2000mm	Fibreglass ramp that has a concertina folding action for storage and transport. 2000mm length.	750mm width 2000mm length 11.7kg	
	DecPac Portable Ramp 2000mm (with EDGE)	Fibreglass ramp that has a concertina folding action for storage and transport. 2000mm length with EDGE.	820mm width 2000mm length 13.7kg	

ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
LIFTERS AND HOISTS				
Lifters				
LHEH01	Patient Lifter - 180kg	<p>The LiftAid 180 is a lightweight, easy to operate mobile lifter that allows patients to be lifted and transferred from floor, bed or chair with ease and safety.</p> <p>3 height settings - tool-less adjustment</p> <p>Quick release hook system</p> <p>Electrically operated base</p> <p>Handheld controller</p> <p>Emergency stop</p> <p>Manual emergency lowering</p> <p>Optional Pivot Frame</p> <p>Compatible with range of slings</p>	<p>Width: 61.5cm</p> <p>Height: 43cm - 181cm (3 height settings)</p> <p>Length: 109cm</p> <p>Product Weight: 35.3kg</p> <p>Warranty: Lifter & Actuator - 2 years Battery - 1 year</p> <p>S.W.L: 180kg</p> <p>Max Spread: 101cm</p> <p>Underbed Clearance: 12.5cm</p>	
LHEH02	Patient Lifter - 320kg	<p>The LiftAid 320 is a heavy-duty mobile lifter ideal for the majority of common lifting situations.</p> <p>Patients can be lifted and transferred from floor, bed or chair. Fold-down armrests assist with gait training during rehabilitation or for general toileting.</p> <p>Adjustable height frame allows overall mast height be set to match environmental and lifting needs. Features the advanced Linak® Jumbo Care Diagnostic Control System for intelligent self-diagnostic system capabilities. Low lifting height and soft start/stop function,</p> <p>3 height settings with tool-less adjustment</p> <p>quick release hook system for easy sling</p>	<p>Length Overall - 1410 mm</p> <p>Width Overall 7950 mm - 1190 mm</p> <p>Height Overall 2180 - mm</p> <p>Height of Lift 650 mm - 1840 mm</p> <p>Clearance Under Bed 16 - mm</p> <p>Diameter of Front Wheel - 105 mm</p> <p>Diameter of Rear Wheel - 127 mm</p> <p>Load Capacity - 320 kg</p>	

LHSLGPS

Lifter/Hoist Slings

Sling General Purpose - Sml

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports
Suits the majority of lifters on the market
Brightly coloured handles and adjusters secured to sling for easy access and visibility
Quick-drying mesh fabric machine washable up to 85°C for microbiological control
Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

840 x 1160 x 490 mm (W x L x D)
Capacity - 205 kg



LHSLGPM

Sling General Purpose - Med

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports

Suits the majority of lifters on the market

Brightly coloured handles and adjusters secured to sling for easy access and visibility

Quick-drying mesh fabric machine washable up to 85°C for microbiological control

Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

890 x 1200 x 545 mm (W x L x D) Capacity - 205kg



LHSLGPL

Sling General Purpose - Lge

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair.

Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports Suits the majority of lifters on the market Brightly coloured handles and adjusters secured to sling for easy access and visibility Quick-drying mesh fabric machine washable up to 85°C for microbiological control Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

1135 x 1325 x 640 mm
(W x L x D)
Capacity - 205kg



LHSLGPXL

Sling General Purpose - Xlge

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports
Suits the majority of lifters on the market
Brightly coloured handles and adjusters secured to sling for easy access and visibility
Quick-drying mesh fabric machine washable up to 85°C for microbiological control
Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort



LHSLGPHS



Sling General Purpose Head Support Sml

General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).

For Yoke frame lifters only

900, 230, 450, 790, 965 mm, 205 kg



LHSLGPHM	Sling General Purpose Head Support Med	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	1010, 230, 490, 890, 1060 mm, 205 kg	
LHSLGPHL	Sling General Purpose Head Support Lge	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	1080, 250, 480, 990, 1180 mm, 300 kg	

LHSLGPHXL

Sling General Purpose Head Support XLge

General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).

For Yoke frame lifters only



LHSLHS

Sling Hygiene Sml

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access
buckle waist strap for added security and increased infection control.
Generous padding under arms
3 loop support strap system controls positioning of head and leg supports.
Suits majority of lifters on the market.
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
manufactured from the highest quality material.



LHSLHM

Sling Hygiene Med

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access
buckle waist strap for added security and increased infection control.
Generous padding under arms
3 loop support strap system controls positioning of head and leg supports.
Suits majority of lifters on the market.
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
manufactured from the highest quality material.

Maximum user weight
205kg



LHSLHL

Sling Hygiene Lge

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access
buckle waist strap for added security and increased infection control.
Generous padding under arms
3 loop support strap system controls positioning of head and leg supports.
Suits majority of lifters on the market.
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
manufactured from the highest quality material.

Maximum user weight
300kg



LHSLHXL

Sling Hygiene XLge

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access
buckle waist strap for added security and increased infection control.
Generous padding under arms
3 loop support strap system controls positioning of head and leg supports.
Suits majority of lifters on the market.
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
manufactured from the highest quality material.

Maximum user weight
300kg



LHSLHHS

Sling Hygiene With Head Support Sml

Open area from mid back to mid thigh for hygiene access
inbuilt head support for patients with less control
buckle waist strap for added security and increased infection control
generous padding under arms
3 loop support strap system controls positioning of head and leg supports
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
suits majority of lifters on the market
manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)

Maximum user weight 205kg



LHSLHHM

Sling Hygiene With Head
Support Med

Open area from mid back to mid thigh
for hygiene access
inbuilt head support for patients with
less control
buckle waist strap for added security
and increased infection control
generous padding under arms
3 loop support strap system controls
positioning of head and leg supports
brightly coloured handles and
adjusters, secured to sling for easy
access and visibility
quick drying mesh fabric; machine
washable up to 85°C for microbiological
control
suits majority of lifters on the market
manufactured from the highest quality
material, in compliance with industry
standards (AS/NZS ISO 10535:2011)

Maximum user weight
205kg



LHSLHHL

Sling Hygiene With Head
Support Lge

Open area from mid back to mid thigh
for hygiene access
inbuilt head support for patients with
less control
buckle waist strap for added security
and increased infection control
generous padding under arms
3 loop support strap system controls
positioning of head and leg supports
brightly coloured handles and
adjusters, secured to sling for easy
access and visibility
quick drying mesh fabric; machine
washable up to 85°C for microbiological
control
suits majority of lifters on the market
manufactured from the highest quality
material, in compliance with industry
standards (AS/NZS ISO 10535:2011)

Maximum user weight
300kg



LHSLHHXL

Sling Hygiene With Head Support XLge

Open area from mid back to mid thigh for hygiene access
inbuilt head support for patients with less control
buckle waist strap for added security and increased infection control
generous padding under arms
3 loop support strap system controls positioning of head and leg supports
brightly coloured handles and adjusters, secured to sling for easy access and visibility
quick drying mesh fabric; machine washable up to 85°C for microbiological control
suits majority of lifters on the market
manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)

Maximum user weight 300kg



Stand Aid Slings




LHSLSATM

Sling - Standard Torso Medium

Fully padded, including soft padding under the arms
Various attachments
Built-in Keyhole Plate Safety Straps (Applicable for the version of this product which uses Keyhole Plates)
A comfortable, well padded, high back

Chest Size 85cm - 110cm
Load capacity 200kg



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
MOBILITY AIDS				
Rollator Frames				
MARFMIN	Rollator Walking Frame - Mini	Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Max height 82cm Min height 72cm. Max user weight. 150kg	
MARFST	Rollator Walking Frame - Standard 5 inch wheels	Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Adjustable height 810 - 910mm. Overall width 670mm. Width between handles 430mm Depth 500mm. Recommended mainly for indoor use Max user weight. 150kg	
MARFBA	Rollator Walking Frame - Bariatric 227kgs	Oversize frame, wider and deeper to accommodate individuals up to 227kgs Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Maximum user weight: 227kg Height adjustable: 32.5"~39.5" with 5 wheels	

Pick Up Walking Frames

MAPWFMN

Pick Up Walking Frame - Mini

Height adjustable
Folds flat for easy storage and transportation
Foam handles
Lightweight aluminium
Tip size 28mm

Handle Height 700 - 805mm
Max Base Size 630 x 500mm
Weight 2.5kg
Max user weight 150kg



MAPWFST

Pick Up Walking Frame -- Standard

Folding aluminium frame with foam rubber hand grips.

Adjustable height 770 - 940mm with rubber ferrules on all 4 legs.
Width 455mm, Depth 400mm.



MAPWFBA

Pick-up Walking Frame - Bariatric 227kgs

Oversize frame design provides wider and deeper to accommodate individuals up to 227kgs
Precision design provides additional strength while adding minimal weight

Maximum user weight 227kgs
Height adjustable: 32.5"~39.5"



Walking Sticks

MAWSTH**Walking Sticks - T Handle**

T shaped handle, single point, adjustable. To provide support people who have difficulty walking due to balance problems, pain or injury.

The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.

Handle height: 79 cm - 102 cm (25 mm increments) Approximate user height: 5' 3" - 6' 1" Rubber tip/ferrule: 19 mm Product weight: 300 gm Maximum user weight capacity: 114 kg




**MAWSB01****Walking Stick - Bariatric - 225 Kg**



This sturdy steel walking stick provides support to a person who has difficulty walking due to balance problems, pain or injury.

The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.

Rubber Hand Grip
Handle Height: 73 cm - 96 cm (25 mm Increments)
Approximate User Height: 5' 2" - 6' 1"
Rubber tip/ferrule: 22 mm
Maximum user weight capacity: 225 kg



MAWSCH	Walking Stick - Crook Handle	To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight. The height of the stick can be adjusted to	Height adjustment range: 760mm - 990mm	
MAWSEL	Walking Stick - Ergonomic - Left Hand	Anatomical Walking StickConstructed from aluminium to allow durability and strength whilst being light in weight. Supplied with anti-slip rubber tip. This walking stick is available in both a left hand and right hand option to ergonomically fit the natural shape of the user's hand for maximum comfort and ease of use.	Height adjustment range: 760mm - 990mm	
MAWSER	Walking Stick - Ergonomic - Right Hand	Anatomical Walking StickConstructed from aluminium to allow durability and strength whilst being light in weight. Supplied with anti-slip rubber tip. This walking stick is available in both a left hand and right hand option to ergonomically fit the natural shape of the user's hand for maximum comfort and ease of use.	:Height adjustment range: 760mm - 990mm	

MAWSF	Walking Stick - Folding	<p>A lightweight, aluminium walking stick with compact folding mechanism for convenient storage. The stick folds and unfolds simply and easily, making it ideal for everyday use or travel.</p> <p>The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>lightweight folding walking stick height adjustable easy to collapse and unfold maximum user weight 113kg black colour</p>	
MAWSSN	Walking Stick - (Swan neck) with strap	<p>To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight.</p> <p>The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Height adjustment range: 810mm - 1035mm Tube size: Upper tube 22.2mm, lower tube (telescopic) 19mm Rubber tip (Outer size Ø38XH48; Inner sizeØ18.5xD34) Maximum User Weight: 113kg Colour: Silver</p>	

MAWSSNS

Walking Stick - (Swan neck)
without strap

To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight. The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.

Height adjustment range:
810mm - 1035mm
Tube size: Upper tube
22.2mm, lower tube
(telescopic) 19mm
Rubber tip (Outer size
Ø38XH48; Inner
sizeØ18.5xD34)
Maximum User Weight:
113kg
Colour: Silver



4 Wheel Walkers**MA4WMI**

4 Wheel Walker - CAREQUIP
MINI

Frame and backrest fold compact for storage and transportation
Height adjustable handles
Narrow frame and low seat height ideal for the smaller client, or where space is limited
Handbrakes can be locked for added stability
Supplied with removable backrest, vinyl bag and padded seat
Walker can be folded with vinyl bag still attached
Curved padded backrest for additional comfort and support
*Note: Product comes with a vinyl bag and not a basket as pictured.

Seat height 470mm Seat
width 310mm Max
user weight 130kg
Unit weight 6.2kg



MA4WSB

4 Wheel Walker - STANDARD
BLUE

Frame and backrest fold compact for
storage and transportation
Height adjustable handles
Soft-lock handbrakes
Curved padded backrest for additional
comfort and support
Supplied with vinyl bag and padded seat
Walker can be folded with vinyl bag still
attached

Seat height 525mm Seat
width 350mm Max user
weight 130kg
Unit weight 5.9kg



MA4WSR

4 Wheel Walker - STANDARD
RED

Frame and backrest fold compact for
storage and transportation
Height adjustable handles
Soft-lock handbrakes
Curved padded backrest for additional
comfort and support
Supplied with vinyl bag and padded seat
Walker can be folded with vinyl bag still
attached

Seat height 525mm Seat
width 355mm Max user
weight 130kg
Unit weight 5.9kg



MA4WRE

4 Wheel Walker - DAYS RED

The Days 100 Series aluminium frame ensures a combination of lightweight design and robust build quality with a maximum user weight of 125kg.

6" wheels

Ergonomic handgrips to make them easy to grip.

Adjustable height handles with 5 settings at 25mm intervals

Rest Seat

Storage Bag

Loop lockable brakes provide added security and convenience. The brake tension can also be easily adjusted to suit requirements

Foldable backrest

Folds easily for storage & transport - simply fold together by pulling the strap on the seat

Max User Weight: 125kg

Seat Height: 500mm

Seat Width: 350mm

Seat Depth: 320mm

Width Between Arms:

460mm

Backrest Height: 750mm

Overall Height: 755-875mm

Overall Width: 620mm

Overall Depth: 640mm

Backrest: Straight, folding

Wheel: 6"

Tyre: PVC

Net Weight: 6kg



MA4WBSR

DAYS Walker Boxed -
ADJUSTABLE Standard Red

Easy to assemble - no tools required
Retail carrying box with handle
Height adjustable seat and handles
Padded seat with storage compartment
bag underneath
Padded backrest
Hand operated brake system

Max User Weight: 136kg
with foldable backrest
Seat Height: 500 -
600mm
Seat Width: 356mm
Seat Depth: 305mm
Overall Height: 785 -
1000mm
Overall Width: 622mm
Overall Depth: 680 -
800mm
Width Between Handles:
460mm
Backrest: Curved, folding
Wheel: 7.5"
Tyre: PVC, Black
Net Weight: 7.4kg



MA4WBSB

DAYS Walker Boxed-
ADJUSTABLE Standard Blue

Easy to assemble - no tools required
Retail carrying box with handle
Height adjustable seat and handles
Padded seat with storage compartment
bag underneath
Padded backrest
Hand operated brake system

Max User Weight: 136kg
with foldable backrest
Seat Height: 500 -
600mm
Seat Width: 356mm
Seat Depth: 305mm
Overall Height: 785 -
1000mm
Overall Width: 622mm
Overall Depth: 680 -
800mm
Width Between Handles:
460mm
Backrest: Curved, folding
Wheel: 7.5"
Tyre: PVC, Black
Net Weight: 7.4kg

**MA4WLB**

4 Wheel Walker - BLUE (LARGE
WHEELS)

8" wheels
Ergonomic handgrips and adjustable
height handles
Rest Seat
Storage Bag
Loop lockable brakes
Foldable backrest
Folds easily for storage & transport

Max User Weight: 125kg
Seat Height: 565mm
Seat Width: 350mm
Seat Depth: 320mm
Width Between Arms:
460mm
Backrest Height: 800mm
Overall Height: 840-
940mm
Overall Width: 620mm
Overall Depth: 645mm
Backrest: Straight,
folding
Wheel: 8"
Tyre: PVC
Net Weight: 7.2kg



MA4WLG**4 Wheel Walker - GOLD (LARGE WHEELS)**

8" wheels
Ergonomic handgrips and adjustable height handles
Rest Seat
Storage Bag
Loop lockable brakes
Foldable backrest
Folds easily for storage & transport

Max User Weight: 125kg
Seat Height: 565mm
Seat Width: 350mm
Seat Depth: 320mm
Width Between Arms: 460mm
Backrest Height: 800mm
Overall Height: 840-940mm
Overall Width: 620mm
Overall Depth: 645mm
Backrest: Straight, folding
Wheel: 8"
Tyre: PVC
Net Weight: 7.2kg

**MA4WHD****4 Wheel Walker - HEAVY DUTY**

To provide support and stability whilst walking. Seat lifts up to allow for storage in basket below
7.5" wheels
Adjustable height handles
Rest seat
Storage bag
Loop lockable brakes
Folds easily for storage & transport
Screen reader support enabled.

Max User Weight: 200kg
Seat Height: 570mm
Seat Width: 495mm
Seat Depth: 355mm
Overall Height: 855 - 975mm
Overall Width: 735mm
Overall Depth: 655mm
Backrest: Curved
Wheel: 7.5"
Tyre: PVC, Black
Net Weight: 8.8kg
Screen reader support enabled.



MA4WPD

4 Wheel Walker -PUSH DOWN
BRAKES **SPECIAL**
ORDER

To provide support and stability whilst walking. Seat lifts up to allow for storage in basket below. The unique brake system is activated by pushing down on the handles or whenever you sit on the seat. The robust yet lightweight alloy frame has tube-in-tube construction for increased strength and longer lifespan. The castors are made from solid rubber, making them free from punctures. They are designed for easy operation and good manoeuvrability especially on indoor surfaces. Their slightly smaller size makes them especially good for negotiating doorways or tight spaces. The orthopedic moulded handgrips are designed to contour entirely to your hand, minimising stress or strain to the wrist. The height of the handles is easy to adjust, with an extension range of 13 cm to ensure that the handles can be set to best match the height of the user. With a compact fold down design, this walker is easy to transport or store when not in use.

Suitable for people up to
130 kilograms in weight.
Height 945mm max -
810mm min
Length 610mm
Width 575mm
Seat Depth 310mm
Seat Width 360mm
Seat height 557mm
Folded width 215mm
Castor size 153mm



MA4WGST

Gutter Frame

Sturdy wheeled walking frame for indoor mobility training and rehabilitation. Designed for users who require the ability to weight bear through the forearms. Available in both non folding and folding models.

Note: Not designed to help users stand from a seated position

Forearm cuff height: 93 cm - 123 cm (25 mm increments)
Width: Front 40 cm / Rear 60 cm
Front swivel wheels and rear ski's for maneuverability
Adjustable hand grips to suit users forearm length
Maximum user weight capacity: 110 kg



Tri Walkers

MA3WW

Tri Walker with Tray and Basket

The Aluminium Tri Wheel Walkers have a lightweight aluminium frame that makes these walkers extremely easy to manoeuvre making them ideal for use indoors or out.

Ergonomic handgrips for optimum comfort
Adjustable height handles with 5 settings at 25mm intervals
Loop lockable brakes provide added security and convenience. The brake tension can also be easily adjusted to suit requirements
Vinyl bag supplied as standard
Folds easily for storage or transport using the simple easy to grab pull handle
Supplied with basket and tray.



Crutches





MACRL01

Crutches Axillary Large

Aluminium, push button adjustable height axilla crutches with rubber ferrules.

Maximum User Weight -
136kg
135 - 155cm



MACRM01	Crutches Axillary Medium	Aluminium, push button adjustable height axilla crutches with rubber ferrules.	Maximum User Weight - 136kg 114 - 135cm	
	MACRS01	Crutches Axillary Youth/Small	Durable and strong, these Aluminium Underarm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Fitted with high quality rubber tips for stability.	Maximum User Weight - 136kg 94 - 114cm 
	MACREM1	Elbow Crutches Medium	Durable and strong, these Forearm Crutches are constructed from high tensile extruded aluminium, making them both lightweight and very strong. Fitted with high quality rubber tips for stability.	Maximum User Weight - 160kg 89 - 120cm 
	MACRERM	Elbow Crutches Medium Ergonomic Grip	Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Ergonomic palm grip offers pressure distribution for maximum comfort. Fitted with high quality rubber tips for stability. Both hand to floor and palm grip to arm cuff sections are adjustable.	Maximum User Weight - 160kg Hand Grip Height - 66 - 94cm, Arm Ring Adjustment- 20.5 - 27.5cm 

MACRERL	Elbow Crutches Large Ergonomic Grip	<p>Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Ergonomic palm grip offers pressure distribution for maximum comfort. Fitted with high quality rubber tips for stability.</p> <p>Both hand to floor and palm grip to arm cuff sections are adjustable.</p>	<p>Maximum User Weight - 160kg Hand Grip Height - 83 - 111cm, Arm Ring Adjustment- 20.5 - 27.5cm</p>	
Quad Sticks				
MAQSS01	Quad Stick Small - 15cm x 20cm	<p>This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided,</p>	<p>Handle height: 74 cm - 97 cm (25 mm increments) Soft foam handgrip Base: 20 cm x 15 cm Rubber tips/ferrule: 13 mm Product weight: 800 gm Maximum user weight: 114 kg</p>	
MAQSL01	Quad Stick Large - 17cm x 28cm	<p>This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided,</p>	<p>Handle height: 76 cm - 98 cm (25 mm increments) Soft foam handgrip Base: 28 cm x 17 cm Rubber tips/ferrule: 16 mm Product weight: 800 gm Maximum user weight capacity: 113 kg</p>	




MAQSOSL

Quad Stick - Off Set - 23 x 39cm

This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided, Adjustable height, larger base.

Height 795-1020mm
Base width 230mm
Base Depth 390mm
Tip size 16mm
SWL 75kg



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
MOBILITY AID ACCESSORIES				
IV POLES				
MAIVMO	IV Stand Mobile	2 Hook design 560mm nylon base Chromed steel pole 50mm castors	Adjustable height 1240-2130mm	
ROLLATOR ACCESSORIES				
MAARSK	Rollator Ski Replacement (x2)	Pair of Rollator Skis to replace lost or damaged skis.	To fit 25mm Rollator frames	
OXYGEN CARRIER				
MAAC02	Oxygen Bottle Carrier For Wheelchair	Oxygen bottle carrier frame to attach to wheelchair. Contact store for availability or suitability to client's chair.		
STUMP SUPPORTS				


MABF01

Fracture Board

Designed to support an extended leg and is normally used in a wheelchair. The Board can be used for either the left or right leg. Made from a strong plastic that is easy to wipe clean, to help Infection Control. This plastic is virtually maintenance free and will not splinter. The webbing straps can be removed for washing or replacing if necessary.

Length 112cm
Width 40cm



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
PRESSURE CARE				
Pressure Cushions - Air				
PCROL01	Cushion ROHO Low Profile 8x9 (fits 14x16 W/C)	<p>The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.</p> <p>Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.</p>	Cell size 8 x 9 Dimensions 38 x 43 cm Profile Low Height 60mm Max User Weight Unlimited	

PCROL02

Cushion ROHO Low Profile 9x9

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 9 x 9
Dimensions 43 x 43 cm
Profile Low
Height 60mm
Max User Weight
Unlimited



PCROL03

Cushion ROHO Low Profile 10 cells x 9 cells (fits 18x16 W/C)

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 10 x 9
Dimensions 47 x 43cm
Profile Low
Height 60mm
Weight Unlimited



PCROL04

Cushion ROHO Low Profile 10
cells x 10 cells

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 10 x 10
Dimensions 47 x 47 cm
Profile Low
Height 60mm
Max User Weight
Unlimited



PCROL05

Cushion ROHO Low Profile 11 cells x 10 cells (fits 20x18 W/C)

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 11 x 10
Dimensions 52 x 47 cm
Profile Low
Height 60mm
Max User Weight
Unlimited



PCROH01

Cushion ROHO High Profile 10 cells x 10 cells

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 10 x 10
Dimensions 47 x 47 cm
Profile High
Height 100mm
Max User Weight Unlimited



PCROH02

Cushion ROHO High Profile 11 cells x 10 cells (fits 20x18 W/C)

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 11 x 10
Dimensions 52 x 47 cm
Profile High
Height 100mm
Max User Weight Unlimited
Screen reader support enabled.



PCROH03

Cushion ROHO High Profile 8x9

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 8 x 9
Dimensions 38 x 43 cm
Profile High
Height 100mm
Max User Weight Unlimited



PCROH04

Cushion ROHO High Profile 9x9

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 9 x 9
Dimensions 43 x 43 cm
Profile High Height
100mm Max User
Weight Unlimited



PCROQ01

Cushion ROHO Quadtro 8 cells
x 9 cells HP (fits 14x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 8 x 9
Dimensions 38 x 43 cm
Profile High
Height 100mm
Max User Weight
Unlimited



PCROQ02

Cushion ROHO HP Quadtro 9 cells x 9 cells (fits 16x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 9 x 9
Dimensions 43 x 43 cm
Profile High
Height 100mm
Max User Weight
Unlimited



PCROQ03

Cushion ROHO Quadtro 10 cells x 9 cells High Profile (fits 18x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 10 x 9
Dimensions 47 x 43 cm
Profile High
Height 100mm
Max User Weight
Unlimited



PCROQH1

Cushion ROHO HP Quadtro
10x10

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 10 x 10
Dimensions 47 x 47 cm
Profile High
Height 100mm
Max User Weight
Unlimited



PCROQH2

Cushion HP ROHO Quadtro
11x10 (fits 20x18 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 11 x 10
Dimensions 52 x 47 cm
Profile High
Height 100mm
Max User Weight
Unlimited



Combination Cushions

PCJC11

Cushion Jay 2 Standard
460mmx410mm (18"x16")

Contoured foam cushion with Jay Flow fluid pad 46 x 41 cm
18"x16"

Designed for the client with moderate to high postural and pressure care needs

Firm, anatomically correct contoured foam base creates a foundation for postural stability

Provides excellent lateral and forward/rear stability

Fluid Tripad in rear seat well maintains even pressure distribution and keeps fluid under areas at high risk of skin breakdown

Supplied with breathable ballistic stretch cover



PCROHY2

Cushion ROHO Hybrid Elite
460mmx460mm (fits 18x16
W/C)

Combines ROHOs Air Flotation
qualities and Jays base technology
Ideal for clients who require moderate
to high pressure care and positional
support needs

Contoured foam base improves
alignment of the pelvis and lower
extremities, enhancing sitting posture
and function

Air cells are interconnected and their
unique shape promote blood circulation
to facilitate pressure prevention,
healing and relief

Flexible and adjustable air cells adapt
to patient contours and reduce shear
forces while transferring or
repositioning

Durable neoprene construction is soft
to suit skin sensitivity, and is easily
cleaned

Jay foam base is made from fluid-
resistant closed cell foam for additional
stability, assisting in safer wheelchair
transfers

Two-way stretch cover with encased
comfort foam is machine washable

Dimensions 46 x 46cm

Height 95mm

Max User Weight 225kg



Foam Cushions

PCCU02

Cushion Basic Foam 405mm x
405mm Black Cover

Max user weight 150kg
405mm x 405mm
50mm thick



PCCU03

Cushion Basic Foam 460mm x
405mm Black Cover

Max user weight 150kg
460mm x 405mm
50mm thick



PCCU04

Cushion Basic Foam 500mm x
450mm Black Cover

Max user weight 150kg
500mm x 450mm
50mm thick



PCCU05

Basic Foam Cushion with Gel Insert. (450 x 400 x 50mm)

This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.

Max user weight: 110Kg
450x400x50mm




**PCCU06**

Basic Foam Cushion with Gel Insert. (450 x 450 x 50mm)

This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.

Max user weight: 110Kg
450x450x50mm



PCCU07	Basic Foam Cushion with Gel Insert. (500 x 450 x 50mm)	This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.	Max user weight: 110Kg 500x450x50mm	
PCWCSC	Contoured Wheelchair Seat Cushion (18" x 16")	An anatomically moulded wheelchair cushion that is used to relieve and prevent back pain, numbness and other discomfort that can be experienced after prolonged sitting in a slung seat wheelchair.	18" x 16" 460 x 410 mm	
Comfort Products				
PCSS01	Sheepskin Medical Grade	Genuine Australian lambswool overlay Provides additional comfort to existing mattress or chair Soft pile allows air circulation to keep skin dry Natural fibre reduces the risk of skin irritation Machine washable Dimensions are calculated on an average	640mm wide x 960mm Long	

PCCFT02	Sheepskin Slippers Small - Pair	<p>Low profile medical sheepskin slipper</p> <p>pressure-redistributing and microclimate controlling properties</p> <p>non-slip suede soles</p> <p>slim and lightweight</p> <p>washable to 80°C</p>	
PCCFT03	Sheepskin Slippers Medium - Pair	<p>Low profile medical sheepskin slipper</p> <p>pressure-redistributing and microclimate controlling properties</p> <p>non-slip suede soles</p> <p>slim and lightweight</p> <p>washable to 80°C</p>	
PCCFT04	Sheepskin Slippers Large - Pair	<p>Low profile medical sheepskin slipper</p> <p>pressure-redistributing and microclimate controlling properties</p> <p>non-slip suede soles</p> <p>slim and lightweight</p> <p>washable to 80°C</p>	
PCCFT06	Sheepskin Heel Protectors	<p>Genuine medical sheepskin</p> <p>Soft pile allows air circulation to keep skin dry</p> <p>Natural fibres reduces the risk of skin irritation</p> <p>Adjustable velcro securing straps</p> <p>Machine washable</p> <p>One size fits all</p> <p>Supplied as a pair.</p>	

PCHBHE	Heel Protector - Bed Heel Elevator	<p>Wedge-shaped contoured foam provides excellent support for the lower leg</p> <p>Rigid foam aids blood flow and reduces pressure</p>
Static Mattresses		
PCFM01	Mattress - Area Care Plus - Static Foam	<p>The AreaCare Plus is available in Standard and King Single size. The three layered, zoned design of this very high risk pressure surface is ideal for patients with poor mobility or uncertain skin sensitivity.</p> <p>Three-layer foam pressure surface offers maximum pressure reduction in a three-layer foam surface for higher-risk patients.</p> <p>A top layer of high density, temperature-sensitive, memory foam conforms to the body's shape, working in combination with the higher resistance support layer to disperse pressure evenly. The top layer also incorporates a special immersion, slow recovery foam foot zone for the ulcer prone heel area. The core layer and edges provide form, strength, and stability. Standard model suitable for users up to 220 kilograms; King Single up to 250 kilograms.</p>



PCPRM01

Mattress - Area Care Pro Static Foam

The AreaCare Pro multi-density foam pressure surface offers maximum pressure reduction in a three layer foam surface for higher risk patients. The three layered, zoned design of this very high risk pressure surface is ideal for heavier patients with poor mobility or uncertain skin sensitivity. Made from premium materials, each foam layer is made of varying density to maximise the redistribution of pressure. The top layer memory foam also incorporates a special immersion, slow recovery foam foot zone for the ulcer prone heel area, improves ventilation to the skins surface and increases the overall weight capacity to 220 kilograms.

Single (1960 x 865 x 150mm)



Alternating Pressure Reduction Mattresses



PCMR01


ProCair Mattress Replacement System

ProCair is designed for individuals at the high-risk pressure injury category. Built for optimal patient therapy and comfort, ProCair is now available in Standard and King Single sizes with sealed base to support your infection control strategy.

Weight Capacity (kg)
40-200



PCMR02	ProCair Plus Mattress Replacement System	<p>The ProCair Plus is our premium solution for high-risk patients optimal patient therapy and comfort; fully automated operational features; advanced safety features including back-up battery; and advanced error and service indicators. Available in Standard and King Single models.</p>	<p>King Single (2000 x 1050 x 200mm) Single (2000 x 880 x 200mm) Weight Capacity 40-200Kg</p>	
		<p>Pressure Reduction Overlays</p> <p>Functional overlay for pressure prevention Automatic pressure setting Used over the existing mattress for basic support Easy to operate and maintain 40-150 kg Ideally suited to home care or aged care environments, this product can be used for people at risk to high risk of pressure injury. Automatic weight sensing technology that delivers fail-safe patient comfort levels. Conform to all relevant Australian standards.</p>	<p>2000 x 880 x 130 mm Weight Capacity (kg) 40-150kg</p>	

ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
REHABILITATION EQUIPMENT				
REEXPED	Exercise Equipment			
	Exercise Pedals	Lightweight and portable, the versatile product can be used for a variety of alternative exercises, for both lower and upper body. The compact design makes it easy to move the exerciser from room to room, to take outside in the fresh air or even to take when travelling. The exerciser is designed to be placed on the floor when exercising with the feet, or placed on a table top when exercising with the hands. It features both forwards and backwards movement, to ensure you are working all muscle groups. It also offers an easily adjustable resistance level - simply turn the knob to either increase or decrease the resistance at will.	This product weighs 2.5 kilograms and measures 280 mm high, 445 mm wide and 530 mm deep.	

RESQBP

Balance Pad

Balance Pad ultrasoft speciality foam is both water- and dirt-resistant. Application: Balance training, co-ordination and reaction training, mobility training, standing stability training, functional retraining of the muscles in the lower extremities, motor-skill training to maintain balance, e.g. for age related physiological changes to the muscles and skeletal frame. The Physiomed Balance Pad - several can be laid on top of each other to increase the 'wobble effect'. The closed-cell robust and hygienic. They are multi-functional in use and are ideally suited for use as swimming floats.

Length 50cm
41cm
0.7kg

Width
Weight



Pain Management

RETU01

Tens Unit Spectrum

This easy to use, digital TENS machine offers natural, non-invasive pain relief without the use of chemicals or drugs. Transcutaneous Electrical Nerve Stimulation (TENS) delivers tiny electrical impulses to the nerve fibres that lie underneath the skin's surface, providing pain relief by blocking pain signals to the brain via the spinal cord and peripheral nervous system.

1 x TENS machine unit
1 x D size 9-volt battery
1 x AC power adaptor
2 x electrode leads
4 x self adhesive electrodes (pack of 4)
1 x instruction manual
1 x hard plastic carrying case

Large liquid crystal display and easy operation press buttons, ideal for people who have difficulty reading smaller panel inscriptions or turning smaller dials. The unit offers 32 stimulation modes for pain relief and muscle training. There are 5 adjustable TENS modes, 3 adjustable EMS modes and 24 professional presets for TENS and EMS available.



REHW0.5

Weights

HAND WEIGHT Plastic
COATED .5KG

These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.

Dumbell, Vinyl Coated,
Purple 0.5kg



REHW1.0	HAND WEIGHT Plastic COATED 1.0KG	These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.	Dumbell, Vinyl Coated, Pink 1.0 kg
	REHW1.5	HAND WEIGHT Plastic COATED 1.5KG	Dumbell, Vinyl Coated, Green 1.5kg
	Exercise Media		
RETBTHN	Theraband Yellow 5.5m		Yellow thin 5.5m
		TheraBand Resistance Exercise Bands are best used before working out for preventative stretches that prime muscles for strenuous exercises. They are an essential part of any workout or athletic activity.	
RETPXSF	Theraputty Tan Extra Soft 85g	Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems. It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.	Tan Extra Soft 85g



RETPSFT

Theraputty Yellow Soft 85g

Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems.

It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.

Yellow Soft 85g

**RETPMDM**

Theraputty Red Medium 85g

Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems.

It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.

Red Medium 85g

**RETPFRM**

Theraputty Green Firm 85g

Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems.

It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.

Green Firm 85g



REAW0.9




Weights






Ankle Weight - Vinyl 0.9kg
(This is a substitute for the
1.0kg ankle weight which is
currently not available)

Vinyl Cuff Weights with a long velcro
closure strap assures a secure and
comfortable fit even during the most
strenuous exercise program.

White 0.9kg Ankle
Weight with velcro
fastening.



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
TRANSFER AIDS				
	Handy Bars			
	TATHB	Transfer Handy Bar	Provides a safe and supportive handhold for transfers in and out of the car Multifunctional device designed to help users get in and out of a car easier.	Maximum User Weight Capacity: 320 kg 
	Leg Lifters			
TALL	Leg Lifter (Single) Pelican Brand	Rigid device with loop for foot at one end and two straps at the other for assisting independence to lift a leg into a car or bed.	Overall length: 102 cm	
TADHLL	Dual Handle Leg Lifter	A single leg lifter with 2 handles for users who cannot lift their leg with one hand.	Length 86.5 cm Weight 110g	
Slide Sheets				

TASLS00	Slide Sheet SMALL - 1000mm x 1450mm	Slippery surface both sides, use folded onto itself to provide low resistance under clients requiring assistance with repositioning. Size: 1m x 1.4m sq for bed transfers	1m x 1.4m vary	Colour may vary		
	TASLS01	Slide Sheet LARGE - 2000mm x 1450mm	Slippery surface both sides, use folded onto itself to provide low resistance under clients requiring assistance with repositioning. Size: 2m x 1.4m sq for bed transfers	2m x 1.4m vary	Colour may vary	
	Slide Boards					
	TACTP01	Car Transfer Pad - Pelican	Car slide transfer pad has 2 slip surfaces that decreases friction allowing the client to slide and pivot into the car seat.	Top Pad size 58cm x 58cm		
	TATB01	The Slider	"Assists patient transfers between wheelchair and car, chair, bed etc Strong perspex construction Rear wheel cutout, covered with a flexible plastic guard, minimises skin damage while transferring over the arch of the wheel Not to be used as a lifting device"	Length 640mm Width 350mm Weight 1.5kg Max User Weight 120kg		
TATB03	Curved Transfer Board PVC	PVC transfer board	Weight: 1.95kg. Max. user weight: 152kg			
	Turntables					


TASC01

Swivel Cushion Rota

This cushioned seat glides you in and out of cars smoothly without struggling. Placed on any seat it will swivel in both directions and may help those who suffer from arthritis, leg or back pain. It's lightweight, portable and consists of a plastic base with a fleece covered cushion. A strap secures the position of the swivel seat.

Diameter: 40cm
Height: 70mm
Weight: 1.6kg
Maximum user weight:
190kg



ITEM CODE	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
WHEELCHAIRS				
WCSP14S	Self Propelled Wheelchair			
	Self-Propelled 14" SWIFT	The Days Swift Self-Propelled 14" Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair.	Light frame 7.8kg Folding backrest (for storage or transport) Quick release wheels (for transport) Fixed arms Swing away and removable leg rests Solid tyres (no maintenance required) Extended brake handle (ease of brake application) Maximum User Weight: 80kg	

WCSP16L

Self-Propelled 16" LINK
SPECIAL ORDER

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
- Solid tyres
- 16" seat width
- Swing in/out legrests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Anti-tipper comes as standard inclusion

**WCSP16S**

Self-Propelled 16" SWIFT

The Days Swift Self-Propelled Wheelchair with Handbrakes has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair.

Light frame 7.8kg

- Folding backrest (for storage or transport)
- Quick release wheels (for transport)
- Fixed arms
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Handbrakes for added safety*
- Extended brake handle (ease of brake application)
- Maximum User Weight: 115kg



WCSP18L**Self-Propelled 18" LINK**

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
 - Solid tyres
 - 18" seat width
 - Swing in/out legrests
 - Swing back removable armrests
 - Height adjustable back upholstery
 - Adjustable castors
 - Multi position rear wheels
- Can fit stump supports
- Drum brakes on the 20"-22"
 - Anti-tipper comes as standard inclusion



WCSP18S**Self-Propelled 18" SWIFT**

The Days Swift Self-Propelled Wheelchair with Handbrakes has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair.

Light frame 7.8kg
Folding backrest (for storage or transport)
Quick release wheels (for transport)
Fixed arms
Swing away and removable leg rests
Solid tyres (no maintenance required)
Handbrakes for added safety*
Extended brake handle (ease of brake application)
Maximum User Weight: 115kg



WCSP18W**Self-Propelled 18" WHIRL**

The Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about. Features: Quick release wheels (for transport) Removable armrests Swing away and removable leg rests Solid tyres (no maintenance required)

Specifications: * Weight: 16.8kg (complete) * User Weight: 120kg * Seat Width: 450mm * Seat Depth: 430mm * Backrest Height: 440mm * Seat Height - Front: 500mm * Seat Height - Rear: 475mm * Front Castor Wheel: 8" x 1" PU * Rear Wheel: 24" x 1 3/8" PU * Overall Width: 640mm * Upholstery Nylon (Flame Retardant) * Colour: Silver Warranty: 12 months on frame & crossbars



WCSP20L

Self-Propelled 20" LINK
SPECIAL ORDER

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
- Solid tyres
- 16", 18", 20" and 22" seat widths available
- Swing in/out legrests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Drum brakes on the 20"-22"
- Anti-tipper comes as standard inclusion

**WCSP20M**

Self-Propelled Bariatric 20"
MERITS

Reinforced double cross-brace for extra strength with weight capacity of 182kg
Detachable full length armrests
Detachable swing away leg rests with oversize footplates, with optional elevating leg rests
24" x 1.3/8" puncture free foam filled tyres

Overall Width 73.5cm(20")
Overall Length 123cm
Overall Height 91 cm
Seat Width 48cm(20")
Seat Depth 46cm
Seat height from ground 49cm
Seat Back Height 40cm
Footrest Length
Adjustment 37cm – 50cm
Armrest Height 23cm
Weight Capacity 182kg
Total weight 15.5kg(20")



WCSP22M

Self-Propelled Bariatric 22"
MERITS

Reinforced double cross-brace
for extra strength with weight
capacity of 182kg
Detachable full length armrests
Detachable swing away leg
rests with oversize footplates,
with optional elevating leg
rests
24" x 1.3/8" puncture free foam
filled tyres

Overall Width 78.5cm(22")
Overall Length 123cm
Overall Height 91 cm Seat
Width 56cm(22") Seat
Depth 46cm Seat height
from ground 49cm Seat
Back Height 40cm Footrest
Length Adjustment 37cm –
50cm Armrest Height 23cm
Weight Capacity 182kg
Total weight 16.5kg(22")

**Transit Wheelchairs****WCTR14A**

Transit Wheelchair 14"
ASPIRE

The Assist 2 – is designed for
an individual or facility who
requires a robust chair that will
meet the rigorous demands of
the community and healthcare
environments.
Swing away and removable arm
rests – this facilitates user
transfers and also makes for
easier transport and storage.

SWL: 140kg
Weight: 17.3kg



WCTR16A**Transit Wheelchair 16"
ASPIRE**

The Assist 2 – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments. Swing away and removable arm rests – this facilitates user transfers and also makes for easier transport and storage.

SWL: 140kg
Weight: 16.4kg

**WCTR16E****Transit Wheelchair 16"
ESCAPE**

With a foldable aluminium frame, the Escape has been specially designed for ease of transportation and storage, making it an ideal chair for occasional use. Suitable for both indoors and out, it gives you flexibility wherever you go.

In addition to being portable, these attendant-propelled chairs offer both comfort and practicality. With a full-length armrest, padded upholstery and detachable swing-away footrests, they provide all the basic comfort necessities. Peace of mind and low maintenance are provided by the puncture-proof tyres.

Seat Height - Front (from the ground): 510mm
Seat Height - Rear (from the ground): 490
Seat Width (on the outside of the tubes): 450mm
Seat Width (on the inside of the tubes): 400mm
Seat Depth: 400mm
Width Between Arms: 400mm
Armrest Height Above Seat: 230mm
Backrest Height (from the seat): 400mm
Footrest Length (from seat to surface of footrest): 350-450mm



WCTR18S**Transit Wheelchair 18inch
SWIFT**

The Days Swift Transit Attendant Propelled Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair. Attendant propelled in one size.

Light frame 7.8kg
Folding backrest (for storage or transport)
Fixed arms
Swing away and removable leg rests
Solid tyres (no maintenance required)
Extended brake handle (ease of brake application)



WCTR18W**Transit Wheelchair 18" WHIRL**

The Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. * This sleek and stylish wheelchair is practical and will look good when you are out and about. Features: * Quick release wheels (for transport) * Removable armrests * Swing away and removable leg rests * Solid tyres (no maintenance required)

Specifications: * Weight: 16.8kg (complete) * User Weight: 120kg * Seat Width: 450mm * Seat Depth: 430mm * Backrest Height: 440mm * Seat Height - Front: 500mm * Seat Height - Rear: 475mm * Front Castor Wheel: 8" x 1" PU * Rear Wheel: 24" x 1 3/8" PU * Overall Width: 640mm * Upholstery Nylon (Flame Retardant) * Colour: Silver Warranty: 12 months on frame & crossbars

**WCTR20A****Transit Wheelchair 20" ASPIRE**

The Assist 2 – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments. Swing away and removable arm rests – this facilitates user transfers and also makes for easier transport and storage.

SWL: 140kg
Weight: 17.9kg

