






# RDNS EQUIPMENT CATALOGUE

## EXTERNAL PRESCRIBERS



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>BATHING AND TOILETING</b>					
<b>Bath Aids</b>					
BTBDWH	Yes	Bathboard with handle	Provides a support surface to place across the bath which enables safe and easy showering in a seated position. Helps when transferring in and out of the bath as this can be performed in a seated position. The included handle assists with transfers in and out of the bath as well as safety whilst seated showering.	Adjustable width clamps: 41 cm - 62 cm Length:69 cm Width: 26 cm Maximum user weight capacity: 150 kg	
BTBDBA	Yes	Bathboard - Bariatric	Designed to provide a stable seat over the bath. Constructed of reinforced PVC Anti - slip rubber pads improves safety and protects the bath surface	Adjustable width clamps: 50 cm - 66 cm Length: 66 cm Width: 34 cm Maximum user weight capacity: 190 kg	
BTBSSW	Yes	Bath Seat - Swivel	Designed for showering over bath and aiding with transfers into the bath by allowing this occur in a safe seated position.	Maximum user weight: 110kg Adjustable width to fit most baths Moulded plastic seat with drainage holes Contoured back rest Side support arms for added security 360° swivel seat Locking lever for set positioning	

BTBTBWB	Yes	Bath Transfer Bench - With Back	<p>Provides a support surface to place across the bath which enables safe and easy showering in a seated position. Helps when transferring in and out of the bath as this can be performed in a seated position.</p> <p>Plastic seat and backrest provides a non slip surface with holes for water drainage Adjustable backrest and handle to suit left or right facing baths without using tools Built-in dish and hand-held shower slot for user convenience</p>	<p>Height adjustable seat: 43 cm - 56 cm (25mm increments) Seat dimension: Width 66 cm Depth 37 cm Maximum user weight: 136kg</p>	
BTSMNS	Yes	Shower Mat / Bath Mat Non Slip	<p>With their raised ridges and unique design, the mats allow you to stand, without slipping, on a wet surface. Use on flat firm wet surface &amp; press down first so the suction caps engage. The mats grip securely and safely due to their high-grade medical rubber content and additional suction cups.</p>	<p>Machine washable at 40°. Dimension: 555mm x 555mm</p>	

**Mobile Shower Commodes**

BTMSTR6



No

Mobile Shower Commode,  
Transit, 457mm/18inch Wide,  
Aluminium

The attendant propelled mobile shower commode is lightweight and maneuverable to allow for easier toileting and showering whilst the user is seated. The chair features a lightweight, silver painted aluminium frame, which provides a rust free, high strength product. Brakes on all four wheels for security, the swing away, detachable footrests offer comfort and support whilst being transported and the armrests can be flipped back to facilitate side transfer. This chair can be used in a shower cubicle or wet room, used as a commode or over a standard toilet bowl. The front of the seat is cut-away to facilitate personal hygiene and a set insert is provided for comfort during transport. Supplied with a square commode pan as standard. The detachable back enables the shower commode chairs to be partially dismantled for more compact storage.

Overall Height: 958 – 1008mm (with castors)  
Overall Length: 550mm (without footrests) / 810mm (with footrests)  
Overall Width: 555mm  
Seat Height: 495 – 570mm  
Seat Width: 457mm  
Seat Depth: 425mm  
Aperture Size: 190mm x 360mm  
Backrest Width: 462mm  
Backrest Height (from seat): 436mm  
Front & Rear Wheel Diameter: 130mm  
Under Seat Clearance: 430 – 505mm  
Weight: 12.4kg  
Maximum User Weight: 125kg  
Warranty: 12 months



<p><b>BTMSTR1</b>      <b>No</b></p>	<p>Mobile Shower Commode, Transit, Stainless Steel, 460mm</p>	<p>Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Swingback arms with padded armpads assist side transfers and can be positioned at two different heights Adjustable safety arms for additional security Sliding footrest is textured and weight-bearing for added safety Footrest slides back under frame to assist patient transfers PVC slatted backrest Mounted on 4 x 125mm stainless steel locking castors Durable stainless steel frame Tested and conforms to Australian Standards AS/NZS 3973: 2009I</p>	<p>Seat Width - 460mm Seat Depth - 430mm Seat Height - 595 Armrest Height - 200/240mm Backrest Height - 340mm Overall Width - 540mm Max User Weight - 200kg</p>	
<p><b>BTMSTR3</b>      <b>No</b></p>	<p>Etac Swift Mobile Shower Commode, Transit, Swing Away Footplates, Adjustable Height</p>	<p>Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Soft, adjustable, non-stick, quick-drying nylon back support Seat and back support easily detachable for cleaning or storage Swing away arm and foot supports Curved foot supports with side heel protection</p>	<p>Fully height adjustable from 450 - 600mm Seat height (at 20mm intervals) Maximum user weight 135kgs</p>	

BTMSTR4

No

Mobile Shower Commode,  
Transit, Bariatric, 560mm wide,  
Stainless Steel

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove.

Wider seat widths and reinforced frame to suit the bariatric client

Swingback arms with padded armpads assist side transfers and can be positioned at two different heights

Adjustable safety arms for additional security

Sliding footrest is textured and weight-bearing for added safety

Footrest slides back under frame to assist patient transfers

PVC slatted backrest

Mounted on 4 x 125mm stainless steel locking castors

Durable stainless steel frame

Tested and conforms to Australian Standards AS/NZS 3973: 2009I

Seat Width - 560mm

Seat Depth - 430mm

Seat Height - 595mm

Armrest Height -

200/240mm

Backrest Height - 340mm

Overall Width - 640mm

Max User Weight - 300kg



BTMSTR5

No

Mobile Shower Commode,  
Transit, Bariatric 610mm wide,  
Stainless Steel

Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Wider seat widths and reinforced frame to suit the bariatric client

Swingback arms with padded armpads assist side transfers and can be positioned at two different heights

Adjustable safety arms for additional security

Sliding footrest is textured and weight-bearing for added safety

Footrest slides back under frame to assist patient transfers

PVC slatted backrest

Mounted on 4 x 125mm stainless steel locking castors

Durable stainless steel frame

Tested and conforms to Australian Standards AS/NZS 2072: 2000

Seat Width - 610mm

Seat Depth - 430mm

Seat Height - 595mm

Armrest Height -



200/240mm

Backrest Height - 340mm



Overall Width - 690mm

Max User Weight - 300kg



<p><b>BTMSSP1</b>      <b>No</b></p>	<p>Mobile Shower Commode, Self Propelled, SS, 460mm</p>	<p>Assists with the transportation of the user for showering and toileting. The chair can be wheeled over a standard toilet and used in a large shower alcove. Swingback arms with padded arm pads for side transfers Adjustable safety arms for additional security Armrests can be positioned at two different heights Removable PVC slatted backrest Adjustable height swingaway legrests or slide out weight-bearing footrest Mounted on 2 x 125mm stainless steel castors and 2 x 22 inch rear wheels Durable stainless steel frame.</p>	<p>Seat Width - 460mm Seat Depth - 430mm Seat Height - 595mm Armrest Height - 200/240mm Backrest Height - 340mm Overall Width - 660mm Max User Weight - 200kg</p>	
<p><b>BTMSSP2</b>      <b>No</b></p>	<p>Mobile Shower Commode, Self Propelled, Adjustable Height</p>	<p>Easy to maneuver, for the independent user, the 24-inch wheels are positioned to make it easy to maneuver the chair without compromising user safety or comfort. It is easy to transfer to and from the chair thanks to the large, unrestricted seat in front of the wheels. The brake handles are also positioned below the seat so that they do not get in the way during transfer, while remaining easy to reach. Front seat recess; fully open at rear Long, removable, swing away arm supports Curved footplate with side protection for heels; slides in under seat Footplate features full height and tilt adjustment.</p>	<p>Width - 69 cm Width between arm supports - 43.5 cm Seat height - 55 cm Seat width - 48 cm Seat depth - 44 cm Space below seat - without pan holder: 49.5 cm, with pan holder: 48 cm Opening - depth x width, 36 cm x 21,5 cm Weight - 20 kg Maximum user weight - 130 kg</p>	



<p><b>BTDMSC</b></p> <p><b>No</b></p>	<p>Denmark Mobile Shower Commode</p>	<p>The Denmark mobile shower/commode chair is attendant propelled and has height adjustable and angle adjustable footrests. The slip resistant armrests can be removed sideways providing an advantage for side way transfers. The Chair also has 100 mm castors, an ergonomic cutout seat for easy hygiene access and a soft flexible plastic backrest. The weight capacity is 150kgs.</p>	<p>Weight (incl. Footrests): 12.6kg  Capacity: 150kg  Seat Depth: 45cm  Seat Height: 44.5 - 56.5cm (adjustable)  Width Between Armrests: 47.5cm  Max Width: 58cm  Length (without the footrests): 60cm</p>	
<p><b>BTMSTR7</b></p> <p><b>No</b></p>	<p>Aquatec Ocean VIP Ergo Tilt Shower Commode</p>	<p>The Aquatec Ocean Ergo Family of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well established Aquatec Ocean range, the Ocean Ergo family comes packed with a range of new features, designed to make personal care activities easier for both the user and carer.</p>	<p>SWL - 150kg.  Seat Height 475 - 600mm  Seat width 480mm  Overall width 565mm  Overall Depth 1000mm (including footrests)</p>	

BTMSEST

No

ETAC Swift Mobil-2 Tilt in Space Mobile Shower Commode with Pan

The Swift Mobil Tilt-2, has a unique tilt feature that gives increased security for the user and an ergonomic working position for the carer. It's -5° to +35° tilt function can be operated from either side of the chair.

SWL - 160kg

Height 500 - 650mm Seat

Width 500mm

Width 640mm

Seat

Overall



This allows the carer to reassure the user when tilting.

### Shower Chairs

BTSCST



Yes



Shower Chair, Standard,  
Adjustable - With Arms



A standard no-fuss, safe and reliable chair which is height adjustable. It features a contoured, plastic moulded seat, which is both comfortable & hygienic. It has been cleverly designed with a gentle slope, so that water can flow away for easy drainage. The steel frame using 22.2mm tubing enables users to gain a firm grip making transferring safer and easier. Tested to Australian standards.





Overall Height : 775mm - 900mm  
Overall Width : 545mm  
Overall Depth : 515mm  
Seat Width: 440mm  
Seat Depth: 405mm  
Seat Height Range: 410mm–535mm  
Width between Armrests: 490mm  
Armrest (Seat to top of handle): 200mm  
Backrest Height: 350mm  
Frame: Steel  
Frame Diameter (armrest and leg): 22.2mm  
Extension Leg Diameter: 25.8mm  
Extension Leg (outer/inner): Outer  
Rubber Tip (White): DIA 40mm x DIA 23.8mm  
Colour: White  
Net Weight: 3.7kg  
Maximum User Weight: 125kg  
Warranty: 12 months



<p><b>BTSCB00</b>      <b>Yes</b></p>	<p>Shower Chair - Bariatric - 525 Wide</p>	<p>Reinforced wide frame – Provides stability for larger users          Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion          Height adjustable legs – Easy push button height adjustment for changing needs          Additional height increments on one leg – Accommodates uneven surfaces          All frames are individually treated for corrosion resistance          Raised armrests – Assist the user transferring on and off the elevated seat          Durable non-slip feet – Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability</p>	<p>SWL 200kg      Overall          Height: 735 - 910 mm          Overall Width: 570 mm          Overall Depth: 570 mm          Width Between Arms: 525 mm</p>	
<p><b>BTSCB03</b>      <b>Yes</b></p>	<p>Shower Chair - Bariatric - 550 Wide</p>	<p>Height adjustable          Wider seat widths and reinforced frame to suit the bariatric client          Padded seat and curved backrest provide comfort and postural support          Non-slip rubber feet          Powder coated steel frame          Electrophoretic surface treatment to prevent rusting</p>	<p>Seat Width 550mm          Seat Height 470 - 640mm          Base Width 650mm          Base Depth 520mm          Max User weight 300kg</p>	

<p>BTSCB04</p> <p>Yes</p>	<p>Shower Chair - Bariatric - 600 Wide</p>	<p>Height adjustable  Wider seat widths and reinforced frame to suit the bariatric client  Padded seat and curved backrest provide comfort and postural support  Non-slip rubber feet  Powder coated steel frame  Electrophoretic surface treatment to prevent rusting</p>	<p>Seat Width 600mm  Seat Height 470 - 640mm  Base Width 700mm  Base Depth 520mm  Max User weight 300kg</p>	
<p>BTSCB01</p> <p>Yes</p>	<p>Shower Chair - Bariatric - 610 Wide</p>	<p>Backrest can be adjusted in depth +/- 60mm and has curved tubing to ensure comfort for any body shape.  Having the right height for the user ensures success in getting in and out of the chair and these chairs are adjustable from 420 to 550mm.  Butterfly armrests ensure the body doesn't get squeezed getting in to and out of the chair whilst providing a good grip for sitting and getting up.  The high-tech, lightweight tubing provides the highest strength at the lowest weight  Removable seat pad ensures these items are easy to clean.</p>	<p>Seat Width: 610mm  Seat Depth: 400 to 520mm  Overall Width: 730mm  Weight: 7kg Maximum  User Weight: 325kg</p>	

<p><b>BTSCB05</b>      <b>Yes</b></p>	<p>Shower Chair - Bariatric - 650 Wide</p>	<p>Height adjustable  Wider seat widths and reinforced frame to suit the bariatric client  Padded seat and curved backrest provide comfort and postural support  Non-slip rubber feet  Powder coated steel frame  Electrophoretic surface treatment to prevent rusting</p>	<p>Seat Width 650mm  Seat Height 470-640mm  Base Width 750mm  Base Depth 520mm  Max User Weight 300kg</p>	
<p><b>BTSCB02</b>      <b>Yes</b></p>	<p>Shower Chair, Commode &amp; Toilet Aid, Barriatric (All in One)</p>	<p>Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.</p>	<p>Total Width: 790mm  Total Depth: 530mm  Seat Depth: 450mm  Seat Height: 465 - 615mm (to over seat top)  Width Between Arms: 600mm  Height Adjustable: 430-570mm  Safe Working Limit: 295kg</p>	

BTSCSW1	Yes	Shower Chair, Swivel	<p>Height adjustable  Swivel motion assists easy transfers from standing to showering position  Seat locks in both transfer and showering positions  Integrated arms assist transfers  Drainage slots drain water away from seating area  Non-slip rubber tips  Powder coated steel frame</p>	<p>Width Between Arms 480mm  Seat Height 480-560mm  Base Width 560mm  Base Depth 540mm  Max User Weight 160kg</p>	
<b>Shower Stools</b>					
BTSSWA	Yes	Shower Stool Adjustable - With Arms	<p>The standard shower stool is height adjustable and features side handles to assist with transfers. The stool has a small footprint and is ideal for smaller shower alcoves.</p>	<p>Overall Width : 510mm  Overall Depth : 410mm  Width Between Arms: 450mm  Weight of Stool: 3.2kg  Height Adjustable: 490 - 610mm  Maximum User Weight: 110kg</p>	
BTSSNA	Yes	Shower Stool Adjustable - Without Arms	<p>Height adjustable  Side handles to assist transfers  Durable plastic moulded seat  Non-slip rubber tips  Electrophoretic surface treatment to prevent rusting  Powder coated steel frame</p>	<p>Width 380mm  Seat height 510 - 660mm  Base width 550mm  Base Depth 380mm  Max user weight 125kg</p>	
BTSSAW2	Yes	Shower Stool Adjustable - With Backrest	<p>Height adjustable  Padded seat and back for added comfort  Integrated arms assist transfers  Non-slip rubber tips  Electrophoretic surface treatment to prevent rusting  Powder coated steel frame</p>	<p>Width Between arms 460mm  Seat Height 520-670mm  Base Width 550mm  Base Depth 380mm  Max User Weight 125kg</p>	

BTSSAHD

Yes

Shower Stool Adjustable -  
Heavy Duty Padded - Aluminium




Constructed from sturdy no rust  
aluminium with easy push button  
adjustment legs  
The internal sliding legs provide smooth  
and effortless height adjustment whilst  
avoiding any chance of water filling the  
legs  
Reinforced frame using 25.4mm tubing  
enables user to gain a firmer grip  
making transferring safer & easier  
The seat is padded for comfort and is  
removable for cleaning





Overall Height (footprint):  
679mm - 829mm  
Overall Width (footprint):  
655mm  
Overall Depth (footprint):  
445mm  
Seat Width: 445mm  
Seat Depth: 275mm  
Seat Height Range:  
539mm–689mm  
Width between Armrests:  
538mm  
Armrest (Seat to top of  
handle): 140mm  
Frame: Aluminium  
Frame Diameter (armrest  
and leg): 25.4mm  
Extension Leg Diameter:  
28.6mm  
Extension Leg  
(outer/inner): Outer  
Rubber Tip (White): DIA  
40mm x DIA 27mm  
Colour: White  
Net Weight: 3.0kg  
Maximum User Weight:  
150kg  
Warranty: 12 months  
ARTG ID: 291268  
Tested to Australian





BTSSAWB	Yes	Shower Stool With Removable Back and Arms - Bariatric	<p>Aluminum construction in light weight and corrosion free</p> <p>Cross-brace to enhance the stability during sitting</p> <p>Blow molded bench and back provide comfort and strength</p> <p>Drainage in seat to reduce slipping</p> <p>Limited Lifetime Warranty</p>	<p>Seat width 51cm</p> <p>Seat depth 36cm</p> <p>Seat height 36-51cm</p> <p>Overall width 58cm</p> <p>Overall depth 38cm</p> <p>Overall height 68-84cm</p> <p>Weight Capacity 227kg</p> <p>Total weight 3.3kg</p>	
<b>Accessories</b>					
BTSH1T	Yes	Shower Hose Single Tap 2 Metre - Push on	<p>Fixes to existing tap spouts Lightweight plastic handpiece Supplied with two rubber reducing collars to suit most tap spouts Easily fitted and removed</p>	Length 2m	
BTSH2T	Yes	Shower Hose Double Tap 2 Metre - Push on	<p>Hand Shower – Dual Tap Clamp-on. Fixes to existing tap spouts Lightweight plastic handpiece Supplied with two rubber reducing collars to suit most tap spouts Easily fitted and removed</p> <p>Length 2m</p>	Length 2m	
<b>Toileting Aids</b>					

BTTAX01	Yes	<p>Toilet Seat Extension - With Lid - Ergonomic contoured one piece seat with lid Tapered flange with front and rear access Stain and odour resistant Easy to secure with two hand adjustable brackets</p>	<p>Height 50mm. Maximum user weight 130kg</p>	
BTTAX02	Yes	<p>Toilet Seat Extension - With Lid - Ergonomic contoured one piece seat with lid Tapered flange with front and rear access Stain and odour resistant Easy to secure with two hand adjustable brackets</p>	<p>Height 100mm. Maximum user weight 130kg</p>	
BTTAX03	Yes	<p>Toilet Seat Extension - 50mm</p> <p>Increases toilet seat height Ideal for people with hip or knee injuries or those who require greater seat height Easily positioned and removed using location lip and adjustable clamps Front and rear cutouts for personal cleaning Contoured seat for greater comfort Moulded plastic is stain and odour resistant</p>	<p>Height: 50mm Max user weight: 190kg</p>	

BTTAX04	Yes	Toilet Seat Extension - 100mm	Increases toilet seat height Ideal for people with hip or knee injuries or those who require greater seat height Easily positioned and removed using location lip and adjustable clamps Front and rear cutouts for personal cleaning Contoured seat for greater comfort Moulded plastic is stain and odour resistant	Height: 100mm Max user weight: 190kg	
BT TSA01	Yes	Toilet Surround - Adjustable	Frame to assist toilet transfers Height adjustable Non-slip rubber tips Powder coated steel	Height - 625-775 Width between Arms - 455 Base Width - 540 Base Depth - 490 Max User Weight - 125kg	
BTT SW01	Yes	Toilet Surround - Wide	Heavy duty fixed height frame to assist toilet transfers Fitted with rear cross bar for additional strength and stability Non-slip rubber tips Powder coated steel	Height - 685 Width between Arms - 620 Base Width - 670 Base Depth - 520 Max User Weight - 130kg	
BTT SRST	Yes	Toilet Seat Raiser - Adjustable with Splash Guard (SWL 125kg)	The Over Toilet Aid assists transfers on & off the toilet. This toilet surround frame is height adjustable, with a push button adjustment. Powder coated steel frame with non slip rubber feet. Comes complete with a seat and splashguard	Seat Height 440-590mm Seat Width 370mm Width Between Arms 450mm Max User Weight 125kg	

BTTSR525	Yes
BTTSR55	Yes

Toilet Seat Raiser - Adjustable with Splash Guard (SWL 175kg)

Toilet Seat Raiser - Bariatric 550mm Wide



Raised armrests – Assist the user transferring on and off the Over Toilet Aid  
 Additional height increments on one leg – Accommodates uneven surfaces  
 Height adjustable legs – Easy push button height adjustment for changing needs

Provides a raised seat to assist in toileting and transfers.  
 Wider seat widths and reinforced frame to suit the bariatric client.  
 Ideal for people with hip or knee injuries or those who require greater seat height  
 Ergonomic angled arms assist transfers  
 Supplied with lift-back toilet seat and removable splashguard  
 Non-slip rubber tips  
 Powder coated steel

Seat Height 450 - 630mm  
 Seat Depth 435mm  
 Width Between Arms 525mm  
 Max User Weight 175kg

Width Between Arms 550mm  
 Seat height 550mm  
 Armrest height 230mm  
 Base width 680mm  
 Base Depth 610mm  
 Maximum user weight 200kg



<p><b>BTTSR60</b></p> <p>Yes</p>	<p>Toilet Seat Raiser - Bariatric 600mm Wide</p>	<p>Provides a raised seat to assist in toileting and transfers Wider seat widths and reinforced frame to suit the bariatric client Ideal for people with hip or knee injuries or those who require greater seat height Ergonomic angled arms assist transfers Supplied with lift-back toilet seat and removable splashguard Non-slip rubber tips Powder coated steel</p>	<p>Width Between Arms 600mm Seat height 550mm Armrest height 230mm Base width 730mm Base Depth 610mm Maximum user weight 250kg</p>	
<p><b>BTTSR65</b></p> <p>Yes</p>	<p>Toilet Seat Raiser - Bariatric 650mm Wide</p>	<p>Provides a raised seat to assist in toileting and transfers Wider seat widths and reinforced frame to suit the bariatric client Ideal for people with hip or knee injuries or those who require greater seat height Ergonomic angled arms assist transfers Supplied with lift-back toilet seat and removable splashguard Non-slip rubber tips Powder coated steel</p>	<p>Width Between Arms 650mm Seat height 550mm Armrest height 230mm Base width 780mm Base Depth 610mm Maximum user weight 300kg</p>	

BTSCB02



Yes



Shower Chair, Commode & Toilet Aid, Barriatric (All in One)

The bariatric All in One provides for being a toilet seat raise, shower chair and commode chair. The legs are height adjustable. It has a square bowl with easy access for removal from the front of the chair. The square bowl does not have a lid. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Total Width: 790mm  
Total Depth: 530mm  
Seat Depth: 450mm  
Seat Height: 465 - 615mm  
(to over seat top)  
Width Between Arms:  
600mm  
Height Adjustable: 430-  
570mm  
Safe Working Limit: 295kg



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>BEDROOM</b>					
<b>Bedside Commode</b>					
BRCMAH1	Yes	Commode - Adjustable Height	<p>The Bedside Commode has a clip on toilet seat, which clips on in seconds. There are no metal fittings on toilet seat. It is also height adjustable and easy to clean. It has a light blue deluxe padded overseat, as well as a removable round bowl with lid and handle. The lid eliminate spillages and reduces risks of cross infection. It also features a powder coated dark blue frame. This bedside commode looks very smart in a bedroom.</p>	<p>Overall Height: 820 - 970mm            Overall Width: 610mm            Overall Length: 550mm            Seat Height: 460 - 610mm            Weight: 8kg            Maximum User Weight: 110kg</p>	
BRCMCL	Yes	Classic Bedside Commode	<p>Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom            Available in Champagne or Slate Vinyl            Raised Padded Armrests            Easy Clean Vinyl Surfaces            Height Adjustable            Removable, Padded Seat            Non-Slip Rubber Tips</p>	<p>Seat Height: 475 - 615 mm</p>	

<p><b>BRCM58</b></p> <p>Yes</p>	<p>Commode - Adjustable Height - Heavy Duty - 580mm</p>	<p>Steel Heavy Duty Commode With removable Backrest. Durable and rigid frame in extra wide for high weight loading Adjustable height.</p>	<p>Seat 58cm width x 46cm depth Max User Weight 200kg.</p>	
<p><b>BRCMB65</b></p> <p>Yes</p>	<p>Commode - Bariatric - 650mm Wide</p>	<p>Wider seat widths and reinforced frame to suit the bariatric client Padded seat, padded armrests, lift-back toilet seat and removable commode bowl Adjustable seat height Padded seat cover slides back easily for access to commode Three sided vanity skirt conceals commode bowl Powder coated steel frame</p>	<p>Width between arms 650mm Seat height: 390 - 540mm Armrest Height: 300mm Max user weight 150kg</p>	



BTSCB02

Yes

Shower Chair, Commode & Toilet Aid, Barriatric (All in One)

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with square bowl with easy front access to remove bowl. The square bowl does not have a lid. Optional splash guard available. The seat lifts up and comes with a seat cover. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring.

Total Width: 790mm  
Total Depth: 530mm  
Seat Depth: 450mm  
Seat Height: 465 - 615mm (to over seat top)  
Width Between Arms: 600mm  
Height Adjustable: 430-570mm  
Safe Working Limit: 295kg



BRHBR01

Discuss with CL

### Bed Sticks



Horizontal Bed Rail




Once the user is in bed, simply slide the Bed Rail under the mattress. The weight of the user secures the Bed Rail. Constructed from sturdy epoxy coated steel.




Height: 500mm  
Width (Base): 580mm  
Length (Rail): 600mm  
Weight: 2.5kg



### Bed Positioning

<p><b>BRBB01</b>      <b>Yes</b></p>	<p><b>Bed Bolster - Mattress Extension</b></p>	<p>200 mm long foam mattress extension for additional pressure relief under heel zone, available in Standard and King Single widths.          Suitable for use in most standard residential care, hospital and home care beds with sleeping platform extension capabilities. Standard model suits any 150 mm high foam mattress; King Single model suits Triple Layer King surface.</p>	<p>King Single (180 high, 1050mm wide)          Single (150 high, 865mm wide)</p> <p>Features          high-density memory foam upper layer with medium density foam core base          double-strength multi-stretch waterproof cover reduces the potential for shear injuries          seams are stitched, welded and taped for maximum infection control          Standard and King Single widths available</p>	
<p><b>BRBRA01</b>      <b>Yes</b></p>	<p><b>Back Rest - Adjustable</b></p>	<p>Multi-angled pillow support          Provides assistance and support in bed          Wide aluminium frame fitted with strong breathable nylon back          Angle can be adjusted via a simple ratchet mechanism</p>	<p>Length 495mm          Width 700mm          Height 490mm</p>	<hr/> 

<p>BRBW01</p> <p>Yes</p>	<p>Foam Bed Wedge</p>	<p>Angled foam cushion provides comfort and support to the upper body Ideal for those who suffer from respiratory and reflux problems Removable quilted cotton cover Wedge can be turned to provide two different recline angles (15 degrees or 75 degrees) Foam insert can be removed to alter the angle by approx 5 degrees</p>	<p>Length 680mm Width 710mm Height 240-600mm</p>	
<p><b>Bed / Chair Blocks</b></p>				
<p>BRBK01</p> <p>Yes</p>	<p>BED Blocks Plastic 40mm</p>	<p>Plastic block raisers designed to increase the height of a bed. Increases height by 40 mm Stackable</p>	<p>Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg</p>	
<p>BRBK01</p> <p>Yes</p>	<p>CHAIR Blocks Plastic 40mm</p>	<p>Plastic block raisers designed to increase the height of a chair. Increases height by 40 mm Stackable</p>	<p>Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg</p>	

BRBK02	Yes	BED Blocks Plastic 100mm	Plastic block raisers designed to increase the height of a bed. Increases height by 100 mm Stackable	Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg	
BRBK02	Yes	CHAIR Blocks Plastic 100mm	Plastic block raisers designed to increase the height of a chair. Increases height by 100 mm Stackable	Internal radius: 65 mm Suits square legs up to 45 mm x 45 mm Weight Limit: Furniture plus user must weigh less than 250 kg	
<b>Overbed Tables</b>					
BROBT01	Yes	Overbed Table	Eoxy coated frame with a vinyl top. Recessed edges, which enables the table to be cleaned easily. The 50mm twin wheel castors allow the table to be easily moved away from the bed when not in use.	Length 800mm Overall Width 400mm Minimum Height 710mm Maximum Height 1050mm Height of Castor Wheels 50mm	

BROBT03

Yes

Overbed Table - Split Top Tilt

Recessed Thermoform top to prevent spills  
Lightweight mobile tilting table positions over bed or chair  
Can be used for eating, drawing or other activities  
Angle adjustable to 45 degrees, ideal for homecare environment  
Side table does not tilt, suitable for drink and food  
Fully height adjustable  
"C" base shape fitted with four castors enables close positioning to the bedside or chair  
Stylish cream colour finish  
Powder coated mild steel frame construction

Height: 780 - 1140mm

Width: 430mm

Length: 800mm



**Self Help Poles/Ladders**

BRBCA01



Yes




Bed Cradle - Adjustable Height




The bed cradle is designed to keep bed linen off sensitive areas of the patient's lower leg or feet by creating clearance between bed linen and the mattress  
The lower frame slides underneath the mattress, securing the cradle to the bed, while the upper frame is positioned horizontally over the lower limbs  
Bedding is then placed over the frame to remove any weight or pressure from the sensitive limbs below  
Height adjustable frame allows the user to customise the space between mattress surface and bedding

Height adjustable 440 - 640 mm





<p><b>BRBL01</b>      <b>Yes</b></p>	<p><b>Bed Ladder</b></p>	<p>A solution for clients who have physical difficulties sitting up in or getting out of bed  Step by step they gradually pull themselves up from a lying to a sitting position  The action involved also strengthens the stomach and arm muscles of the user  Can also be used to assist in transfers from wheelchairs to beds  Plastic ladder rungs give a firm and safe grip  Cord ends tie around the bed legs  Suitable for single or double beds</p>	<p>Rung Width 195mm  Rung Diameter 25mm  Cord Length 3m (to first rung)</p>	
<p><b>BRSHPO1</b>      <b>Yes</b></p>	<p><b>Self Help Pole - Freestanding</b></p>	<p>Free standing support handle provides assistance for the user to reposition or transfer in and out of bed  Handle height is easily adjusted to suit the user  Sits behind or beside the bed</p>	<p>Safe work load 150kg</p>	



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>MOBILITY AIDS</b>					
<b>Rollator Frames</b>					
MARFMIN	Yes	Rollator Walking Frame - Mini	Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Max height 82cm Min height 72cm. Max user weight. 150kg	
MARFST	Yes	Rollator Walking Frame - Standard 5 inch wheels	Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Adjustable height 810 - 910mm. Overall width 670mm. Width between handles 430mm Depth 500mm. Recommended mainly for indoor use Max user weight. 150kg	
MARFBA	Yes	Rollator Walking Frame - Bariatric 227kgs	Oversize frame, wider and deeper to accommodate individuals up to 227kgs Provides support to assist a person to walk who has difficulty due to balance, pain or injury. Frames comes set up with 2 wheels at the front and 2 skis at the back for easier maneuvering. Frame folds up flat for transporting.	Maximum user weight: 227kg Height adjustable: 32.5"~39.5" with 5	
<b>Pick Up Walking Frames</b>					



MAPWFMN	Yes	Pick Up Walking Frame - Mini	Height adjustable Folds flat for easy storage and transportation Foam handles Lightweight aluminium Tip size 28mm	Handle Height 700 - 805mm Max Base Size 630 x 500mm Weight 2.5kg Max user weight 150kg	
MAPWFST	Yes	Pick Up Walking Frame -- Standard	Folding aluminium frame with foam rubber hand grips.	Adjustable height 770 - 940mm with rubber ferrules on all 4 legs. Width 455mm, Depth 400mm.	
MAPWFBA	Yes	Pick-up Walking Frame - Bariatric 227kgs	Oversize frame design provides wider and deeper to accommodate individuals up to 227kgs Precision design provides additional strength while adding minimal weight	Maximum user weight 227kgs Height adjustable: 32.5"~39.5"	



**Walking Sticks**






<p><b>MAWSTH</b>      <b>Yes</b></p>	<p>Walking Sticks - T Handle</p>	<p>T shaped handle, single point, adjustable. To provide support people who have difficulty walking due to balance problems, pain or injury.</p> <p>The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Handle height: 79 cm - 102 cm (25 mm increments) Approximate user height: 5' 3" - 6' 1" Rubber tip/ferrule: 19 mm Product weight: 300 gm Maximum user weight capacity: 114 kg</p>	
<p><b>MAWSB01</b>      <b>Yes</b></p>	<p>Walking Stick - Bariatric - 225 Kg</p>	<p>This sturdy steel walking stick provides support to a person who has difficulty walking due to balance problems, pain or injury.</p> <p>The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Rubber Hand Grip Handle Height: 73 cm - 96 cm (25 mm Increments) Approximate User Height: 5' 2" - 6' 1" Rubber tip/ferrule: 22 mm Maximum user weight capacity: 225 kg</p>	



<p><b>MAWSCH</b>      <b>Yes</b></p>	<p>Walking Stick - Crook Handle</p>	<p>To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight. The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Height adjustment range: 760mm - 990mm</p>	
<p><b>MAWSEL</b>      <b>Yes</b></p>	<p>Walking Stick - Ergonomic - Left Hand</p>	<p>Anatomical Walking StickConstructed from aluminium to allow durability and strength whilst being light in weight. Supplied with anti-slip rubber tip. This walking stick is available in both a left hand and right hand option to ergonomically fit the natural shape of the user's hand for maximum comfort and ease of use.</p>	<p>Height adjustment range: 760mm - 990mm</p>	



<p><b>MAWSER</b></p> <p>Yes</p>	<p>Walking Stick - Ergonomic - Right Hand</p>	<p>Anatomical Walking Stick Constructed from aluminium to allow durability and strength whilst being light in weight. Supplied with anti-slip rubber tip. This walking stick is available in both a left hand and right hand option to ergonomically fit the natural shape of the user's hand for maximum comfort and ease of use.</p>	<p>:Height adjustment range: 760mm - 990mm</p>	
<p><b>MAWSF</b></p> <p>Yes</p>	<p>Walking Stick - Folding</p>	<p>A lightweight, aluminium walking stick with compact folding mechanism for convenient storage. The stick folds and unfolds simply and easily, making it ideal for everyday use or travel.</p> <p>The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>lightweight folding walking stick height adjustable easy to collapse and unfold maximum user weight 113kg black colour</p>	

<p><b>MAWSSN</b></p> <p>Yes</p>	<p>Walking Stick - (Swan neck) with strap</p>	<p>To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight. The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Height adjustment range: 810mm - 1035mm          Tube size: Upper tube 22.2mm, lower tube (telescopic) 19mm          Rubber tip (Outer size Ø38XH48; Inner sizeØ18.5xD34)          Maximum User Weight: 113kg          Colour: Silver</p>	
<p><b>MAWSSNS</b></p> <p>Yes</p>	<p>Walking Stick - (Swan neck) without strap</p>	<p>To provide support people who have difficulty walking due to balance problems, pain or injury. Constructed from aluminium to allow durability and strength whilst being light in weight. The height of the stick can be adjusted to match the height of the user an important feature when selecting the stick that is right for you. Adjust the height of the telescopic leg with the push of a button. The tip of the stick is covered in a solid rubber tip, which can be replaced should it show signs of wear after prolonged use.</p>	<p>Height adjustment range: 810mm - 1035mm          Tube size: Upper tube 22.2mm, lower tube (telescopic) 19mm          Rubber tip (Outer size Ø38XH48; Inner sizeØ18.5xD34)          Maximum User Weight: 113kg          Colour: Silver</p>	
<p><b>Quad Sticks</b></p>				

MAQSS01	Yes	Quad Stick Small - 15cm x 20cm	This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided,	Handle height: 74 cm - 97 cm ( 25 mm increments) Soft foam handgrip Base: 20 cm x 15 cm Rubber tips/ferrule: 13 mm Product weight: 800 gm Maximum user weight: 114 kg	
MAQSL01	Yes	Quad Stick Large - 17cm x 28cm	This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided,	Handle height: 76 cm - 98 cm ( 25 mm increments) Soft foam handgrip Base: 28 cm x 17 cm Rubber tips/ferrule: 16 mm Product weight: 800 gm Maximum user weight capacity: 113 kg	
MAQSOSL	Yes	Quad Stick - Off Set - 23 x 39cm	This lightweight aluminium quad base walking stick is suitable for a user with limited mobility who requires a four point base for extra stability. The offset foam handle can rotate for left or right hand use. It is adjustable in height, wrist strap provided, Adjustable height, larger base.	Height 795-1020mm Base width 230mm Base Depth 390mm Tip size 16mm SWL 75kg	

**4 Wheel Walkers**

<p><b>MA4WMI</b>      <b>Yes</b></p>	<p>4 Wheel Walker - CAREQUIP MINI</p>	<p>Frame and backrest fold compact for storage and transportation          Height adjustable handles          Narrow frame and low seat height ideal for the smaller client, or where space is limited          Handbrakes can be locked for added stability          Supplied with removable backrest, vinyl bag and padded seat          Walker can be folded with vinyl bag still attached          Curved padded backrest for additional comfort and support          *Note: Product comes with a vinyl bag and not a basket as pictured.</p>	<p>Seat height 470mm    Seat width 310mm    Max user weight 130kg          Unit weight 6.2kg</p>	
<p><b>MA4WSB</b>      <b>Yes</b></p>	<p>4 Wheel Walker - STANDARD BLUE</p>	<p>Frame and backrest fold compact for storage and transportation          Height adjustable handles          Soft-lock handbrakes          Curved padded backrest for additional comfort and support          Supplied with vinyl bag and padded seat          Walker can be folded with vinyl bag still attached</p>	<p>Seat height 525mm    Seat width 350mm    Max user weight 130kg          Unit weight 5.9kg</p>	

MA4WSR	Yes	4 Wheel Walker - STANDARD RED	<p>Frame and backrest fold compact for storage and transportation</p> <p>Height adjustable handles</p> <p>Soft-lock handbrakes</p> <p>Curved padded backrest for additional comfort and support</p> <p>Supplied with vinyl bag and padded seat</p> <p>Walker can be folded with vinyl bag still attached</p>	<p>Seat height 525mm</p> <p>Seat width 355mm</p> <p>Max user weight 130kg</p> <p>Unit weight 5.9kg</p>	
MA4WRE	Yes	4 Wheel Walker - DAYS RED	<p>The Days 100 Series aluminium frame ensures a combination of lightweight design and robust build quality with a maximum user weight of 125kg.</p> <p>6" wheels</p> <p>Ergonomic handgrips to make them easy to grip.</p> <p>Adjustable height handles with 5 settings at 25mm intervals</p> <p>Rest Seat</p> <p>Storage Bag</p> <p>Loop lockable brakes provide added security and convenience. The brake tension can also be easily adjusted to suit requirements</p> <p>Foldable backrest</p> <p>Folds easily for storage &amp; transport - simply fold together by pulling the strap on the seat</p>	<p>Max User Weight: 125kg</p> <p>Seat Height: 500mm</p> <p>Seat Width: 350mm</p> <p>Seat Depth: 320mm</p> <p>Width Between Arms: 460mm</p> <p>Backrest Height: 750mm</p> <p>Overall Height: 755-875mm</p> <p>Overall Width: 620mm</p> <p>Overall Depth: 640mm</p> <p>Backrest: Straight, folding</p> <p>Wheel: 6"</p> <p>Tyre: PVC</p> <p>Net Weight: 6kg</p>	

MA4WBSR

Yes

DAYS Walker Boxed -  
ADJUSTABLE Standard Red

Easy to assemble - no tools required  
Retail carrying box with handle  
Height adjustable seat and handles  
Padded seat with storage compartment  
bag underneath  
Padded backrest  
Hand operated brake system

Max User Weight: 136kg  
with foldable backrest  
Seat Height: 500 -  
600mm  
Seat Width: 356mm  
Seat Depth: 305mm  
Overall Height: 785 -  
1000mm  
Overall Width: 622mm  
Overall Depth: 680 -  
800mm  
Width Between Handles:  
460mm  
Backrest: Curved, folding  
Wheel: 7.5"  
Tyre: PVC, Black  
Net Weight: 7.4kg





MA4WBSB

Yes

DAYS Walker Boxed-  
ADJUSTABLE Standard Blue

Easy to assemble - no tools required  
Retail carrying box with handle  
Height adjustable seat and handles  
Padded seat with storage compartment  
bag underneath  
Padded backrest  
Hand operated brake system

Max User Weight: 136kg  
with foldable backrest  
Seat Height: 500 -  
600mm  
Seat Width: 356mm  
Seat Depth: 305mm  
Overall Height: 785 -  
1000mm  
Overall Width: 622mm  
Overall Depth: 680 -  
800mm  
Width Between Handles:  
460mm  
Backrest: Curved, folding  
Wheel: 7.5"  
Tyre: PVC, Black  
Net Weight: 7.4kg



MA4WLB



Yes

4 Wheel Walker - BLUE (LARGE WHEELS)

8" wheels  
Ergonomic handgrips and adjustable height handles  
Rest Seat  
Storage Bag  
Loop lockable brakes  
Foldable backrest  
Folds easily for storage & transport

Max User Weight: 125kg  
Seat Height: 565mm  
Seat Width: 350mm  
Seat Depth: 320mm  
Width Between Arms: 460mm  
Backrest Height: 800mm  
Overall Height: 840-940mm  
Overall Width: 620mm  
Overall Depth: 645mm  
Backrest: Straight, folding  
Wheel: 8"  
Tyre: PVC  
Net Weight: 7.2kg



<p><b>MA4WLG</b></p>	<p><b>Yes</b></p>	<p>4 Wheel Walker - GOLD (LARGE WHEELS)</p>	<p>8" wheels Ergonomic handgrips and adjustable height handles Rest Seat Storage Bag Loop lockable brakes Foldable backrest Folds easily for storage &amp; transport</p>	<p>Max User Weight: 125kg Seat Height: 565mm Seat Width: 350mm Seat Depth: 320mm Width Between Arms: 460mm Backrest Height: 800mm Overall Height: 840-940mm Overall Width: 620mm Overall Depth: 645mm Backrest: Straight, folding Wheel: 8" Tyre: PVC Net Weight: 7.2kg</p>	
<p><b>MA4WHD</b></p>	<p><b>Yes</b></p>	<p>4 Wheel Walker - HEAVY DUTY</p>	<p>To provide support and stability whilst walking. Seat lifts up to allow for storage in basket below 7.5" wheels Adjustable height handles Rest seat Storage bag Loop lockable brakes Folds easily for storage &amp; transport Screen reader support enabled.</p>	<p>Max User Weight: 200kg Seat Height: 570mm Seat Width: 495mm Seat Depth: 355mm Overall Height: 855 - 975mm Overall Width: 735mm Overall Depth: 655mm Backrest: Curved Wheel: 7.5" Tyre: PVC, Black Net Weight: 8.8kg Screen reader support enabled.</p>	

MA4WPD

Yes

4 Wheel Walker -PUSH DOWN BRAKES

To provide support and stability whilst walking. Seat lifts up to allow for storage in basket below. The unique brake system is activated by pushing down on the handles or whenever you sit on the seat. The robust yet lightweight alloy frame has tube-in-tube construction for increased strength and longer lifespan. The castors are made from solid rubber, making them free from punctures. They are designed for easy operation and good manoeuvrability especially on indoor surfaces. Their slightly smaller size makes them especially good for negotiating doorways or tight spaces. The orthopedic moulded handgrips are designed to contour entirely to your hand, minimising stress or strain to the wrist. The height of the handles is easy to adjust, with an extension range of 13 cm to ensure that the handles can be set to best match the height of the user. With a compact fold down design, this walker is easy to transport or store when not in use.

Suitable for people up to 130 kilograms in weight.  
Height 945mm max - 810mm min  
Length 610mm  
Width 575mm  
Seat Depth 310mm  
Seat Width 360mm  
Seat height 557mm  
Folded width 215mm  
Castor size 153mm



MA4WGST

Yes

Gutter Frame

Sturdy wheeled walking frame for indoor mobility training and rehabilitation. Designed for users who require the ability to weight bear through the forearms. Available in both non folding and folding models.



Note: Not designed to help users stand from a seated position




Forearm cuff height: 93 cm - 123 cm (25 mm increments)  
Width: Front 40 cm / Rear 60 cm  
Front swivel wheels and rear ski's for maneuverability  
Adjustable hand grips to suit users forearm length  
Maximum user weight capacity: 110 kg





---

**Tri Walkers**

<p>MA3WW</p> <p>Yes</p>	<p>Tri Walker with Tray and Basket</p>	<p>The Aluminium Tri Wheel Walkers have a lightweight aluminium frame that makes these walkers extremely easy to manoeuvre making them ideal for use indoors or out.</p>	<p>Ergonomic handgrips for optimum comfort Adjustable height handles with 5 settings at 25mm intervals Loop lockable brakes provide added security and convenience. The brake tension can also be easily adjusted to suit requirements Vinyl bag supplied as standard Folds easily for storage or transport using the simple easy to grab pull handle Supplied with basket and tray.</p>	
<p><b>Crutches</b></p>				
<p>MACRL01</p> <p>Yes</p>	<p>Crutches Axillary Large</p>	<p>Aluminium, push button adjustable height axilla crutches with rubber ferrules.</p>	<p>Maximum User Weight - 136kg 135 - 155cm</p>	

MACRM01	Yes	Crutches Axillary Medium	Aluminium, push button adjustable height axilla crutches with rubber ferrules.	Maximum User Weight - 136kg 114 - 135cm	
MACRS01	Yes	Crutches Axillary Youth/Small	Durable and strong, these Underarm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Fitted with high quality rubber tips for stability. Colour Grey. Maximum User Weight - 136kg	Maximum User Weight - 136kg 94 - 114cm	
MACREM1	Yes	Elbow Crutches Medium	Durable and strong, these Forearm Crutches are constructed from high tensile extruded aluminium, making them both lightweight and very strong. Fitted with high quality rubber tips for stability.	Maximum User Weight - 160kg 89 - 120cm	

<p><b>MACRERM</b>      <b>Yes</b></p>	<p>Elbow Crutches Medium Ergonomic Grip</p>	<p>Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Ergonomic palm grip offers pressure distribution for maximum comfort. Fitted with high quality rubber tips for stability. Both hand to floor and palm grip to arm cuff sections are adjustable.</p>	<p>Maximum User Weight - 160kg Hand Grip Height - 66 - 94cm, Arm Ring Adjustment- 20.5 - 27.5cm</p>	
<p><b>MACRERL</b>      <b>Yes</b></p>	<p>Elbow Crutches Large Ergonomic Grip</p>	<p>Durable and strong, these Ergonomic Grip Forearm Crutches are constructed from high tensile extruded aluminium making them both lightweight and very strong. Ergonomic palm grip offers pressure distribution for maximum comfort. Fitted with high quality rubber tips for stability. Both hand to floor and palm grip to arm cuff sections are adjustable.</p>	<p>Maximum User Weight - 160kg Hand Grip Height - 83 - 111cm, Arm Ring Adjustment- 20.5 - 27.5cm</p>	
<p style="text-align: center;"><b>Electric Scooters</b></p>				



PMAES745

No

745 Plus

Stylish eye catching design with attractive modern colour range

This new design offers generous foot space & leg room, also with large battery capacity

Loop handlebars offer comfortable finger or thumb throttle operation

Adjustable high back contoured seat with simple operation of slide and swivel mechanisms to allow maximum user comfort

Infinite adjustment of steering tiller allowing everyone to easily find a comfortable position

High mount charger socket means easy access when charging

11" wheels which provide extra ground clearance

Adjustable suspension to give maximum comfort

4 pole motor delivering more power

Adjustable power saving LED headlight, indicators front and rear, rear tail lights and brake lights

Built in speed reduction switch which is activated when turning offering exceptional safety.

Max speed: 10kph

Range: 30-40km

Weight: 98kg

SWL: 150kg



PMAESF

No

Fende

Stylish new modern design with a modesty panel, giving this scooter a smart and attractive look

Adjustable prestige car-like front & rear suspension delivers a smooth ride over varied terrain

Comfortable and contoured leather feel high back seat with slide, swivel, recline and adjustable armrests

Unique battery storage design, allowing for greater floor space and a secure lockable storage area for personal items

New modern LCD dashboard display, making reading the battery, speed, time & temperature gauges clear and easy

LED lighting system delivers brighter illumination with less power consumption

Built-in speed reduction switch, which is activated when turning offering exceptional safety

With new easy to access 12 volt socket (for phone or GPS)

Max speed: 10kph

Range: 30km

Weight: 114.5kg

SWL: 136kg



MACRM01


No

Pegasus Pro

Comfort, safety and reliability with added performance. Ideal for longer trips into town or afternoons driving around the park with your grandchildren, the new ergonomically designed Pegasus PRO is a powerful upgrade to our standard Pegasus METRO model.

Max speed: 10kph  
Range: 52km  
Weight: 136kg  
SWL: 160kg



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>WHEELCHAIRS</b>					
WCSP14S	Yes	Self-Propelled 14" SWIFT	<p data-bbox="535 320 902 352" style="color: red;"><b>Self Propelled Wheelchair</b></p> <p>The Days Swift Self-Propelled 14" Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors &amp; the light construction and ease to fold for storage or transport make it a very versatile chair.</p>	<p>Light frame 7.8kg</p> <p>Folding backrest (for storage or transport)</p> <p>Quick release wheels (for transport)</p> <p>Fixed arms</p> <p>Swing away and removable leg rests</p> <p>Solid tyres (no maintenance required)</p> <p>Extended brake handle (ease of brake application)</p> <p>Maximum User Weight: 80kg</p>	

WCSP16L

Yes

Self-Propelled 16" LINK  
**SPECIAL ORDER**

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
- Solid tyres
- 16" seat width
- Swing in/out legrests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Anti-tipper comes as standard inclusion



WCSP16S

Yes

Self-Propelled 16" SWIFT

The Days Swift Self-Propelled Wheelchair with Handbrakes has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair.

Light frame 7.8kg

- Folding backrest (for storage or transport)
- Quick release wheels (for transport)
- Fixed arms
- Swing away and removable leg rests
- Solid tyres (no maintenance required)
- Handbrakes for added safety\*
- Extended brake handle (ease of brake application)
- Maximum User Weight: 115kg



WCSP18L

Yes

Self-Propelled 18" LINK

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
  - Solid tyres
  - 18" seat width
  - Swing in/out legrests
  - Swing back removable armrests
  - Height adjustable back upholstery
  - Adjustable castors
  - Multi position rear wheels
- Can fit stump supports
- Drum brakes on the 20"-22"
  - Anti-tipper comes as standard inclusion



WCSP18S

Yes

Self-Propelled 18" SWIFT

The Days Swift Self-Propelled Wheelchair with Handbrakes has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair.

Light frame 7.8kg  
Folding backrest (for storage or transport)  
Quick release wheels (for transport)  
Fixed arms  
Swing away and removable leg rests  
Solid tyres (no maintenance required)  
Handbrakes for added safety\*  
Extended brake handle (ease of brake application)  
Maximum User Weight: 115kg



WCSP18W

Yes

Self-Propelled 18" WHIRL

The Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. This sleek and stylish wheelchair is practical and will look good when you are out and about. Features: Quick release wheels (for transport) Removable armrests Swing away and removable leg rests Solid tyres (no maintenance required)

Specifications: \* Weight: 16.8kg (complete) \* User Weight: 120kg \* Seat Width: 450mm \* Seat Depth: 430mm \* Backrest Height: 440mm \* Seat Height - Front: 500mm \* Seat Height - Rear: 475mm \* Front Castor Wheel: 8" x 1" PU \* Rear Wheel: 24" x 1 3/8" PU \* Overall Width: 640mm \* Upholstery Nylon (Flame Retardant) \* Colour: Silver Warranty: 12 months on frame & crossbars





WCSP20L

Yes

Self-Propelled 20" LINK  
**SPECIAL ORDER**

Days Premium Link Wheelchair is a lightweight adaptable wheelchair suitable for full-time users and is designed for maximum practicality and adjustability.

Light weight aluminium - 17.2kg

- Quick release wheels
- Solid tyres
- 16", 18", 20" and 22" seat widths available
- Swing in/out legrests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Drum brakes on the 20"-22"
- Anti-tipper comes as standard inclusion



WCSP20M

Yes

Self-Propelled Bariatric 20"  
MERITS

Reinforced double cross-brace for extra strength with weight capacity of 182kg  
Detachable full length armrests  
Detachable swing away leg rests with oversize footplates, with optional elevating leg rests  
24" x 1.3/8" puncture free foam filled tyres

Overall Width 73.5cm(20")  
Overall Length 123cm  
Overall Height 91 cm  
Seat Width 48cm(20")  
Seat Depth 46cm  
Seat height from ground 49cm  
Seat Back Height 40cm  
Footrest Length Adjustment 37cm – 50cm  
Armrest Height 23cm  
Weight Capacity 182kg  
Total weight 15.5kg(20")



WCSP22M

Yes

Self-Propelled Bariatric 22"  
MERITS

Reinforced double cross-brace for extra strength with weight capacity of 182kg  
Detachable full length armrests  
Detachable swing away leg rests with oversize footplates, with optional elevating leg rests  
24" x 1.3/8" puncture free foam filled tyres

Overall Width 78.5cm(22")  
Overall Length 123cm  
Overall Height 91 cm  
Seat Width 56cm(22")  
Seat Depth 46cm  
Seat height from ground 49cm  
Seat Back Height 40cm  
Footrest Length Adjustment 37cm – 50cm  
Armrest Height 23cm  
Weight Capacity 182kg  
Total weight 16.5kg(22")



WCSP16B1

Yes

Self-Propelled Breezy 16" with solid Tyres and Fixed armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

Seat Width: 38 - 52cm  
Seat Depth: 41 - 46cm  
Seat Height: 42.5 - 51cm  
Seat Angle: 0 - 6 degrees  
Backrest Height: 41 - 46cm  
SWL: 125kg  
Total Weight: 14.9Kg



Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

WCSP16B2

Yes

Self-Propelled Breezy 16" with pneumatic Tyres and Fixed armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Seat Width: 38 - 52cm  
Seat Depth: 41 - 46cm  
Seat Height: 42.5 - 51cm  
Seat Angle: 0 - 6 degrees  
Backrest Height: 41 - 46cm  
SWL: 125kg  
Total Weight: 14.9Kg



WCSP16B3

Yes

Self-Propelled Breezy 16" with solid Tyres and Adjustable armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Seat Width: 38 - 52cm  
Seat Depth: 41 - 46cm  
Seat Height: 42.5 - 51cm  
Seat Angle: 0 - 6 degrees  
Backrest Height: 41 - 46cm  
SWL: 125kg  
Total Weight: 14.9Kg



WCSP18B1

Yes

Self-Propelled Breezy 18" with solid Tyres and Fixed armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Seat Width: 38 - 52cm  
Seat Depth: 41 - 46cm  
Seat Height: 42.5 - 51cm  
Seat Angle: 0 - 6 degrees  
Backrest Height: 41 - 46cm  
SWL: 125kg  
Total Weight: 14.9Kg



WCSP18B2

Yes

Self-Propelled Breezy 18" with pneumatic Tyres and Fixed armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

Seat Width: 38 - 52cm  
Seat Depth: 41 - 46cm  
Seat Height: 42.5 - 51cm  
Seat Angle: 0 - 6 degrees  
Backrest Height: 41 - 46cm  
SWL: 125kg  
Total Weight: 14.9Kg



WCSP18B3

Yes

Self-Propelled Breezy 18" with solid Tyres and Adjustable armrests.

The Breezy BasiX<sup>2</sup> is a lightweight folding wheelchair that features amazing adjustability and rigidity, combined with a modern and clean look. Adjustments and stump supports available for amputee set up.

SWL: 125kg

Weight: 14.9kg

Seat Width: 38 - 52cm

Seat Height: 42.5 - 51cm

Seat Depth: 41 - 46cm



Recent Improvements include greater seat depth adjustment (41cm to 46cm), a tension adjustable backrest which is extendable from 41cm to 46cm and a more robust, rigid castor attachment (0, 3 and 6 degrees).

**Transit Wheelchairs**

WCTR14A

Yes

Transit Wheelchair 14"  
ASPIRE

The Transit 2 – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments.  
Swing away and removable arm rests – this facilitates user transfers and also makes for easier transport and storage.

SWL: 140kg  
Weight: 16.4kg



WCTR16A

Yes

Transit Wheelchair 16"  
ASPIRE

The Transit 2 – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments.  
Swing away and removable arm rests – this facilitates user transfers and also makes for easier transport and storage.

SWL: 140kg  
Weight: 16.6kg





WCTR16E

Yes

Transit Wheelchair 16"  
ESCAPE

With a foldable aluminium frame, the Escape has been specially designed for ease of transportation and storage, making it an ideal chair for occasional use. Suitable for both indoors and out, it gives you flexibility wherever you go.

In addition to being portable, these attendant-propelled chairs offer both comfort and practicality. With a full-length armrest, padded upholstery and detachable swing-away footrests, they provide all the basic comfort necessities. Peace of mind and low maintenance are provided by the puncture-proof tyres.

Seat Height - Front (from the ground): 510mm  
Seat Height - Rear (from the ground): 490  
Seat Width (on the outside of the tubes): 450mm  
Seat Width (on the inside of the tubes): 400mm  
Seat Depth: 400mm  
Width Between Arms: 400mm  
Armrest Height Above Seat: 230mm  
Backrest Height (from the seat): 400mm  
Footrest Length (from seat to surface of footrest): 350-450mm



WCTR18S

Yes

Transit Wheelchair 18inch  
SWIFT

The Days Swift Transit Attendant Propelled Wheelchair has been designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Its compact dimensions allow good manoeuvrability indoors & the light construction and ease to fold for storage or transport make it a very versatile chair. Attendant propelled in one size.

Light frame 7.8kg  
Folding backrest (for storage or transport)  
Fixed arms  
Swing away and removable leg rests  
Solid tyres (no maintenance required)  
Extended brake handle (ease of brake application)



WCTR18W

Yes

Transit Wheelchair 18" WHIRL

The Healthcare Whirl is a versatile wheelchair with a range of features to accommodate different user needs. \* This sleek and stylish wheelchair is practical and will look good when you are out and about. Features: \* Quick release wheels (for transport) \* Removable armrests \* Swing away and removable leg rests \* Solid tyres (no maintenance required)

Specifications: \* Weight: 16.8kg (complete) \* User Weight: 120kg \* Seat Width: 450mm \* Seat Depth: 430mm \* Backrest Height: 440mm \* Seat Height - Front: 500mm \* Seat Height - Rear: 475mm \* Front Castor Wheel: 8" x 1" PU \* Rear Wheel: 24" x 1 3/8" PU \* Overall Width: 640mm \* Upholstery Nylon (Flame Retardant) \* Colour: Silver Warranty: 12 months on frame & crossbars



WCTR20A

Yes

Transit Wheelchair 20" ASPIRE

The Transit 2 – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments.

Swing away and removable arm rests – this facilitates user transfers and also makes for easier transport and storage.

SWL: 140kg  
Weight: 17.9kg



**Stump Supports**

WCAMSSRB

Yes

BREEZY Stump Support -  
RIGHT

This wheelchair Accessory  
stump support replaces the  
footrest in the Breezy Basix  
wheelchair.






WCAMSSLB



Yes

BREEZY Stump Support -  
LEFT

This wheelchair Accessory  
stump support replaces the  
footrest in the Breezy Basix  
wheelchair.



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>CHAIRS AND SEATING</b>					
<b>Foot Stools</b>					
CSFS01	Yes	Foot Stool / Leg Rest	The Legrest Stool is height adjustable, designed with a 10 degree angle for comfort and support and should be placed with the higher end against the chair.	Seat Pad 400mm x 500mm x 40mm Height Adjustable 420mm - 570mm Maximum User Weight 160kg	
<b>Kitchen Stools</b>					
CSKPSWA	Yes	Kitchen Perching Stool - With Arms	Ideal for use in the kitchen, or any activity where a perched seating position is of assistance Height adjustable Padded seat and backrest Integrated arm supports assist transfers Vinyl upholstery	Width Between Arms: 455mm Seat Depth: 300mm Seat to Top of Backrest: 300mm Leg Height Adjustment: 500mm Floor to Seat: 635mm Maximum User Weight: Up to 100kg	
<b>Utility Chairs</b>					
CSLBCWA	Yes	Low Back Utility Chair	The Standard Utility Chair features a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder coated steel frame. The adjustable height assists users who struggle to get out of standard chairs.	Overall Width 580mm Seat Width 455mm Seat Depth 460mm Seat Height 450 - 600mm Back Height 465mm Maximum User Weight 160kg	

<p>CSHBCWA</p> <p>Yes</p>	<p>High Back Utility Chair</p>	<p>This chair can be used in a lounge room instead of a couch or regular armchair that is too low to get in and out of with ease. Ffeatures a contoured backrest, padded armrests, vinyl upholstery, adjustable height, and a powder coated steel frame. Available in Teal/Blue and Mocha.</p>	<p>Overall Width 640mm            Seat Width 490mm            Seat Depth 480mm            Seat Height 455 - 605mm            Back Height 730mm            Maximum User Weight 160kg</p>	
<p>CSHBCB</p> <p>Yes</p>	<p>High Back Utility Chair - Bariatric</p>	<p>This chair can be used in a lounge room instead of a couch or regular armchair that is too low to get in and out of with ease.            Wider seat widths and reinforced frames to suit the bariatric client            Adjustable height and seat depth            Contoured backrest            Powder coated steel frame</p>	<p>Upholstery Champagne Vinyl            Seat Width 600mm            Seat Depth 490-550mm            Seat Height 430-580mm            Arm Height 250mm            Back Height 740mm            Overall Width 800mm            Max User Weight 300kg</p>	
<p><b>Adult Postural Chairs</b></p>				

CSAPCL

No

Chair Air Deluxe Bed Adult with Footrest - Large

Cushioned adjustable air cell support has been incorporated into the design. The four lateral support wings fold away for easy transferring of the patient. This chair has a reclinable seat along with an adjustable backrest and calf rest, tray, neck support and footplate. The frame is completely constructed from sturdy stainless steel.

Features:



- Adjustable calf rest
- Reclinable seat
- Air system for pressure relief
- Full stainless steel frame
- Three lockable castors and one steering castor (non-marking)
- Breathable polyurethane cover (fire retardant)
- Replaceable and easy to clean cover
- Four folding wings
- Standard dining tray
- Magazine pocket

Seat Height: 63cm  
Seat Width: 51cm  
Seat Depth: 50cm  
Overall Width: 73cm  
Overall Length: 135 - 188cm  
Table to Floor: 88cm  
Table Slide: 23cm  
Table Size: 60 x 40.5 x 2cm  
Back Reclining: 102° - 136°  
Seat Tilting: 0 - 18°  
Leg Elevating: 70° down to 5° up  
Weight of Chair: 43.5kg  
Maximum User Weight: 120kg  
Warranty: 1 year



---

**Electric Recliner Chairs**

<p>CSPLRA1</p> <p>No</p>	<p>Monarch Chair ( PETITE SINGLE MOTOR)</p>	<p>Seamed backrest provides additional lumbar support  Padded legrest is synchronised with the backrest  Full chaise pad provides additional legrest support and comfort  Supplied with magazine pocket  Hand held controller is simple to use  Rise or recline motion stops at any point when button is released</p>	<p>Petite  Seat Width 480mm  Seat Depth 400mm  Seat Height 460mm  Arm Height 140mm  Back Height 620mm  Overall Width 780mm  Max User Weight 130Kg</p>	
<p>CSPLRA3</p> <p>No</p>	<p>Monarch Chair (SINGLE MOTOR)</p>	<p>Seamed backrest provides additional lumbar support  Padded legrest is synchronised with the backrest  Full chaise pad provides additional legrest support and comfort  Supplied with magazine pocket  Hand held controller is simple to use  Rise or recline motion stops at any point when button is released</p>	<p>Standard  Seat Width 550mm  Seat Depth 460mm  Seat Height 500mm  Arm Height 190mm  Back Height 680mm  Overall Width 820mm  Max User Weight 130Kg</p>	



CSPLRA2

No

Monarch Chair (DUAL MOTOR)

Seamed backrest provides additional lumbar support  
Padded legrest is synchronised with the backrest  
Full chaise pad provides additional legrest support and comfort  
Supplied with magazine pocket  
Hand held controller is simple to use  
The dual motor provides independent control of backrest and legrest.

Standard  
Seat Width 550mm  
Seat Depth 460mm  
Seat Height 500mm  
Arm Height 190mm  
Back Height 680mm  
Overall Width 820mm  
Max User Weight 130Kg



CSSMMFB

No

Mayfair Chair - Black (SINGLE MOTOR - SPACE SAVER)

The Select Mayfair Chair features a very comfortable 3 Section Waterfall Backrest Design for Optimum Comfort and Upper, Middle and Lower back support.

Quality Foams and no Sag springs provide supreme Postural Support.

Backup Power Supply.

Featuring quality Actuators and American Branded Mechanisms means quiet motors for a Smooth & Discreet Operation.

Standard  
Seat Width 520mm  
Seat Depth 480mm  
Seat Height 480mm  
Overall Width 850mm  
Max User Weight 150Kg



CSSMMFV

No

Mayfair Chair - Vanilla (SINGLE MOTOR - SPACE SAVER)

The Select Mayfair Chair features a very comfortable 3 Section Waterfall Backrest Design for Optimum Comfort and Upper, Middle and Lower back support.

Quality Foams and no Sag springs provide supreme Postural Support.

Backup Power Supply.

Featuring quality Actuators and American Branded Mechanisms means quiet motors for a Smooth & Discreet Operation.

Standard  
Seat Width 520mm  
Seat Depth 480mm  
Seat Height 480mm  
Overall Width 850mm  
Max User Weight 150Kg



CSDPRPT

No

Drive Portland Recliner PETITE  
(DUAL MOTOR)

The Portland Riser Recliner Lift Chair is designed with comfort in mind. It features an exclusive memory foam seat and its dual motor mechanism allows the user to find their perfect comfort position in minutes.

Height: 101cm  
Width: 83cm  
Depth: 90cm  
Seat Height: 43cm  
Seat Width: 47cm  
Seat Depth: 48cm  
SWL: 150kg



The Portland Riser Lift Chair Armchair features luxurious Sofia fabric and is available in two colours- Mushroom and Graphite.

Features of the Portland Riser Lift Chair include:

Reclines to an excellent snooze position.

High quality upholstery.

Excellent quiet and smooth actions.

Built in convenient storage pockets.

Battery back up facility

New 2 year warranty

Weight capacity 150kg

CSDPRTM

No

Drive Portland Recliner (DUAL MOTOR)

The Portland Riser Recliner Lift Chair is designed with comfort in mind. It features an exclusive memory foam seat and its dual motor mechanism allows the user to find their perfect comfort position in minutes.

Height: 110cm  
Width: 86cm  
Depth: 94cm  
Seat Height: 438cm  
Seat Width: 51cm  
Seat Depth: 54cm  
SWL: 150kg



The Portland Riser Lift Chair Armchair features luxurious Sofia fabric and is available in two colours- Mushroom and Graphite.

Features of the Portland Riser Lift Chair include:

Reclines to an excellent snooze position.

High quality upholstery.

Excellent quiet and smooth actions.

Built in convenient storage pockets.

Battery back up facility

New 2 year warranty

Weight capacity 150kg

**Step Stools**

CSSS01

No

Step Stool with optional hand rail.

Chrome plated steel frame step stool with hand rail. It has a slip-resistant fluted rubber covered step and a tall hand rail can be fitted to one side. It also has non-marking, slip resistant rubber tips.

Height of Step: 23cm

- Height of Handrail (from floor): 104cm

- Overall Width: 44cm

- Overall Depth: 35.5cm




- Step Width: 36cm





- Step Depth: 29cm

- Max User Weight:


160kg



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>CONTINENCE</b>					
<b>Bed Pans</b>					
CNBPSL1	Yes	Bed Pan - Slipper	Assists toileting of user whilst confined to bed Wedge shape assists in correct positioning Poly Propylene construction		
BTBPP	Yes	Bed Pan Plastic	Light weight polypropylene bed pan for toileting while remaining in bed  Wide opening cut-out area to allow for correct positioning Durable moulded plastic design with autoclavable construction Compatible with most Mobile Shower Commodes		
<b>Urinals</b>					
COMUNS1	Yes	Urinal Male - Non spill	Autoclavable for sanitising Integrated handle Patented non-spill spout prevents spillage in any position Spout unscrews easily for disposal and cleaning 100ml graduation markings Moulded plastic	Weight 240g 1 litre capacity 500ml non-spill capacity	

COMUST1	Yes	Urinal Male - Standard	Standard urinal bottle for males. Ideal for bed confinement or travelling	Graduated up to 1000 ml Secured lid to reduce spills and odours	
COFUST	Yes	Urinal Female	Assists the user or carer in personal toileting  Ideal for bed confinement or travelling Anatomically shaped for females Secured lid to reduce spills and odours	Graduated up to 1000 ml	
COASTI	Yes	<b>Absorbent Sheets</b> Absorbent Sheet - With Bed Tuck Ins	Highly absorbent draw sheet with side tuck-ins. This item comes with built in waterproofing.	Width (not incl. tuck-ins): 1000mm Length: 1000mm Absorption Capacity: 2.5L	
COWPTI	Yes	<b>McIntoshes</b> Waterproof Sheeting / McIntosh PVC	Mack-E Mattress Protector with wings		



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>ELECTRIC BEDS</b>					
<b>Electric Adjustable Beds</b>					
BREBA01	No	Community Electric Bed	<p>A versatile and adaptable nursing bed which is easy to use and built for strength and durability for the home care market. Electric bed, with electric backrest/headraise, knee break and manual leg lift. Trendelenburg and Reverse Trendelenburg. 4 locking castors. Bed folds for transportation.</p> <ul style="list-style-type: none"> <li>• Easy to fold and superior clearance when folded for transport or storage</li> <li>• Multiple positioning to suit a range of care needs</li> <li>• Easy to use hand controller with the ability to lockout individual functions</li> <li>• Linak actuators providing superior reliability</li> <li>• Strong and durable frame with a 5-year warranty</li> </ul>	<ul style="list-style-type: none"> <li>• Safe Working Load - 220kg</li> <li>• Dimensions - 2324x1050mm</li> <li>• Height Range - 330 to 730mm</li> <li>• Mattress size - 1960</li> <li>• Bed Frame Weight - 74 kg</li> <li>• Warranty - 5 year bed frame / 2 year electrical / 1 year castor</li> </ul>	

BREBA02

No

Electric Bed, Aged Care

A fully-featured, easy to operate, electric bed.

Four-part, high low profiling bed

Adjustable bed height, backrest, knee break and heel lift

Quick-release head and foot panels

Choice of panel designs

Central locking castors

In-built 200mm base extension

Electric hand controller with auto contour and a lockout function

4 Detachable mattress retainers

Drainage bag holder

- Safe working load - 220 kg

- Maximum patient weight - 185 kg

- Warranty - 3 year

- Electrical / 10 year Frame / 1 year Back-Up Battery

- Frame material - Powder Coated Steel

- Model - Standard, King

- Overall dimensions - W 920 mm† L 2036 mm± W 1090 mm† L 2036 mm±

- Height base - 330 - 780 mm 330 - 780 mm

- Mattress dimensions

- recommended Single - W 860 - 900 mm

- Mattress dimensions

- recommended Single - L 1960 - 2200\* mm

- Mattress dimensions

- recommended King - W 1050 -1070 mm

- Mattress



BREBICPO

No

iCare King Single Complete Bed Package ONYX

Safe working load 200kg on king single.

Hand Controller Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.

Hard wearing commercial grade fabric. Onyx range (BLACK) 82% Polyester, 18% Olefin. Stone range (GREY) 70% Polyester, 30% Acrylic.

Strong metal mesh under mattress allowing increased airflow

8x heavy duty swivel castors with one lockable castor on each corner.

Head lift angle 62° and knee break angle 42°

Max height 66cm from floor to top of base

Min height 22cm from floor to top of base.

Includes the IC20 Mattress a mid range mattress with a soft feel but excellent support and a matching headboard in ONYX fabric.

Safe working load - 200kg.

Minimum height from floor to base is 220mm.



BREBILP

No

iCare Long Single Complete  
Bed Package STONE

Safe working load 200kg on long  
single.  
Hand Controller Backlit for night time  
use. Easy press buttons. Wired stretch  
cord. Hook attachment for hanging to  
side of bed. Wipeable.  
Hard wearing commercial grade fabric.  
Onyx range (BLACK) 82% Polyester,  
18% Olefin. Stone range (GREY) 70%  
Polyester, 30% Acrylic.  
Strong metal mesh under mattress  
allowing increased airflow  
8x heavy duty swivel castors with one  
lockable castor on each corner.  
Head lift angle 62° and knee break  
angle 42°  
Max height 66cm from floor to top of  
base  
Min height 22cm from floor to top of  
base.  
Includes the IC20 Mattress a mid range  
mattress with a soft feel but excellent  
support and a matching headboard in  
STONE fabric.

Safe working load -  
200kg.  
Minimum height from  
floor to base is 220mm.



BREBICP

No

iCare King Single Complete Bed Package STONE

Safe working load 200kg on king single.

Hand Controller Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable.

Hard wearing commercial grade fabric. Onyx range (BLACK) 82% Polyester, 18% Olefin. Stone range (GREY) 70% Polyester, 30% Acrylic.

Strong metal mesh under mattress allowing increased airflow

8x heavy duty swivel castors with one lockable castor on each corner.

Head lift angle 62° and knee break angle 42°

Max height 66cm from floor to top of base

Min height 22cm from floor to top of base.



Includes the IC20 Mattress a mid range mattress with a soft feel but excellent support and a matching headboard in STONE fabric.





Safe working load - 200kg.


Minimum height from floor to base is 220mm.



**Accessories**

BREBBS1	No	Bedstick Low - Townsend for Electric Beds	<p>Clamped firmly to bed base railing Acts as a support to roll over, to sit from lying position and to assist with getting in and out of bed.</p> <p>For use for Hospital Style Beds</p> <p>Manufactured in Adelaide, Powdercoated mild steel. Clean using mild detergent and water and dry with soft cloth.</p>	<p>Low profile 525 x 175 2.1kg Max User Weight 136kg</p>	
BREBBS2	No	BedstickStandard - Townsend for Electric Beds	<p>Clamped firmly to bed base railing Acts as a support to roll over, to sit from lying position and to assist with getting in and out of bed.</p> <p>For use for Hospital Style Beds</p> <p>Manufactured in Adelaide, Powdercoated mild steel. Clean using mild detergent and water and dry with soft cloth.</p>	<p>Standard Height 715x175 2.4kg Max User Weight 136kg</p>	

BREBBR3	No	Bedrail for Elec Bed - Fold Down - 3/4 length	Clamp on, fold down metal side rail system attaches firmly to bed frame. Three quarter length rails for increased access.	
BREBBR1	No	Bed Rail - Full Length for Electric Bed	Drop down metal side rail system attaches firmly to bed frame. Full length rails for added safety and security.	
BREBBR2	No	Bed Rail Cover - Padded for Electric Bed	Soft cushions to encase metal bed rails protects limbs from hitting rails or getting caught. Soft, medical grade foam cushions, stretch cover, easy secure Velcro closure, easy wipe clean. Two sizes available, fits majority of bed rails.	
BREBSHP	No	Pole Self Help For Community Bed	Support pole to aid patient repositioning while lying in bed. Mounts into inbuilt brackets at either side of bed head.	

ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>PRESSURE CARE</b>					
PCROL01	No	<b>Pressure Cushions - Air</b> Cushion ROHO Low Profile 8x9 (fits 14x16 W/C)	<p>The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.</p> <p>Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.</p>	Cell size 8 x 9 Dimensions 38 x 43 cm Profile Low Height 60mm Max User Weight Unlimited	



PCROL02

No

Cushion ROHO Low Profile 9x9

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 9 x 9  
Dimensions 43 x 43 cm  
Profile Low  
Height 60mm  
Max User Weight  
Unlimited



PCROL03

No

Cushion ROHO Low Profile 10 cells x 9 cells (fits 18x16 W/C)

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 10 x 9  
Dimensions 47 x 43cm  
Profile Low  
Height 60mm  
Weight Unlimited



PCROL04

No

Cushion ROHO Low Profile 10 cells x 10 cells

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 10 x 10  
Dimensions 47 x 47 cm  
Profile Low  
Height 60mm  
Max User Weight  
Unlimited



PCROL05

No

Cushion ROHO Low Profile 11 cells x 10 cells (fits 20x18 W/C)

The Low-Profile Cushions effectively redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO Low-Profile Cushion maximizes stability while reducing the risk of potential tissue breakdown.

Available with either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use. The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

Cell size 11 x 10  
Dimensions 52 x 47 cm  
Profile Low  
Height 60mm  
Max User Weight Unlimited



PCROH01

No

Cushion ROHO High Profile 10 cells x 10 cells

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Cell size 10 x 10  
Dimensions 47 x 47 cm  
Profile High  
Height 100mm  
Max User Weight Unlimited



Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

PCROH02

No

Cushion ROHO High Profile 11 cells x 10 cells (fits 20x18 W/C)

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Cell size 11 x 10  
Dimensions 52 x 47 cm  
Profile High  
Height 100mm  
Max User Weight Unlimited  
Screen reader support enabled.

Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.  
Easy Clean Wheelchair Cushion  
The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.



PCROH03

No

Cushion ROHO High Profile 8x9

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Cell size 8 x 9  
Dimensions 38 x 43 cm  
Profile High  
Height 100mm  
Max User Weight Unlimited



Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.

PCROH04

No

Cushion ROHO High Profile 9x9

The High Profile Cushions redistribute pressure and conform to the user's body in order to prevent against ulcers and tissue breakdown. Its soft, air-filled cells permit blood flow while simultaneously providing users with the additional support they require while seated. The ROHO High-Profile Cushion assists those who are unable to independently shift their own weight.

Cell size 9 x 9  
Dimensions 43 x 43 cm  
Profile High Height  
100mm Max User  
Weight Unlimited



Available with your choice of either a nylon lycra or incontinent resistant cover. Both are machine washable, and they protect the cushion's surface from wear and tear as well as provide aid during slide transfers. The bottom of the cover has a non-skid coating to prevent sliding while in use.

Easy Clean Wheelchair Cushion

The sturdy neoprene rubber cells resist liquid/moisture and clean easily with soap and water.



PCROQ01

No

Cushion ROHO Quadtro 8 cells  
x 9 cells HP (fits 14x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 8 x 9  
Dimensions 38 x 43 cm  
Profile High  
Height 100mm  
Max User Weight  
Unlimited



PCROQ02

No

Cushion ROHO HP Quadtro 9 cells x 9 cells (fits 16x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 9 x 9  
Dimensions 43 x 43 cm  
Profile High  
Height 100mm  
Max User Weight  
Unlimited



PCROQ03

No

Cushion ROHO Quadtro 10 cells x 9 cells High Profile (fits 18x16 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 10 x 9  
Dimensions 47 x 43 cm  
Profile High  
Height 100mm  
Max User Weight  
Unlimited



PCROQH1

No

Cushion ROHO HP Quadtro  
10x10

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 10 x 10  
Dimensions 47 x 47 cm  
Profile High  
Height 100mm  
Max User Weight  
Unlimited



PCROQH2

No

Cushion HP ROHO Quadtro  
11x10 (fits 20x18 W/C)

The Quadtro Select cushion is made for individuals who require superior pressure relief and protection from soft tissue breakdown. Its high-profile nature facilitates blood flow via advanced shape-conforming technology.

The ROHO Quadtro Select is lightweight, easy to transport, and easy to clean (just use soap and water). It features four individual sections that are uniformly inflated through a single valve. A simple push of the ISOFLO Memory Control knob locks air into each section, therefore filling them to maximum capacity. The ROHO cushion is also available in varying dimensions to accommodate for users' individual seating needs.

Cushion Cover Included

This seat cushion comes with a durable cover that reduces unnecessary friction and shear. Spill something on your seat? Machine wash the cover in cold water for an easy clean up.

Cell size 11 x 10  
Dimensions 52 x 47 cm  
Profile High  
Height 100mm  
Max User Weight  
Unlimited



**Combination Cushions**

PCJC11

No

Cushion Jay 2 Standard  
460mmx410mm (18"x16")

Contoured foam cushion with Jay Flow fluid pad 46 x 41 cm  
18"x16"

Designed for the client with moderate to high postural and pressure care needs

Firm, anatomically correct contoured foam base creates a foundation for postural stability

Provides excellent lateral and forward/rear stability

Fluid Tripad in rear seat well maintains even pressure distribution and keeps fluid under areas at high risk of skin breakdown

Supplied with breathable ballistic stretch cover



PCROHY2

No

Cushion ROHO Hybrid Elite  
460mmx460mm (fits 18x16  
W/C)

Combines ROHOs Air Flotation  
qualities and Jays base technology  
Ideal for clients who require moderate  
to high pressure care and positional  
support needs

Contoured foam base improves  
alignment of the pelvis and lower  
extremities, enhancing sitting posture  
and function

Air cells are interconnected and their  
unique shape promote blood circulation  
to facilitate pressure prevention,  
healing and relief

Flexible and adjustable air cells adapt  
to patient contours and reduce shear  
forces while transferring or  
repositioning

Durable neoprene construction is soft  
to suit skin sensitivity, and is easily  
cleaned

Jay foam base is made from fluid-  
resistant closed cell foam for additional  
stability, assisting in safer wheelchair  
transfers

Two-way stretch cover with encased  
comfort foam is machine washable  
Supplied with pump

Cushions is available with fluid resistant  
cover

Dimensions 46 x 46cm

Height 95mm

Max User Weight 225kg



**Foam Cushions**

PCCU08

Yes

iCare Basic Foam Cushion 400  
x 400 x 50

Range of sizes and thicknesses

- » Zip off washable covers
- » Range of sizes available
- » Pressure and heat sensitive  
ActiveX™ core

400 x 400 x 50

The standard seat cushion range  
features high density visco cores which  
respond to pressure in a different way.

Instead of pushing back like standard  
foam cushions the Icare cushions will  
mould around the shape of your body  
taking away points of pressure.





PCCU11

Yes

iCare Basic Foam Cushion 400  
x 400 x 100


Range of sizes and thicknesses

- » Zip off washable covers
- » Range of sizes available
- » Pressure and heat sensitive  
ActiveX™ core

400 x 400 x 100

The standard seat cushion range  
features high density visco cores which  
respond to pressure in a different way.

Instead of pushing back like standard  
foam cushions the Icare cushions will  
mould around the shape of your body  
taking away points of pressure.



PCCU09

Yes

iCare Basic Foam Cushion 450  
x 450 x 50

Range of sizes and thicknesses

- » Zip off washable covers
- » Range of sizes available
- » Pressure and heat sensitive  
ActiveX™ core

450 x 450 x 50

The standard seat cushion range  
features high density visco cores which  
respond to pressure in a different way.

Instead of pushing back like standard  
foam cushions the Icare cushions will  
mould around the shape of your body  
taking away points of pressure.



PCCU12

Yes

iCare Basic Foam Cushion 450  
x 450 x 50

Range of sizes and thicknesses

- » Zip off washable covers
- » Range of sizes available
- » Pressure and heat sensitive  
ActiveX™ core

450 x 450 x 100

The standard seat cushion range  
features high density visco cores which  
respond to pressure in a different way.

Instead of pushing back like standard  
foam cushions the Icare cushions will  
mould around the shape of your body  
taking away points of pressure.



PCCU10

Yes

iCare Basic Foam Cushion 450  
x 500 x 50

Range of sizes and thicknesses



- » Zip off washable covers
- » Range of sizes available
- » Pressure and heat sensitive  
ActiveX™ core

450 x 500 x 50



The standard seat cushion range features high density visco cores which respond to pressure in a different way.

Instead of pushing back like standard foam cushions the Icare cushions will mould around the shape of your body taking away points of pressure.

<p>PCCU13</p> <p>Yes</p>	<p>iCare Basic Foam Cushion 450 x 500 x 100</p>	<p>Range of sizes and thicknesses</p> <ul style="list-style-type: none"> <li>» Zip off washable covers</li> <li>» Range of sizes available</li> <li>» Pressure and heat sensitive ActiveX™ core</li> </ul> <p>The tandard seat cushion range features high density visco cores which respond to pressure in a different way.</p> <p>Instead of pushing back like standard foam cushions the lcare cushions will mould around the shape of your body taking away points of pressure.</p>	<p>450 x 500 x 100</p>	
<p>PCCU05</p> <p>No</p>	<p>Basic Foam Cushion with Gel Insert. (450 x 400 x 50mm)</p>	<p>This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.</p>	<p>Max user weight: 150Kg 450x400x50mm</p>	

PCCU06

No

Basic Foam Cushion with Gel Insert. (450 x 450 x 50mm)

This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.

Max user weight: 150Kg  
450x450x50mm




PCCU07




No





Basic Foam Cushion with Gel Insert. (500 x 450 x 50mm)

This foam pressure care cushion is suited to people at low to medium risk of pressure injury that is, damage to the skin and underlying tissue as a result of prolonged sitting without moving. Helps relieve chronic pain and pressure for users up to 110kg. Made from a patented gel bladder encased in resilient PU comfort foam cushion for superior weight equalisation and pressure distribution.

Max user weight: 150Kg  
500x450x50mm



<p>PCWCSC</p> <p>Yes</p>	<p>Contoured Wheelchair Seat Cushion (18" x 16")</p>	<p>An anatomically moulded wheelchair cushion that is used to relieve and prevent back pain, numbness and other discomfort that can be experienced after prolonged sitting in a slung seat wheelchair.</p>	<p>18" x 16" 460 x 410 mm</p>	
<p><b>Comfort Products</b></p>				
<p>PCSS01</p> <p>Yes</p>	<p>Sheepskin Medical Grade</p>	<p>Genuine Australian lambswool overlay Provides additional comfort to existing mattress or chair Soft pile allows air circulation to keep skin dry Natural fibre reduces the risk of skin irritation Machine washable Dimensions are calculated on an average</p>	<p>640mm wide x 960mm Long</p>	
<p>PCCFT02</p> <p>Yes</p>	<p>Sheepskin Slippers Small - Pair</p>	<p>Low profile medical sheepskin slipper pressure-redistributing and microclimate controlling properties non-slip suede soles slim and lightweight washable to 80°C</p>		

PCCFT03	Yes	Sheepskin Slippers Medium - Pair	Low profile medical sheepskin slipper pressure-redistributing and microclimate controlling properties non-slip suede soles slim and lightweight washable to 80°C	
PCCFT04	Yes	Sheepskin Slippers Large - Pair	Low profile medical sheepskin slipper pressure-redistributing and microclimate controlling properties non-slip suede soles slim and lightweight washable to 80°C	
PCCFT06	Yes	Sheepskin Heel Protectors	Genuine medical sheepskin Soft pile allows air circulation to keep skin dry Natural fibres reduces the risk of skin irritation Adjustable velcro securing straps Machine washable One size fits all Supplied as a pair.	
PCHBHE	Yes	Heel Protector - Bed Heel Elevator	Wedge-shaped contoured foam provides excellent support for the lower leg Rigid foam aids blood flow and reduces pressure	

---

**Static Mattresses**



PCFM01

No

Mattress - Area Care Plus -  
Static Foam

The AreaCare Plus is available in Standard and King Single size. The three layered, zoned design of this very high risk pressure surface is ideal for patients with poor mobility or uncertain skin sensitivity.

Three-layer foam pressure surface offers maximum pressure reduction in a three-layer foam surface for higher-risk patients.

A top layer of high density, temperature-sensitive, memory foam conforms to the body's shape, working in combination with the higher resistance support layer to disperse pressure evenly. The top layer also incorporates a special immersion, slow recovery foam foot zone for the ulcer prone heel area. The core layer and edges provide form, strength, and stability. Standard model suitable for users up to 220 kilograms; King Single up to 250 kilograms.

Single (1960 x 865 x  
150mm)



PCPRM01

No

Mattress - Area Care Pro Static Foam

The AreaCare Pro multi-density foam pressure surface offers maximum pressure reduction in a three layer foam surface for higher risk patients. The three layered, zoned design of this very high risk pressure surface is ideal for heavier patients with poor mobility or uncertain skin sensitivity. Made from premium materials, each foam layer is made of varying density to maximise the redistribution of pressure. The top layer memory foam also incorporates a special immersion, slow recovery foam foot zone for the ulcer prone heel area, improves ventilation to the skins surface and increases the overall weight capacity to 220 kilograms.

Single (1960 x 865 x 150mm)



PCICSL5

No

iCare Long Single Mattress - SOFT

For clients who like a soft mattress, or who spend a long time in bed. The feel of the IC25 soft Mattress is best suited to aged or petite patients.

The IC25 Mattresses have almost no partner disturbance or roll together.

Long Single (2030 x 900 x 225)



PCICMLS

No

iCare Long Single Mattress - MEDIUM

Middle of the range, the IC20 proves to be popular for most. The medium feel provides a soft feel but excellent support.

Long Single (2030 x 900 x 200)

How pressure care works with the IC20: This IC20 medium mattress is heat and pressure sensitive which means that as it gets warmer it becomes softer. So when there is a pressure point developing which usually gets hot and enflamed, this area of the mattress will become softer which allows the blood to flow freely again.



PCICFLS



No




iCare Long Single Mattress - FIRM

The IC15 firm mattress offers support and comfort. Firm feel makes it easier for turning.

Long Single (2030 x 900 x 150)



<p>PCICSKS</p> <p>No</p>	<p>iCare King Single Mattress - SOFT</p>	<p>For clients who like a soft mattress, or who spend a long time in bed. The feel of the IC25 soft Mattress is best suited to aged or petite patients.</p> <p>The IC25 Mattresses have almost no partner disturbance or roll together.</p>	<p>King Single (2030 x 1070 x 225)</p>	
<p>PCICMKS</p> <p>No</p>	<p>iCare King Single Mattress - MEDIUM</p>	<p>Middle of the range, the IC20 proves to be popular for most. The medium feel provides a soft feel but excellent support.</p> <p>How pressure care works with the IC20: This IC20 medium mattress is heat and pressure sensitive which means that as it gets warmer it becomes softer. So when there is a pressure point developing which usually gets hot and enflamed, this area of the mattress will become softer which allows the blood to flow freely again.</p>	<p>King Single (2030 x 1070 x 200)</p>	

<p>PCICFKS</p> <p>No</p>	<p>iCare King Single Mattress - FIRM</p>	<p>The IC15 firm mattress offers support and comfort. Firm feel makes it easier for turning.</p>	<p>King Single (2030 x 1070 x 150)</p>	
<p><b>Alternating Pressure Reduction Mattresses</b></p>				
<p>PCMR01</p> <p>No</p>	<p>ProCair Mattress Replacement System</p>	<p>ProCair is designed for individuals at the high-risk pressure injury category. Built for optimal patient therapy and comfort, ProCair is now available in Standard and King Single sizes with sealed base to support your infection control strategy.</p>	<p>Weight Capacity (kg) 40-200</p>	
<p>PCMR02</p> <p>No</p>	<p>ProCair Plus Mattress Replacement System</p>	<p>The ProCair Plus is our premium solution for high-risk patients optimal patient therapy and comfort; fully automated operational features; advanced safety features including back-up battery; and advanced error and service indicators. Available in Standard and King Single models.</p>	<p>King Single (2000 x 1050 x 200mm) Single (2000 x 880 x 200mm) Weight Capacity 40-200Kg</p>	

### Pressure Reduction Overlays

PCMO02




No






Procair Mattress Overlay System

Functional overlay for pressure prevention  
Automatic pressure setting  
Used over the existing mattress for basic support  
Easy to operate and maintain  
40-150 kg  
Ideally suited to home care or aged care environments, this product can be used for people at risk to high risk of pressure injury.  
Automatic weight sensing technology that delivers fail-safe patient comfort levels.  
Conform to all relevant Australian standards.




2000 x 880 x 130 mm  
Weight Capacity (kg) 40-150kg



ITEM CODE	CHSP EQUIP IN SCOPE LIST	HCP	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>TRANSFER AIDS</b>						
<b>Handy Bars</b>						
TATHB	Yes	Yes	Transfer Handy Bar	Provides a safe and supportive handhold for transfers in and out of the car Multifunctional device designed to help users get in and out of a car easier.	Maximum User Weight Capacity: 320 kg	
<b>Leg Lifters</b>						
TALL	Yes	Yes	Leg Lifter (Single) Pelican Brand	Rigid device with loop for foot at one end and two straps at the other for assisting independence to lift a leg into a car or bed.	Overall length: 102 cm	
TADHLL	Yes	Yes	Dual Handle Leg Lifter	A single leg lifter with 2 handles for users who cannot lift their leg with one hand.	Length 86.5 cm Weight 110g	
<b>Slide Sheets</b>						

TASLS00	Yes	Yes	Slide Sheet SMALL - 1000mm x 1450mm	Slippery surface both sides, use folded onto itself to provide low resistance under clients requiring assistance with repositioning. Size: 1m x 1.4m sq for bed transfers	1m x 1.4m Colour may vary	
TASLS01	Yes	Yes	Slide Sheet LARGE - 2000mm x 1450mm	Slippery surface both sides, use folded onto itself to provide low resistance under clients requiring assistance with repositioning. Size: 2m x 1.4m sq for bed transfers	2m x 1.4m Colour may vary	
<b>Slide Boards</b>						
TACTP01	Yes	Yes	Car Transfer Pad - Pelican	Car slide transfer pad has 2 slip surfaces that decreases friction allowing the client to slide and pivot into the car seat.	Top Pad size 58cm x 58cm	
TATB01	Yes	Yes	The Slider	"Assists patient transfers between wheelchair and car, chair, bed etc Strong perspex construction Rear wheel cutout, covered with a flexible plastic guard, minimises skin damage while transferring over the arch of the wheel Not to be used as a lifting device"	Length 640mm Width 350mm Weight 1.5kg Max User Weight 120kg	
WCAATB	Yes	Yes	PVC Angled Transfer Board	A smaller transfer board angled at one end to allow better positioning around the wheel of a wheelchair for side transfers. 15833 Angled Transfer Board ( P.V.C.) flexible construction for slight height variation.	SWL 110kg 600mm x250mm x 6mm	



TATB03	Yes	Yes	Curved Transfer Board PVC	PVC transfer board	Weight: 1.95kg. Max. user weight: 152kg	
			<b>Turntables</b>			
TASC01	Yes	Yes	Swivel Cushion Rota	This cushioned seat glides you in and out of cars smoothly without struggling. Placed on any seat it will swivel in both directions and may help those who suffer from arthritis, leg or back pain. It's lightweight, portable and consists of a plastic base with a fleece covered cushion. A strap secures the position of the swivel seat.	Diameter: 40cm Height: 70mm Weight: 1.6kg Maximum user weight: 190kg	
TASC02	Yes	Yes	Patient Swivel Seat	Reduces the risk of strain or back injury for the carer and user Assists the user in transfers from a seated position Ideal for transfers in a variety of situations such as car, wheelchair and household areas Swivels 360? Lightweight and portable Smooth and safe turning action Covered in soft fabric upholstery Nonslip underside	Diameter: 39cm Height: 50mm Weight: 0.9kg Maximum user weight: 100kg	

TASC03

Yes

Yes

Handy Soft Turn Pad

Reduces the risk of strain or back injury for the carer and user

Assists the user in transfers from a seated position

Ideal for transfers in a variety of situations such as car, wheelchair and household areas

Swivels 360°

Lightweight and portable

Cotton upholstered seating surface

Nylon swivel mechanism provides smooth and safe turning action

Nonslip underside





Diameter: 45cm




Height: 40mm



Weight: 0.35kg



Maximum user weight:  
125kg



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>DAILY LIVING AIDS</b>					
<b>Crockery</b>					
DLACU05	Yes	Feeding Cup with Long Spout and Handles	Homecraft Feeding Cup, Standard 8mm Spout, 250ml, 2/pack		
DLAPL01	Yes	Plate Incurve Surround Homecraft	This flexible guard can be fitted to ordinary plates to assist with one-handed eating. Three clips attach it to the rim and create a high inward sloping face to help with food collection and prevent spillage.	Plate size 190 to 254mm external diameter. Dishwasher safe. Melamine material dishwasher safe to 180 degrees F (82°C).	
DLAPL03	Yes	Plate Manoy Large	Manoy Contoured Plate, Large, 279 x 197mm, White		
<b>Cutlery</b>					
DLACT08	Yes	SURE GRIP - Teaspoon	Comfortable ribbed built-up handles help people with weak grasp to hold utensils securely. Stainless steel with rubber handles. Institutional dishwasher safe up to 125 F (52°C)	SURE GRIP Non-bendable Cutlery Dishwasher safe.	

DLACT07	Yes	SURE GRIP - Tablespoon	Comfortable ribbed built-up handles help people with weak grasp to hold utensils securely. Stainless steel with rubber handles. Institutional dishwasher safe up to 125 F (52°C)	SURE GRIP Non-bendable Cutlery Dishwasher safe.	
DLACT05	Yes	SURE GRIP - Fork	Comfortable ribbed built-up handles help people with weak grasp to hold utensils securely. Stainless steel with rubber handles. Institutional dishwasher safe up to 125 F (52°C)	SURE GRIP Non-bendable Cutlery Dishwasher safe.	
DLACT06	Yes	SURE GRIP - Knife	Comfortable ribbed built-up handles help people with weak grasp to hold utensils securely. Stainless steel with rubber handles. Institutional dishwasher safe up to 125 F (52°C)	SURE GRIP Non-bendable Cutlery Dishwasher safe.	

<p><b>DLACT04</b></p> <p>Yes</p>	<p>SURE GRIP BENDABLE - Teaspoon</p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.</p>	<p>SureGrips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.</p>	
<p><b>DLACT03</b></p> <p>Yes</p>	<p>SURE GRIP BENDABLE - Tablespoon</p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.</p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.</p>	

<p><b>DLACT01</b>      <b>Yes</b></p>	<p><b>SURE GRIP BENDABLE - Fork</b></p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle.</p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.</p>	
<p><b>DLACT02</b>      <b>Yes</b></p>	<p><b>SURE GRIP BENDABLE - Knife</b></p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. Please note: The knife is not bendable but features a unique rocking action.</p>	<p>SURE GRIP Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The forks and spoons have a special twist built into the metal shaft that allows them to be bent in any angle for either left or right handed use. Handle length 165mm. Dishwasher safe.</p>	

DLACT18

Yes

GOOD GRIPS BENDABLE -  
Teaspoon

Good Grips Bendable Utensils offer forks and spoons that can be adjusted to different angles. The shaft of the utensils features a special twisting mechanism so that users can customize the reach of their forks and spoons. This allows for more comfortable use that puts less stress on the hand. Note: Rocker Knife is not bendable.

6.5 inch stainless steel bendable utensils that feature 1- $\frac{3}{8}$  inch built up rubber handle  
Latex-free, flexible handle features ribbing for easy and adjustable gripping  
Ideal for individuals that have a limited range of motion, especially with hand to mouth movements  
Available as Fork, Teaspoon, Tablespoon, and Rocker Knife  
Utensils are dishwasher safe  
Handle length 165mm



DLACT26

Yes

GOOD GRIPS BENDABLE -  
Tablespoon

Good Grips Bendable Utensils offer forks and spoons that can be adjusted to different angles. The shaft of the utensils features a special twisting mechanism so that users can customize the reach of their forks and spoons. This allows for more comfortable use that puts less stress on the hand. Note: Rocker Knife is not bendable.

6.5 inch stainless steel bendable utensils that feature 1- $\frac{3}{8}$  inch built up rubber handle  
Latex-free, flexible handle features ribbing for easy and adjustable gripping  
Ideal for individuals that have a limited range of motion, especially with hand to mouth movements  
Available as Fork, Teaspoon, Tablespoon, and Rocker Knife  
Utensils are dishwasher safe  
Handle length 165mm





DLACT25

Yes

GOOD GRIPS BENDABLE -  
Fork

Good Grips Bendable Utensils offer forks and spoons that can be adjusted to different angles. The shaft of the utensils features a special twisting mechanism so that users can customize the reach of their forks and spoons. This allows for more comfortable use that puts less stress on the hand. Note: Rocker Knife is not bendable.

6.5 inch stainless steel bendable utensils that feature 1- $\frac{3}{8}$  inch built up rubber handle  
Latex-free, flexible handle features ribbing for easy and adjustable gripping  
Ideal for individuals that have a limited range of motion, especially with hand to mouth movements  
Available as Fork, Teaspoon, Tablespoon, and Rocker Knife  
Utensils are dishwasher safe  
Handle length 165mm



DLACT27

Yes

GOOD GRIPS NON BENDABLE - Knife

Good Grips Bendable Utensils offer forks and spoons that can be adjusted to different angles. The shaft of the utensils features a special twisting mechanism so that users can customize the reach of their forks and spoons. This allows for more comfortable use that puts less stress on the hand. Note: Rocker Knife is not bendable.

6.5 inch stainless steel bendable utensils that feature 1-3/8 inch built up rubber handle  
Latex-free, flexible handle features ribbing for easy and adjustable gripping  
Ideal for individuals that have a limited range of motion, especially with hand to mouth movements  
Available as Fork, Teaspoon, Tablespoon, and Rocker Knife  
Utensils are dishwasher safe  
Handle length 165mm



DLACT21




Yes





GOOD GRIPS WEIGHTED and BENDABLE - Teaspoon





Good Grips Weighted & Bendable Utensils have 170g of added weight in the built-up handle. This provides more control so that the utensil reaches the mouth easier. The utensil head is bendable. The utensils are ideal for persons with limited hand control, Parkinsons disease, or spasticity.





Available as Fork, Teaspoon, Tablespoon, Souper Spoon, and Rocker Knife  
Built-up handle  
170g  
Institutional dishwasher safe



DLACT21	Yes	GOOD GRIPS WEIGHTED and BENDABLE - Tablespoon	Good Grips Weighted & Bendable Utensils have 170g of added weight in the built-up handle. This provides more control so that the utensil reaches the mouth easier. The utensil head is bendable. The utensils are ideal for persons with limited hand control, Parkinsons disease, or spasticity.	Available as Fork, Teaspoon, Tablespoon, Souper Spoon, and Rocker Knife Built-up handle 170g Institutional dishwasher safe	
DLACT21	Yes	GOOD GRIPS WEIGHTED and BENDABLE -Souper spoon	Good Grips Weighted & Bendable Utensils have 170g of added weight in the built-up handle. This provides more control so that the utensil reaches the mouth easier. The utensil head is bendable. The utensils are ideal for persons with limited hand control, Parkinsons disease, or spasticity.	Available as Fork, Teaspoon, Tablespoon, Souper Spoon, and Rocker Knife Built-up handle 170g Institutional dishwasher safe	
DLACT21	Yes	GOOD GRIPS WEIGHTED and BENDABLE -Fork	Good Grips Weighted & Bendable Utensils have 170g of added weight in the built-up handle. This provides more control so that the utensil reaches the mouth easier. The utensil head is bendable. The utensils are ideal for persons with limited hand control, Parkinsons disease, or spasticity.	Available as Fork, Teaspoon, Tablespoon, Souper Spoon, and Rocker Knife Built-up handle 170g Institutional dishwasher safe	

DLACT22	Yes	GOODGRIPS Knife - Weighted	<p>Good Grips Cutlery Bendable Utensils are made of stainless steel featuring a 35mm diameter built up handle. The Weighted Utensils range are available with 170g of extra weight in the built up handle to provide more control for people with a tremor or limited hand control, Parkinson's disease or spasticity. The knife is weighted but not bendable.</p>	<p>Available as Fork, Teaspoon, Tablespoon, Souper Spoon, and Rocker Knife Built-up handle 170g Institutional dishwasher safe</p>	
<b>Dressing Aids</b>					
DLADA01	Yes	Button Hook Modular Kings Complete	<p>Assists one-handed dressing Stainless steel wire passes through button hole, gripping the button. It is pulled through with a slight twisting motion Contoured handle</p>	<p>Length 230mm Weight 70g</p>	
DLADA02	Yes	Dressing Stick	<p>Wooden stick with epoxy coated hook at one end and rubber tip at the other. Assists in applying and removing clothing. Can also be used to push or pull objects.</p>	<p>66cm length 5/8"-diameter wooden dowel. Finished with lacquer with reinforced hook.</p>	
DLADA03	Yes	Shoe Horn Long Handled Plastic	<p>Plastic handle designed to assist in putting on shoes especially for people who have difficulty bending and reaching.</p>	<p>Length 60cm. Blue Colour. Includes wrist strap.</p>	

<p>DLADA04</p> <p>Yes</p>	<p>Stocking Aid</p>	<p>Stocking aid allows the user to don socks and stockings independently. Side notches hold the sock securely, while the strap ends help pull stocking into position efficiently. Built with a flexible plastic gutter for improved maneuverability.</p>	<p>Available in contoured or notched. Straps 36". Gutter length 9".</p>	
<p>DLADA05</p> <p>Yes</p>	<p>Sock Aid - Single Plastic</p>	<p>Simple aid designed to assist with putting on a stocking especially if bending and reaching the foot is difficult. Cut out side notches designed to hold stocking 1 metre length pull straps</p>		
<p>DLADA07</p> <p>Yes</p>	<p>Stocking Compression Applicator Ezy-As Small</p>	<p>The unique design of the Ezy-As compression garment applicator makes it possible to assist in the fitting of a wide range of compression garments.</p>	<p>Small 28 - 35 cm</p>	
<p>DLADA08</p> <p>Yes</p>	<p>Stocking Compression Applicator Ezy-As Medium</p>	<p>The unique design of the Ezy-As compression garment applicator makes it possible to assist in the fitting of a wide range of compression garments.</p>	<p>Medium 32 - 40cm</p>	

DLADA09	Yes	Stocking Compression Applicator Ezy-As Large	The unique design of the Ezy-As compression garment applicator makes it possible to assist in the fitting of a wide range of compression garments.	Large 38cm - 47cm	
DLADA10	Yes	Compression Stocking / Garment Applicator - Small		Small	
DLADA11	Yes	Sock Aid - Terry Towelling	Simple aid designed to assist with putting on a sock especially if bending and reaching the foot is difficult. Terry towelling outer fabric for greater hold with smooth nylon inner surface to reduce friction.	Nylon pull straps - 700mm	
<b>Household Aids</b>					
DLAHH3	Yes	Tap Turner - Crystal	This Crystal Tap Turner will assist you to use the tap on any tap that does not have an 'X' style knob. The Turner clamps onto the tap and then you can use the handle to control the tap. This is a fantastic tool for those who struggle to grasp a tap and close their hand. Instead, this tool gives you a larger area handle to grip on to.	Length 190mm. Weight 112g	

DLAHHA4

Yes

Tap Turner - Red/Blue

The Homecraft Tapturn Turners are strong moulded turners which simply hook over the top of a 'X' type tap, to make it easy to turn. They are supplied as a red and blue pair for hot and cold taps.

Length 135mm  
Width 47mm



DLAHHA8

Yes

Non Slip Matting




A roll of versatile soft material that can be used on any surface to help stop items from moving, or to protect against scratches.

Comes in a roll (size 1000 mm long by 300 mm wide) and can be cut to size with household scissors.

This non-slip matting has a slightly spongy texture which prevents items moving or slipping on almost any surface.

Place a piece of matting on your table or bench, under trays or when preparing foods - you can even use it under your floor mat to help stop it sliding. Also great to prevent unwanted scratches on any surface. Try wrapping some around a jar lids to improve your grip when opening. The unique plastic material won't fray and and simply wipes clean with a damp cloth - or throw it in the washing machine.



		<b>Tray Mobile</b>			
DLATM01	Yes	Trolley Walker - Adjustable Height	<p>Includes clip-on plastic trays that can be removed easily for cleaning.</p> <p>Can be used with or without trays or it can be used as a table.</p> <p>Rails on tray edges prevent items from slipping off.</p> <p>Lockable hand brakes mean this trolley can also be used as a walking aid.</p>	<p>Overall Width: 585mm</p> <p>Overall Depth: 560mm</p> <p>Handle Height Range: 755mm - 905mm</p> <p>Wheel Diameter: 150mm (6")</p> <p>Maximum User Weight: 150kg</p>	
DLATM02	Yes	Traymobile (no brakes)	<p>Ideal for indoor use</p> <p>Two removable trays to assist transporting smaller household items</p> <p>Height adjustable</p> <p>Four swivel castors provide greater manoeuvrability</p>	<p>Frame Height 800-950mm</p> <p>Frame Width 490mm</p> <p>Frame Depth 540mm</p> <p>Max User Weight 125kg</p>	
		<b>Kitchen Aids</b>			
DLAKA13	Yes	Plastic Spreading Board	no info available		



DLAKA12

Yes

Bread Board

Ideal for people with tremors or impaired sight, the clever design of this cutting aid makes kitchen work much easier. Food is firmly held in the device while a series of rubber stopper feet keep the board from moving around on your kitchen bench or table. Great for bread, cheese, vegetables or other food items, your knife slides between the two guides, making it easy to cut straight and safe, with no risk of injury to fingers or hands. The board can be adjusted for a range of different cutting thickness



DLAKA14

Yes

ETAC Food Prep System

This durable cutting board is ideal for those who have difficulty in gripping, or only have the use of one hand. The clamp can fix items such as tins or bowls, and the stainless steel spikes keep vegetables still when peeling or cutting. It has four non-slip rubber feet to prevent the board from moving. Suction pads are also provided for use on wet surfaces.

Size: 315mm x 295mm.  
Weight: 1.66kg.



DLAKA15

Yes

Cordless Kettle Tipper

Suits most cordless kettles.

Weight: 1kg

There is a metal plate that secures the base of the cordless kettle and a frame with a Velcro strap to keep the upper part of the kettle in place.

A cut out at the foot of the tipper allows cups to be positioned as close to the kettle as required.

The frame has a built-in back stop to prevent spillage if the kettle is accidentally released during use.



DLAKA16

Yes

Uccello Kettle Tipper

The Uccello Tipping Kettle was designed using the innovative and unique PowerPour technology which means the user can safely pour boiling water without removing the kettle from its cradle by just tilting the kettle.

1.5L capacity



Its ergonomic handle makes it comfortable to use for people with arthritis, grip weakness, or hand injury. The kettle features a wide spout and large hinged lid to easily fill with water. The water level indicator assists the user to ensure it isn't overfilled. The base is non-slip and weighted.

The Kettle has a stainless steel heating element, power-on light and auto shut-off with overheating protection.

This product is available in Black & White or Red & White.

DLAKA05




Yes




Jar & Bottle Opener - Mighty Lever




The Jar and Bottle Opener has a flexible stainless steel band which is fitted around a lid, and then clamped tight using the screw at the end of the handle. The handle can then be used


Suitable for tops 20-110mm diameter. Weight 88g.



DLAKA06	Yes	Jar Twister	A cone shaped rubber moulding with fluted finger grips on the outside, and ribbing on the inside to give a good grip for opening twist-off jar lids.	Suitable for tops 20 to 85mm diameter. Weight: 100g	
DLAKA10	Yes	Zyliss 5 in 1 Can Opener	This versatile opener is the ideal gadget for opening pull tabs on cans, plastic and metal twist-off caps, metal bottle caps, pop-tops and jar lids. No more struggling to get the caps and lids off bottles and jars, this handy gadget has all your needs covered. The soft touch grip makes it comfortable and simple to use.	Dimensions: 12Lx15Wx4H	
<b>Personal Care Aids</b>					
DLAPA05	Yes	Long Handled Toe Wiper	Helps people who have limited access to their feet toes to wash and dry between their toes.	Long bendable wire stem with a plastic built-up easy-grip handle Flat wire frame covered in a white towelling pad for washing between the toes Includes 2 pads Length: 710mm	

<p>DLAPA09</p> <p>Yes</p>	<p>Reaching Aid Economy - 67cm</p>	<p>Trigger style mechanism pick up stick Slip resistant jaw provides a secure hold Hook trigger Magnetic tip</p>	<p>67cm</p>	
<p>DLAPA10</p> <p>Yes</p>	<p>Reaching Aid Economy - 82cm</p>	<p>Trigger style mechanism pick up stick Slip resistant jaw provides a secure hold Hook trigger Magnetic tip</p>	<p>82cm</p>	
<p>DLAPA11</p> <p>Yes</p>	<p>Reaching Aid Squeeze - short</p>	<p>Slip resistant rotating jaw designed to pick up objects at different angles Simple hand trigger operation Magnetic tip Wheelchair/walking frame clip supplied.</p>	<p>10204/65 (65 cm)</p>	

<p>DLAPA12</p> <p>Yes</p>	<p>Reaching Aid Squeeze - long</p>	<p>Slip resistant rotating jaw designed to pick up objects at different angles Simple hand trigger operation Magnetic tip Wheelchair/walking frame clip supplied.</p>	<p>10204/81 (81 cm)</p>	
<p>DLAPA13</p> <p>Yes</p>	<p>Bottom Wiper</p>	<p>The Homecraft Toilet Aid provides assistance with independent hygiene. The aid is designed with a curved metal handle and a molded-plastic head that gives the user the proper reach. The handle's recessed serrations allow it to securely grip paper.</p>	<p>Length 25.4cm. Weight 71g.</p>	
<p>DLAPA14</p> <p>Yes</p>	<p>Long Handled Sponge</p>	<p>As you get older, it can be difficult to bend down to wash your feet or to lift them up so you can wash them. Doing so might even put you at risk for a slip and fall in the shower. This long-handled bath sponge has a wire that is strong, but also easily bendable to be rotated, twisted, and moved in any position that helps you accomplish optimal cleaning without bending! It can be used by seniors and people with disabilities to bathe independently, by caretakers to easily clean residents, or for household cleaning tasks.</p>	<p>Length 24". Sponge head 3.75" x 5".</p>	

ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>LIFTERS AND HOISTS</b>					
LHEH01	No	Patient Lifter - 180kg	<p data-bbox="696 331 781 355" style="color: red;"><b>Lifters</b></p> <p>The LiftAid 180 is a lightweight, easy to operate mobile lifter that allows patients to be lifted and transferred from floor, bed or chair with ease and safety.</p> <ul style="list-style-type: none"> <li>3 height settings - tool-less adjustment</li> <li>Quick release hook system</li> <li>Electrically operated base</li> <li>Handheld controller</li> <li>Emergency stop</li> <li>Manual emergency lowering</li> <li>Optional Pivot Frame</li> <li>Compatible with range of slings</li> </ul>	<ul style="list-style-type: none"> <li>Width: 61.5cm</li> <li>Height: 43cm - 181cm (3 height settings)</li> <li>Length: 109cm</li> <li>Product Weight: 35.3kg</li> <li>Warranty: Lifter &amp; Actuator - 2 years   Battery - 1 year</li> <li>S.W.L: 180kg</li> <li>Max Spread: 101cm</li> <li>Underbed Clearance: 12.5cm</li> </ul>	

LHEH02

No

Patient Lifter - 320kg

The LiftAid 320 is a heavy-duty mobile lifter ideal for the majority of common lifting situations.

Patients can be lifted and transferred from floor, bed or chair. Fold-down armrests assist with gait training during rehabilitation or for general toileting.

Adjustable height frame allows overall mast height be set to match environmental and lifting needs.

Features the advanced Linak® Jumbo Care Diagnostic Control System for intelligent self-diagnostic system capabilities. Low lifting height and soft start/stop function,

3 height settings with tool-less adjustment

quick release hook system for easy sling bar installation

swivel-style front and rear castors; foot-operated brake locks on rear

3 height settings with tool-less adjustment

Quick release hook system for easy sling bar installation

Swivel style front and rear castors; foot operated brake locks on rear

Electrically operated 'U' shaped base with wide flanged front legs for enhanced stability and access

Length Overall - 1410 mm

Width Overall 7950 mm - 1190 mm

Height Overall 2180 - mm

Height of Lift 650 mm - 1840 mm

Clearance Under Bed 16 - mm

Diameter of Front Wheel - 105 mm

Diameter of Rear Wheel - 127 mm

Load Capacity - 320 kg



**Lifter/Hoist Slings**



LHSLGPS

No

Sling General Purpose - Sml

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair.

Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports

Suits the majority of lifters on the market

Brightly coloured handles and adjusters secured to sling for easy access and visibility

Quick-drying mesh fabric machine washable up to 85°C for microbiological control

Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

840 x 1160 x 490 mm ( W x L x D )  
Capacity - 205 kg



LHSLGPM

No

Sling General Purpose - Med

An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair.

Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports

Suits the majority of lifters on the market

Brightly coloured handles and adjusters secured to sling for easy access and visibility

Quick-drying mesh fabric machine washable up to 85°C for microbiological control

Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

890 x 1200 x 545 mm (W x L x D) Capacity - 205kg



LHSLGPL



No



Sling General Purpose - Lge



An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports Suits the majority of lifters on the market Brightly coloured handles and adjusters secured to sling for easy access and visibility Quick-drying mesh fabric machine washable up to 85°C for microbiological control Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort

1135 x 1325 x 640 mm  
(W x L x D)  
Capacity - 205kg



<p>LHSLGPXL</p> <p>No</p>	<p>Sling General Purpose - Xlge</p>	<p>An easy fitting general purpose sling providing good support for patients with basic upper body control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. 3 loop support strap system controls the positioning of the head and leg supports Suits the majority of lifters on the market Brightly coloured handles and adjusters secured to sling for easy access and visibility Quick-drying mesh fabric machine washable up to 85°C for microbiological control Heavy-duty, 3 mm quilted cushioning foam in leg sections for greater comfort</p>	
<p>LHSLGPHS</p> <p>No</p>	<p>Sling General Purpose Head Support Sml</p>	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	<p>900, 230, 450, 790, 965 mm, 205 kg</p> 

<p>LHSLGPHM</p> <p>No</p>	<p>Sling General Purpose Head Support Med</p>	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	<p>1010, 230, 490, 890, 1060 mm, 205 kg</p>	
<p>LHSLGPHL</p> <p>No</p>	<p>Sling General Purpose Head Support Lge</p>	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	<p>1080, 250, 480, 990, 1180 mm, 300 kg</p>	

<p>LHSLGPHXL</p> <p>No</p>	<p>Sling General Purpose Head Support XLge</p>	<p>General Purpose with Head Support is an easy fitting general purpose sling ideally suited for patients with less control. Suits general lifting from the floor, bed or chair. Versatile enough for general patient transfers, toileting or hygiene care. Suitable for users up to 205 kilograms (Small and Medium) or 300 kilograms (Large, Extra Large and Extra, Extra Large).</p> <p>For Yoke frame lifters only</p>	
<p>LHSLHS</p> <p>No</p>	<p>Sling Hygiene Sml</p>	<p>An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access  buckle waist strap for added security and increased infection control.  Generous padding under arms  3 loop support strap system controls positioning of head and leg supports.  Suits majority of lifters on the market.  brightly coloured handles and adjusters, secured to sling for easy access and visibility  quick drying mesh fabric; machine washable up to 85°C for microbiological control  manufactured from the highest quality material.</p>	

LHSLHM

No

Sling Hygiene Med

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access  
buckle waist strap for added security and increased infection control.  
Generous padding under arms  
3 loop support strap system controls positioning of head and leg supports.  
Suits majority of lifters on the market.  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
manufactured from the highest quality material.

Maximum user weight  
205kg



LHSLHL

No

Sling Hygiene Lge

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access  
buckle waist strap for added security and increased infection control.  
Generous padding under arms  
3 loop support strap system controls positioning of head and leg supports.  
Suits majority of lifters on the market.  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
manufactured from the highest quality material.

Maximum user weight  
300kg





LHSLHXL

No

Sling Hygiene XLge

An easy fitting sling specifically designed for toilet transfer with open area from mid back to mid thigh. open area from mid back to mid thigh for hygiene access  
buckle waist strap for added security and increased infection control.  
Generous padding under arms  
3 loop support strap system controls positioning of head and leg supports.  
Suits majority of lifters on the market.  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
manufactured from the highest quality material.

Maximum user weight  
300kg



LHSLHHS

No

Sling Hygiene With Head Support Sml

Open area from mid back to mid thigh for hygiene access  
inbuilt head support for patients with less control  
buckle waist strap for added security and increased infection control  
generous padding under arms  
3 loop support strap system controls positioning of head and leg supports  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
suits majority of lifters on the market  
manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)

Maximum user weight 205kg



LHSLHHM

No

Sling Hygiene With Head Support Med

Open area from mid back to mid thigh for hygiene access  
inbuilt head support for patients with less control  
buckle waist strap for added security and increased infection control  
generous padding under arms  
3 loop support strap system controls positioning of head and leg supports  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
suits majority of lifters on the market  
manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)

Maximum user weight 205kg



LHSLHHL



No




Sling Hygiene With Head Support Lge

Open area from mid back to mid thigh for hygiene access  
inbuilt head support for patients with less control  
buckle waist strap for added security and increased infection control  
generous padding under arms  
3 loop support strap system controls positioning of head and leg supports  
brightly coloured handles and adjusters, secured to sling for easy access and visibility  
quick drying mesh fabric; machine washable up to 85°C for microbiological control  
suits majority of lifters on the market  
manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)

Maximum user weight 300kg



<p>LHSLHHXL</p> <p>No</p>	<p>Sling Hygiene With Head Support XLge</p>	<p>Open area from mid back to mid thigh for hygiene access          inbuilt head support for patients with less control          buckle waist strap for added security and increased infection control          generous padding under arms          3 loop support strap system controls positioning of head and leg supports          brightly coloured handles and adjusters, secured to sling for easy access and visibility          quick drying mesh fabric; machine washable up to 85°C for microbiological control          suits majority of lifters on the market          manufactured from the highest quality material, in compliance with industry standards (AS/NZS ISO 10535:2011)</p>	<p>Maximum user weight 300kg</p>	
<p><b>Stand Aid Slings</b></p>				
<p>LHSLSATM</p> <p>No</p>	<p>Sling - Standaid Torso Medium</p>	<p>Fully padded, including soft padding under the arms          Various attachments          Built-in Keyhole Plate Safety Straps (Applicable for the version of this product which uses Keyhole Plates)          A comfortable, well padded, high back Stand Up Sling with a strong wide</p>	<p>Chest Size 85cm - 110cm          Load capacity 200kg</p>	

ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>MOBILITY AID ACCESSORIES</b>					
<b>IV POLES</b>					
MAIVMO	Yes	IV Stand Mobile	2 Hook design 560mm nylon base Chromed steel pole 50mm castors	Adjustable height 1240- 2130mm	
<b>ROLLATOR ACCESSORIES</b>					
MAARSK	Yes	Rollator Ski Replacement (x2)	Pair of Rollator Skis to replace lost or damaged skis.	To fit 25mm Rollator frames	
<b>OXYGEN CARRIER</b>					
MAAC02	Yes	Oxygen Bottle Carrier For Wheelchair	Oxygen bottle carrier frame to attach to wheelchair. Contact store for availabilty or suitability to client's chair.		

MAOC4WW

Yes

Oxygen Carrier for 4 Wheel Walker

Wire frame designed to fit under the seat of the Rollator Deluxe range.

Suits oxygen bottles up to 125 mm diameter (C size)

Allows the under seat bag or basket to remain on Rollator

Angled for optimal oxygen supply & security



MABF01

Yes


Fracture Board

**STUMP SUPPORTS**


Designed to support an extended leg and is normally used in a wheelchair. The Board can be used for either the left or right leg. Made from a strong plastic that is easy to wipe clean, to help Infection Control. This plastic is virtually maintenance free and will not splinter. The webbing straps can be removed for washing or replacing if necessary.

Length 112cm  
Width 40cm







ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
REEXPED	Yes	<b>Exercise Equipment</b> Exercise Pedals	Lightweight and portable, the versatile product can be used for a variety of alternative exercises, for both lower and upper body. The compact design makes it easy to move the exerciser from room to room, to take outside in the fresh air or even to take when travelling. The exerciser is designed to be placed on the floor when exercising with the feet, or placed on a table top when exercising with the hands. It features both forwards and backwards movement, to ensure you are working all muscle groups. It also offers an easily adjustable resistance level - simply turn the knob to either increase or decrease the resistance at will.	This product weighs 2.5 kilograms and measures 280 mm high, 445 mm wide and 530 mm deep.	






<p>RESQBP</p> <p>Yes</p>	<p>Balance Pad</p>	<p>Balance Pad ultrasoft speciality foam is both water- and dirt-resistant</p> <p>Application: Balance training, co-ordination and reaction training, mobility training, standing stability training, functional retraining of the muscles in the lower extremities, motor-skill training to maintain balance, e.g. for age related physiological changes to the muscles and skeletal frame. The Physiomed Balance Pad - several can be laid on top of each other to increase the "wobble effect". The closed-cell robust and hygienic. They are multi-functional in use and are ideally suited for use as swimming floats.</p>	<p>Length 50cm</p> <p>41cm</p> <p>0.7kg</p>	<p>Width</p> <p>Weight</p>	
<p><b>Pain Management</b></p>					

<p>RETU01</p> <p>Yes</p>	<p>Tens Unit Spectrum</p>	<p>This easy to use, digital TENS machine offers natural, non-invasive pain relief without the use of chemicals or drugs. Transcutaneous Electrical Nerve Stimulation (TENS) delivers tiny electrical impulses to the nerve fibres that lie underneath the skin's surface, providing pain relief by blocking pain signals to the brain via the spinal cord and peripheral nervous system.</p> <p>Large liquid crystal display and easy operation press buttons, ideal for people who have difficulty reading smaller panel inscriptions or turning smaller dials. The unit offers 32 stimulation modes for pain relief and muscle training. There are 5 adjustable TENS modes, 3 adjustable EMS modes and 24 professional presets for TENS and EMS available.</p>	<p>1 x TENS machine unit  1 x D size 9-volt battery  1 x AC power adaptor  2 x electrode leads  4 x self adhesive electrodes (pack of 4)  1 x instruction manual  1 x hard plastic carrying case</p>
<p><b>Weights</b></p> <p>REHW0.5</p> <p>Yes</p> <p>HAND WEIGHT Plastic COATED .5KG</p>		<p>These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.</p>	<p>Dumbbell, Vinyl Coated, Purple 0.5kg</p>



REHW1.0	Yes	HAND WEIGHT Plastic COATED 1.0KG	These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.	Dumbell, Vinyl Coated, Pink 1.0 kg	
REHW1.5	Yes	HAND WEIGHT Plastic COATED 1.5KG	These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.	Dumbell, Vinyl Coated, Green 1.5kg	
REHW2.0	Yes	HAND WEIGHT Plastic COATED 2.0KG	These dumbbells are made of steel but coated in a soft vinyl (PVC) easy-clean coating. Can be used during work-outs to increase strength training, or used alone as part of a weight training program.	Dumbell, Vinyl Coated, Mauve 2.0kg	
		<b>Exercise Media</b>			
RETBTHN	Yes	Theraband Yellow 5.5m	TheraBand Resistance Exercise Bands are best used before working out for preventative stretches that prime muscles for strenuous exercises. They are an essential part of any workout or athletic activity.	Yellow thin 5.5m	

<p>RETPXSF</p> <p>Yes</p>	<p>Theraputty Tan Extra Soft 85g</p>	<p>Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems. It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.</p>	<p>Tan Extra Soft 85g</p>	
<p>RETPSFT</p> <p>Yes</p>	<p>Theraputty Yellow Soft 85g</p>	<p>Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems. It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.</p>	<p>Yellow Soft 85g</p>	
<p>RETPMDM</p> <p>Yes</p>	<p>Theraputty Red Medium 85g</p>	<p>Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems. It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.</p>	<p>Red Medium 85g</p>	

RETPFRM

Yes

Theraputty Green Firm 85g

Metron Therapeutic Putty works to strengthen your hands and can help prevent future hand and wrist problems. It is also perfect for hand grip exercises, stress relief exercises, wrist exercises, and upper arm exercises. It is superior to squeeze balls, hand weights, hand springs, or any other exercise/stress relief product.

Green Firm 85g



**Weights**

REAW0.9


Yes

Ankle Weight - Vinyl 0.9kg  
(This is a substitute for the 1.0kg ankle weight which is currently not available)

Vinyl Cuff Weights with a long velcro closure strap assures a secure and comfortable fit even during the most strenuous exercise program.

White 0.9kg Ankle Weight with velcro fastening.



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>BRACES AND SUPPORTS</b>					
<b>BRACES</b>					
BRSKB01	Yes	Knee Brace - Futuro Sport Adjustable	<p>Walk. Jog. Run. The FUTURO Sport Adjustable Knee Support keeps you moving with compression, strength and support where you need it. It also helps provide protection on the track, court, or green. With this adjustable knee wrap, there is no excuse not to be active anymore.</p> <p>Special Features: Adjustable straps for customised fit and support Use For: Sprains, strains, arthritis Product Type: Wrap Around Features and Benefits:</p> <p>Open knee cap design for added support Side ventilation for breathability and reduced bunching behind knee Durable neoprene-blend material for strong support, soothing warmth, and comfort Easy to put on and take off</p>	One Size Fits Most, Measure knee circumference: 33cm - 45.7cm	

BRSW01

Yes

Wrist Braces Roylan D Ring Small  
Left

Wrist Circumference  
14.6cm to 16.5cm  
Length 16.5cm



The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.

---

BRSW02

Yes

Wrist Braces Roylan D Ring Small  
Right

Wrist Circumference  
14.6cm to 16.5cm  
Length 16.5cm



The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.

---



BRSW03

Yes

Wrist Braces Roylan D Ring  
Medium Left

Wrist Circumference  
17.1cm to 19.1cm Length  
17.8cm



The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.

BRSW04

Wrist Braces Roylan D Ring  
Medium Right

Wrist Circumference  
17.1cm to 19.1cm Length  
17.8cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW05

Yes

Wrist Braces Roylan D Ring Large  
Left

Wrist Circumference  
19.7cm to 21.6 cm  
Length 19cm

The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.



BRSW06

Yes

Wrist Braces Roylan D Ring Large  
Right

Wrist Circumference  
19.7cm to 21.6 cm  
Length 19cm



The right choice for high-compliance treatment of arthritis, carpal tunnel syndrome, tendonitis and other wrist disorders. Prevents unwanted flexion or extension, yet allows unobstructed motion of the fingers, thumb. Non-elastic circumferential design for easy application and for firm, evenly distributed pressure without restricting circulation. Protective padding lined with cotton stockinette (except for white) D-ring closure system for simple application, removal and customized tension control. A double layer of stockinette under closures cushions straps. Forearm length is proportional to brace size. Removable, adjustable aluminum bar for additional support. Machine wash and air dry. Latex free.

**ORTHOSES**

---

BRSAF01a

Yes

Rebound Foot Up Ankle Cuff S/M

Rebound Foot-Up is a lightweight ankle foot orthosis that offers dynamic and discreet support for drop foot or related disorders requiring dorsiflexion support

Anke Cuff 18-25cm  
Foot wrap, strap and bolt 20-28cm

Contoured to the shape of the leg, it is comfortable to wear and easy to use

Rebound Foot-Up provides visible improvement in gait by supporting the foot the moment it is raised



BRSAF02a

Yes

Rebound Foot Up Ankle Cuff L/XL

Rebound Foot-Up is a lightweight ankle foot orthosis that offers dynamic and discreet support for drop foot or related disorders requiring dorsiflexion support

Anke Cuff 26-33cm Foot wrap, strap and bolt 20-28cm

Contoured to the shape of the leg, it is comfortable to wear and easy to use

Rebound Foot-Up provides visible improvement in gait by supporting the foot the moment it is raised



BRSAF06

Yes

Ankle Foot Orthosis Foot-Up  
Shoeless Wrap Universal







Rebound Foot-Up Foot Wrap can be used when barefoot or in open shoes or sandals. The adjustable elastic strap of the Rebound Foot Up is attached to the Foot Wrap using a nut and bolt fastened through the eyelets.

One size cut to size (20-28cm)


The Rebound Foot-Up Ankle Cuff and Rebound Foot-Up Foot Wrap are made from breathable material that can be worn comfortably for long periods of time.

The Rebound Foot-Up is hand-washable.



ITEM CODE	CHSP EQUIP IN SCOPE LIST	PRODUCT	PRODUCT DESCRIPTION	SPECIFICATIONS	IMAGE
<b>HOME ACCESS</b>					
<b>Portable Ramps</b>					
HARMP1	Yes	Multipurpose Ramp 870mm	Fibreglass ramp that has a concertina folding action for storage and transport. 875 mm length.	835mm width 875mm length 5.6kg	
HARMP2	Yes	DecPac Portable Ramp 700mm	Fibreglass ramp that has a concertina folding action for storage and transport. 700mm length.	745mm width 700mm length 4.7kg	
HARMP3	Yes	DecPac Portable Ramp 1200mm	Fibreglass ramp that has a concertina folding action for storage and transport. 1200mm length.	750mm width 1200mm length 8.0kg	
HARMP3E	Yes	DecPac Portable Ramp 1200mm (with EDGE)	Fibreglass ramp that has a concertina folding action for storage and transport. 1200mm length with EDGE.	820mm width 1200mm length 9.7kg	
HARMP5	Yes	DecPac Portable Ramp 1650mm	Fibreglass ramp that has a concertina folding action for storage and transport. 1650mm length.	750mm width 1650mm length 10.8kg	
HARMP5E	Yes	DecPac Portable Ramp 1650mm (with EDGE)	Fibreglass ramp that has a concertina folding action for storage and transport. 1650mm length with EDGE.	820mm width 1650mm length 13kg	



<p><b>HARMP6</b>      <b>Yes</b></p>	<p>DecPac Portable Ramp 2000mm</p>	<p>Fibreglass ramp that has a concertina folding action for storage and transport. 2000mm length.</p>	<p>750mm width 2000mm length 11.7kg</p>	
<p><b>HARMP7</b>      <b>Yes</b></p>	<p>DecPac Portable Ramp 2000mm (with EDGE)</p>	<p>Fibreglass ramp that has a concertina folding action for storage and transport. 2000mm length with EDGE.</p>	<p>820mm width 2000mm length 13.7kg</p>	