

# Shower Chair



## Purpose


To provide a safe seated position whilst showering.

## Features

- Plastic moulded seat which is comfortable and hygienic
- Arm rests assist with getting on and off the chair
- The height of the shower chair can be adjusted to make it easier to sit and stand
- Non slip feet

## Installation

- Adjust equipment to the appropriate height by depressing the silver button on the leg and sliding the leg pieces up and down until the button clicks into the appropriate hole. Repeat on all legs. Ensure all four legs are in contact with the floor to avoid rocking or unsteadiness




## Maintenance

- The frame should be checked monthly for signs of rust, cracks, changes in the appearance of the rubber stoppers, or damage/bending in the metal. If any of these are present, stop use and contact your supplier
- This product should be cleaned regularly with mild detergent unless otherwise stated by the manufacturer

## Safe Usage

- Use two hands when transferring on and off the shower chair to avoid tipping
- This equipment is for sitting only. Do not stand on this equipment
- You should only use the equipment in the environment in which it is prescribed
- 125kg weight limit unless otherwise stated by the manufacturer
- Ensure the silver buttons on the height adjustable legs have slotted into the correct holes before sitting



These instructions are to be used in conjunction with the information that your Health Professional has discussed with you.

P 1800 684 277

E [info@maxhealthcare.com.au](mailto:info@maxhealthcare.com.au)

W [maxhealthcareequipment.com.au](http://maxhealthcareequipment.com.au)

