by **REDGUM**





B580

FEATURES

- Allows for toileting away from the toilet for when users cannot make it to the bathroom
- Soft padded seat / backrest / armrests
- Removable padded seat overlay
- Easy to clean vinyl upholstery
- Adjustable seat height
- Flip up seat to access bowl
- Removable bowl with carry handle and lid
- Non slip rubber tips
- Lightweight easy to move

SPECIFICATIONS All dimensions shown are nominal

610mm (W) x 505mm (D) x 785mm (H)

(lowest height setting)

430mm (W) x 450mm (D) x 410 - 560mm (H)

(7 height settings at 25mm increments)

Backrest:

420mm (W) x 240mm (H) x 380mm

(from seat to top of backrest)

Armrests:

40mm (W) x 215mm (D) x 210mm (H)

(from seat to top of armrest)

510mm between armrests

Product Weight: 7kg

Maximum User Weight: 159kg

Please contact REDGUM Brand if you require further information on this product

Removable Bowl (with carry handle) + Lid

P: 08 9248 4180

E: orders@for-de.com.au

1 Business way Malaga Western Australia 6090 Part of the For-de Group Pty Ltd

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

redgumbrand.com.au