

# Lightweight Rollators



# For who you are & what you need.

At any age, we need to feel good and stay active. Knowing that we can get up and go anywhere, at any time, is vital to your independence and enjoying life to its full potential.

With this in mind, we developed our Airgo® rollators. Offering more than just support, they provide convenience, comfort, stability and security. We take great care to create distinct

models that offer unique features so that you can choose one that's right for you.

The joy of walking is that it is a safe, easy and effective exercise.

An Airgo® rollator can help you improve your health and may even increase your energy



### Our exclusive features make a difference!



### Handles

Only Airgo® features Ergo-form Handgrips that fit much more comfortably and place your hand at a natural angle for better brake access. They allow you to reach and grip the brakes without shifting, allowing for better control and added safety.

Feather-touch brakes offer easy manoeuvrability.

Available on: eXcursion™, Fusion™, Lightweight Rollator, Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



### Backrests

Our large Contoured Backrest provides full back support thanks to its unique profile design. It comfortably encompasses your entire back while providing improved support for your posture.

Available on: eXcursion™, Fusion™ and Vista™.



The Ergo-Back is simply the most comfortable rigid rollator backrest available today. It features an ergonomically shaped soft padded cushion made from high density foam.

Available on: Comfort-Plus and Comfort-Plus XWD.



### Seats

The Ergo-Seat's ergonomic shape and high density foam cushion with laminated memory foam on top provides unrivalled seating comfort. Its unique design distributes weight effectively across the entire surface, reducing average pressure measurements by as much as 70% over standard rollator seats.

Available on: Comfort-Plus and Comfort-Plus XWD.





Ergo-Seat <sup>™</sup> Normal Hard Seat Pressure map with Piezo resistive technology



### Footvests

Our new Swing n' Lock Footrests secure to the frame when folded, allowing for a hassle free experience. Converting from a rollator to transport chair has never been easier!

Available on: Fusion™.



### Basket

Soft Basket (underseat) that does not need to be removed when folding the rollator. With convenient compartments and zipper closure, it is opaque and located under the seat for stability, keeping valuables organized and accessible to the user yet safely out of sight.

Available on: Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



Our covered Soft Basket (side-fold) with reflective stitching doesn't need to be removed when folding your rollator. In addition, it is machine washable and dryable for easy cleaning.

Available on: eXcursion™. Fusion™ and Vista™.



### Wheels

Comfort-grip™ Indoor/Outdoor Wheels are non-marring and offer better grip on carpet along with superior braking action. Solid mags for easy maintenance.

Available on: Comfort-Plus and Comfort-Plus XWD.



Our new Low-Profile Comfort-Grip™ Wheels feature a sport tire design that offers better stability, with low rolling resistance on smooth indoor surfaces. Once outdoors, the outer treads grip in to provide superior control.

Available on: eXcursion™ and Fusion™



### Brakes

Our new Tool-Free Brake Adjustment makes adjusting your brake cables a breeze. The large knob is senior friendly and FAQ videos are available on our website.

Available on: eXcursion™. Fusion™ and Vista™.

### Height options

We're not all created the same. Our Multi-Height Frames allows comfort and safety to fit even more people by offering our rollators at the seat height that's right for you.

Available on: Adventure 8, Comfort-Plus, Comfort-Plus XWD, Ultra-Light 6, eXcursion™ and Fusion™.





#### **Accessories**



Oxygen Bag 703-120

Accommodates size C and M6 oxygen tanks.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus and Comfort-Plus XWD.



Saddle Bag 703-100

Snaps conveniently onto your rollator's backrest.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus, Comfort-Plus XWD and Ultra-Light 6.



Side-fold Cane Holder

730-946

Allows you to securely and easily attach a cane to your side-folding rollator.

Fits on: eXcursion™ and Fusion™.



Cane Holder 730-947

Allows you to securely and easily attach a cane to your rollator.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus and Comfort-Plus XWD.



Side-fold Cup Holder

730-944

Stay hydrated while walking with this convenient rotatable cup holder.

Fits on: eXcursion™ and Fusion™



Apollo LED Safety Light

730-950

Rechargeable LED light featuring 3 usage modes, and 20 Lumens of power.

Fits on: All rollators.



Side-fold Curb Climber

A01-551

Makes it easier to lift your side-folding rollator over obstacles such as sidewalks.

Fits on: eXcursion and Fusion.



Curb Climber 703-050

Makes it easier to lift your rollator over

obstacles such as sidewalks.

Fits on: Lightweight Rollator, Adventure 8, Comfort-Plus, and Comfort-Plus XWD.

### There is an Airgo rollator for your needs!

	Colour	Ergo-Form Handgrips	Ergo-Seat and Back	Contoured Backrest	Basket
Lightweight	Ooloui	Tidilagripo	and Baok	Buomoot	Wire
Adventure8	•	1			Soft closed, underseat
Ultra-Light6	•	<b>✓</b>			Soft closed, underseat
Comfort-Plus Comfort-Plus XWD	•	<b>✓</b>	<b>✓</b>		Soft closed, underseat extra large
eXcursion X23	•				
eXcursion X20	• •				Covered,
eXcursion X18	• •	<b>/</b>		<b>/</b>	soft, in front of rollator
eXcursion XWD	• •				Tollatol
eXcursion6	•				
Fusion F23	•				Carrand
Fusion F20	•				Covered, soft, in front of rollator
Fusion F18	•				
Vista	• •			<b>/</b>	Covered, soft, in front

Seat Adjustment	Ultra-light Frame	Wheels	Transport Feature	MAX. weight capacity	NET weight	
Fixed 22" (56 cm)		Standard 8" (20 cm)		300 lb (136 kg)	16.2 lb (7.3 kg)	
Dual height		Standard 8" (20 cm)		300 lb (136 kg)	16.2 lb (7.3 kg)	
Multi-height	<b>✓</b>	Standard 6" (15 cm)		250 lb (113 kg)	10.5 lb (7.8 kg)	
Dual baiabt		Comfort-Grip		300 lb (136 kg)	17.2 lb (7.8 kg)	
Duai neight	Dual height	8" (20 cm)		400 lb (181 kg)	17.8 lb (8 kg)	
Fixed 23" (58 cm)				300 lb (136 kg)	17.1 lb (7.7 kg)	
Fixed 20" (51 cm)		Low-Profile Comfort-Grip 8" (20 cm)		300 lb (136 kg)	15.4 lb (6.9 kg)	
Fixed 18" (46 cm)	<b>/</b>			300 lb (136 kg)	15.4 lb (6.9 kg)	
Fixed 20" (51 cm)				400 lb (181 kg)	16.4 lb (7.4 kg)	
Fixed 20" (51 cm)		Standard 6" (15 cm)		300 lb (136 kg)	13.5 lb (6 kg)	
Fixed 20" (51 cm)				300 lb (136 kg)	17.1 lb	
Fixed 18" (46 cm)	<b>/</b>	Low-Profile Comfort-Grip 8" (20 cm)	Full 2 in 1	300 lb (136 kg)	(7.7 kg)	
Fixed 23" (58 cm)		,		300 lb (136 kg)	18.9 lb (8.6 kg)	
Fixed 20" (51 cm)	<b>/</b>	Low-Profile 8" (20 cm)		300 lb (136 kg)	15 lb (6.8 kg)	

## eXcursion\*





Ergo-form handgrips



Ultra-light frame



Easy to fold



Recognized by the CAOT

### The ultra-light side-fold with superior comfort & stability!

- Full contoured backrest provides superior back support (3).
- Fold by simply pulling the strap (4), the soft covered basket folds with the frame (5).
- A folding lock keeps the unit closed for secure lifting & storing.
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (1).
- Comfort-grip wheels positioned outside the frame for increased stability (6).
- NEW tool-free easy to adjust brake system.
- Weight capacity of 300 lb / XWD weight capacity of 400 lb.


	<b>!</b>	<b>/</b>	<u></u>	H		<b>*</b>		
X23 (23" seat)	5'9" - 6'11"	33.5" - 39.5"	23"	17.5"	25"	28.5"		
Colour: Pacific I	Blue (700-917)	Net weight (with	nout acc.) 17.1 lb					
X20 (20" seat)	5" - 6'2"	30.5" - 37"	20"	17.5"	25"	28.5"		
	Colours: Pearl Black (700-920) Pacific Blue (700-921) Cranberry (700-922) Lime (700-923) Net weight (without acc.) 15.4 lb							
X18 (18" seat)	4'3" - 5'5"	29.5" - 35"	18"	17.5"	25"	28.5"		
Colours: Pacific Blue (700-926) Pearl Black (700-927) Net weight (without acc.) 15.4 lb								
XWD	5" - 6'2"	30" - 36.5"	20"	19.5"	27"	28.5"		
Colours: Pacific Blue (700-928) Pearl Black (700-929) Net weight (without acc.) 16.4 lb								

## eXcursion\*(\*)



Our ultra-light excursion frame, now even lighter with 6" wheels!

- Full contoured backrest provides superior back support (3).
- Fold by simply pulling the strap (4), the soft covered basket folds with the frame (5).
- A folding lock keeps the unit closed for secure lifting & storing.
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (1).
- 6" wheels are non-marring with low rolling resistance, and ideal for indoor use. (6).
- NEW tool-free easy to adjust brake system.
- · Weight capacity of 300 lb.



Patents 7828305, 2587989, 200710103964.9 • ADP type 2 AA0001355 • AADL C225.



### Converts effortlessly from a side-folding rollator to a transport chair.

- Padded seat and full contoured backrest (1).
- Soft basket doesn't need to be removed when folding (2).
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (3).
- Easy to activate brakes (4).
- Locks in open position and when folded for safer lifting.
- 8" non-marring wheels are ideal for indoor/outdoor use (7).
- Weight capacity: 300 lb (136 kg).
- Seat belt available separately, item 700-863.

	<b>†</b> ‡	<b>/</b>	₹1	<b>+</b>		<b>*</b>		
Fusion F23	5'9" - 6.11"	33.5" - 39.5"	23"	19.5"	27"	28.5"		
Colour: Pearl B	Colour: Pearl Black (700-932) Net weight (without acc.) 19.9 lb							
Fusion F20	5' - 6'2"	30" - 36.5"	20"	19.5"	27"	28.5"		
Colour: Pacific Blue (700-935) Net weight (without acc.) 17.1 lb								
Fusion F18	4'3" - 5.5"	29.5" - 35"	18"	19.5"	27"	28.5"		
Colour: Cranberry (700-938) Net weight (without acc.) 17.1 lb								


Patents 7828305, 2587989, 200710103964.9, 7540527 and patents pending. ADP type 3 AA0000G94.





Our travel-friendly side-folding rollator, with pop-in wheels and a lightweight design.

- Padded seat and adjustable ergo backrest (1).
- Extra-large soft basket with a secure cover, that folds with the frame (2).
- Adjustable ergonomic handles, with integrated parking brakes (3).
- Brakes lock for security when seated (4).
- Easy to fold, and stands when folded. Locks both open and closed for safety.
- Travel friendly pop-in wheels allow it to be stored almost anywhere (5).
- Treaded 8" front wheels for superior steering and rolling comfort (6).



Colours: Rocket red (700-906) Royal blue (700-907) Onyx black (700-908) Stellar silver (700-909) Net weight (without acc.) 15 lb


ADP Type 3 AA0001338, AADL C225.

## Ultra-Light 6



### The lightest rollator available: a breeze to lift up and take anywhere!

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 6" non-marring wheels are ideal for indoor use (6).
- Reflective trim for added visibility (7).
- Weight capacity: 250 lb (113 kg).




Patents 7828305, 2587989, 200710103964.9 ADP type 2 MW2GM0050 AADL C225.

### Lightweight Rollator



An economical choice with all the basic features.

- · Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Removable wire basket (4).
- 8" non-marring wheels are ideal for indoor/outdoor use (5).
- Reflective trim for added visibility (6).
- Weight capacity: 300 lb (136 kg).



Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9 AADL C225.

## Adventure<sup>™</sup>8









#### Light and dependable!

- Easy-to activate brakes (2).
- · Padded seat and backrest (3).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 8" non-marring wheels (6).
- Reflective trim for added visibility (7).
- Weight capacity: 300 lb (136 kg).

	<b>!</b>	<b>/</b>	<u></u>	H	1	<b>*</b>	
Standard	5' - 6'2"	33" - 39"	21.5"	17"	24.5"	27.5"	
HEMI	4'.3" - 6'2"	30" - 36"	18.5"	17"	24.5"	25.5"	
Colours: Titanium (700-942) Net weight (without acc.) 16.2 lb							

Colours: Titanium (700-942) Net weight (without acc.) 16.2 lb

Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9 ADP type 3 MW3GM0060 (Adventure 8) AADL C225.

# COMFORT-PLUS™ and COMFORT-PLUS™





Ergo-form handgrips



Dual height frame technology



Ergo-Seat & Back



Extra-Large shopping basket

#### Simply the most comfortable rollator on the market.

- Easy-to activate brakes (2).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Comfort-Grip indoor/outdoor 8" wheels (6).
- Reflective trim for added visibility (7).
- Comfort-PlusTM weight capacity: 300 lb (136 kg).
   Comfort-Plus XWDTM weight capacity: 400 lb (181 kg).

	<b>!</b>	<b>/</b>	<b>1</b>	<b>+</b>		<b>*</b>	
Standard	5' - 6'2"	33" - 39"	22 .5"	17"	24.5"	27.5"	
HEMI	4'3" - 6'2"	30" - 36"	19.5"	17"	24.5"	25.5"	
Colour: Iridescent Blue (700-946) Net weight (without acc.) 17.2 lb							
XWD	5' - 6'2"	33" - 39"	22 .5"	19"	26.5"	27.5"	
	4'3" - 6'2"	30" - 36"	19.5"	19"	26.5"	25.5"	
Colour: Iridescent Blue (700-949) Net weight (without acc.) 17.8 lb							

Patents 7066484, 2471692, 7828305, 2587989, 200710103964.9 ADP type 3 MW3GM0050 (Comfort-PlusTM) AADL C225. ADP type 3 MW3GM0055 (Comfort-PlusTM XWD) AADL C064.

### **Real Life Engineering**

Your Airgo® Rollator is guaranteed to be free from defects in materials and workmanship. This commitment demands that our products are engineered to last, but more importantly, it allows you the confidence and freedom to get up and go, anywhere, anytime!

We also appreciate your feedback so we can learn how to continuously improve our rollators. We call this "Real Life Engineering".



www.caot.ca/proudsupporter



www.airgomobility.com
Scan to visit our webpage featuring
our NEW educational videos.



www.walker-facts.com
Please visit to learn more about rollators & walkers.





Drive Medical Canada, Inc. 256 Aviva Park Drive, Unit #2 Vaughan, Ontario L4L 9C7

To learn more, call: (866) 752-1363 or visit us at

www.airgomobility.com

P/N 915-213 rev 11 0516