


**NEW PRODUCT RELEASE !**



300mm  
Grab Rail  
shown

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.



**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

## SUCTION CUP GRAB RAILS

RG5903: 300mm Grab Rail (White)

RG5904: 400mm Grab Rail (White)

### FEATURES

- Assist individual to position / balance themselves.
- Visual indicators to assist with rail installation.
- Improve safety and reduce the risk of falls at home.
- Can be attached to any flat glass / tiled surface.
- No fixings required.
- Quick and easy to fit / reposition when required.
- Ergonomic ribbed hand grip on rail.
- Strong lightweight construction.

### SPECIFICATIONS

Overall Dimensions:

RG5903: 300mm Grab Rail

**410mm (L) x 110mm (W) x 85mm (H)**

Hand Grip Length:

**120mm** (45mm gap from wall to inside of rail)

Product Weight:

**0.4kg**

Overall Dimensions:

RG5904: 400mm Grab Rail

**510mm (L) x 110mm (W) x 85mm (H)**

Hand Grip Length:

**220mm** (45mm gap from wall to inside of rail)

Product Weight:

**0.6kg**

Maximum ASSIST load (both rails):

**65kg**

Product Materials:

**Contact Redgum brand for details**

Please contact **REDGUM Brand** if you require further information on this product.

**P: 08 9248 4180**

**F: 08 9248 9458**

**E: [orders@shoprider.com.au](mailto:orders@shoprider.com.au)**

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

**[redgumbrand.com.au](http://redgumbrand.com.au)**