

**STRONG LIGHTWEIGHT
CONSTRUCTION!**

**HEIGHT
ADJUSTABLE**
contour padded kneepad



**LARGE
PNEUMATIC
WHEELS**
for excellent mobility

Folds down for
**TRANSPORT
AND STORAGE**



KNEE WALKER

RG30KWB

FEATURES

- Height adjustable contour padded kneepad
- Height adjustable handlebars
 - Anti slide hand grips
 - Left + Right hand dual locking brake system
- Large pneumatic wheels for excellent mobility
- Strong lightweight construction
- Easy to assemble
- Removable wire carry basket
- Folds down for storage and transport
- Walking stick holders
- Offset knee rest position - can be left or right

SPECIFICATIONS

(All dimensions shown are nominal and at lowest setting)

- Overall: **450mm (W) x 950mm (L)**
790mm - 950mm (H) (from ground to top of the handlebar)
- Kneepad Size: **370mm (L) x 185mm (W)**
- Kneepad Height: **400mm - 515mm**
(from ground to top of the kneepad)
- Handlebar: **450mm (W) x 790mm - 950mm (H)**
(from ground to top of the handlebar)
- Wheel Diameter: **200mm/8"**
- Product Weight: **13kg**
- Max. User Weight: **130kg**

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180
E: orders@for-de.com.au

ASK FOR THE BRAND

redgumbrandmobility.com.au