

**STRONG LIGHTWEIGHT  
CONSTRUCTION!**

**HEIGHT  
ADJUSTABLE**  
contour padded kneepad



**Baricare**

# HEAVY DUTY KNEE WALKER

RG40KWHD

## FEATURES

- Height adjustable contour padded kneepad
- Height adjustable handlebars
- Ergonomic rubber hand grips
- Large pneumatic wheels for excellent mobility
- Dual hand brakes
- Strong lightweight construction
- Easy to assemble
- Removable wire basket
- Folds down for transport and storage
- Walking stick holders

## SPECIFICATIONS

(All dimensions shown are nominal and at lowest setting)

- Overall: **500mm (W) x 1050mm (L)**  
**890mm - 1040mm (H)** (from ground to top of the handlebar)
- Kneepad Size: **380mm (L) x 190mm (W)**
- Kneepad Height: **510mm - 620mm**  
(from ground to top of the kneepad)
- Handlebar: **490mm (W) x 890mm - 1040mm (H)**  
(from ground to top of the handlebar)
- Wheel Diameter: **200mm/8"**
- Product Weight: **15.5kg**
- Max. User Weight: **180kg**

Please contact **REDGUM Brand** if you require further information on this product

Folds down for  
**TRANSPORT  
AND STORAGE**



**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

**ASK FOR THE BRAND**

**[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)**