

# BARIATRIC SIDE FOLDING SEAT WALKER

ADJUSTABLE  
Handle Height



STANDS  
When  
Folded



WIDER SEAT  
FOR MAX COMFORT



Please contact **REDGUM Brand** if you require further information on this product

RG441BDG - Dark Grey

## FEATURES

- Wider seat for maximum comfort
- 200kg weight capacity
- Ergonomic hand grips to reduce pressure
- Brake and park brakes
- Padded sling seat
- Removable storage bag included
- Height adjustable handles and backrest
- Large flat free wheels for excellent mobility
- Folds up for easy storage
- Strong lightweight construction

## SPECIFICATIONS

(All dimensions shown are nominal and are in mm)

Overall Dimensions (at lowest setting):

**760mm(W) x 690mm(L) x 935mm(H)**

Seat: **560mm(W) x 420mm(L) x 580mm(H)**

Folded: **320mm(W) x 690mm(L) x 935mm(H)**

Bag: **495mm(W) x 100mm(L) x 300mm(H)**

**5kg capacity**

Handle Height: **910mm - 1010mm (5 settings)**

Wheels: **8"/200mm (front), 10"/250mm (rear)**

Product weight: **10kg**

Max. User Weight: **200kg**

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)